

RNAscope® 4-plex Ancillary Kit for Multiplex Fluorescent Reagent Kit v2 Technical Note

Introduction

This Technical Note provides instructions for performing 4-plex *in situ* hybridization (ISH) on pretreated formalin-fixed paraffin-embedded (FFPE) tissue sections, fresh frozen tissues, and other sample types using the RNAscope® Multiplex Fluorescent Kit v2 (Cat. No. 323100) and 4-Plex Ancillary Kit [Cat. No. 323120]. PerkinElmer Opal™ fluorophores and multiplexed biomarker imaging systems (Vectra® or Mantra™) are required for detection of fluorescent signals. For detailed sample preparation

procedures and safety guidelines, refer to the *RNAscope® Multiplex Fluorescent Reagent Kit v2 User Manual* (Doc. No. 323100-USM). Consult our Technical Notes available at www.acdbio.com/technical-support/user-manuals to prepare other sample types. For every chemical, read the Material Safety Data Sheet (MSDS) and follow handling instructions. Wear appropriate protective eyewear, clothing, and gloves. For the latest services and support information, go to: www.acdbio.com/support.

Workflow

Part 1: Prepare Tissue Samples

Prepare your samples following the instructions for sample preparation and pretreatment in the *RNAscope® Multiplex Fluorescent Reagent Kit v2 User Manual* (Doc. No. 323100-USM), available at www.acdbio.com/technical-support/user-manuals.

Part 2: Run the RNAscope® Assay

Prepare Materials

1. Warm 50X Wash Buffer for **10–20 MIN** to remove any precipitation.
2. Prepare 3 L of 1X Wash Buffer by adding 2.94 L distilled water and 1 bottle (60 mL) of 50X Wash Buffer to a large carboy. Mix well.

Prepare Probes

1. Warm probes for **10 MIN** at **40°C**, then cool to **RT**.
2. Briefly spin the C2, C3, and C4 probes.
3. Pipette 1 volume of C2, 1 volume of C3, and 1 volume of C4 probes to 50 volumes of C1 probe into a tube. Invert the tube several times to mix.

Note: Do not mix probes of the same channel. Store mixed probes at **2–8°C** for up to six months.

Hybridize Probe

1. Remove excess liquid from the slides, place in the HybEZ™ or EZ-Batch™ Slide Rack, and add 4–6 drops of the probe mix to entirely cover each slide.
2. Insert slide rack containing the slides into the HybEZ™ Oven for **2 HRS** at **40°C**.
3. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh 1X Wash Buffer.

OPTIONAL STOPPING POINT. The slides can be stored in 5X SSC overnight at **RT** (not provided in the kit).

For the following steps, use reagents from the RNAscope® Multiplex Fluorescent Kit v2 (Cat. No. 323100).

Hybridize Amp 1

1. Remove excess liquid from slides, place in the HybEZ™ or EZ-Batch™ Slide Rack, and add 4–6 drops RNAscope® Multiplex FL v2 Amp 1 to entirely cover each slide.
2. Insert slides into the HybEZ™ Oven for **30 MIN** at **40°C**.

3. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh 1X Wash Buffer.

Hybridize Amp 2

1. Remove excess liquid from slides, place in the HybEZ™ or EZ-Batch™ Slide Rack, and add 4–6 drops RNAscope® Multiplex FL v2 Amp 2 to entirely cover each slide.
2. Insert slides into the HybEZ™ Oven for **30 MIN** at **40°C**.
3. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.

Hybridize Amp 3

1. Remove excess liquid from slides, place in the HybEZ™ or EZ-Batch™ Slide Rack, and add 4–6 drops RNAscope® Multiplex FL v2 Amp 3 to entirely cover each slide.
2. Insert slides into the HybEZ™ Oven for **15 MIN** at **40°C**.

Note: Prepare TSA Plus Fluorophores during this step. See the following section.

3. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.

Prepare Opal™ fluorophores

1. Determine the volume of Opal™ fluorophore needed (150–200 µL per slide).
2. Dilute the Opal™ fluorophore stocks using the RNAscope® Multiplex TSA buffer provided in the RNAscope® Multiplex Fluorescent Kit v2. Follow these recommendations:

Opal™ fluorophore	PerkinElmer Reagent Kit	Recommended dilution range*
Opal 520	FP1487001KT: Opal 520 Reagent Pack	1:750–1:3000
Opal 570	FP1488001KT: Opal 570 Reagent Pack	1:750–1:3000
Opal 620	FP1495001KT: Opal 620 Reagent Pack	1:750–1:3000
Opal 690	FP1497001KT: Opal 690 Reagent Pack	1:750–1:3000

*Start with a dilution of 1:1500 and adjust based on signal intensity.

Note: Keep the diluted Opal™ fluorophore in the dark prior to applying to slides.

Develop HRP-C1 Signal

1. Remove excess liquid from slides, place in the HybEZ™ or EZ-Batch™ Slide Rack, and add 4–6 drops RNAscope® Multiplex FL v2 HRP-C1 to entirely cover each slide.

2. Insert slides into the HybEZ™ Oven for **15 MIN** at **40°C**.
3. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.
4. Remove excess liquid from slides, place in the HybEZ™ or EZ-Batch™ Slide Rack, and add 150–200 µL diluted Opal 520 to each slide, and incubate for **30 MIN** at **40°C**.

Note: You can mix and match channels and fluorophores. For example, you may assign Opal 570 to the C1 channel instead of Opal 520. Do not assign the same fluorophore to more than one channel.

5. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.
6. Remove excess liquid from slides, place in the HybEZ™ or EZ-Batch™ Slide Rack, and add 4–6 drops RNAscope® Multiplex FL v2 HRP blocker to entirely cover each slide.
7. Insert slides into the HybEZ™ Oven for for **15 MIN** at **40°C**.
8. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.

Develop HRP- C2 Signal

1. Remove excess liquid from slides, place in the HybEZ™ Slide Rack or EZ-Batch™ Slide Rack, and add 4–6 drops RNAscope® Multiplex FL v2 HRP-C2 to entirely cover each slide.
2. Insert slides into the HybEZ™ Oven for **15 MIN** at **40°C**.
3. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.
4. Remove excess liquid from slides, place in the HybEZ™ or EZ-Batch™ Slide Rack, and add 150–200 µL diluted Opal 570 to each slide, and incubate for **30 MIN** at **40°C**.

Note: You can mix and match channels and fluorophores. For example, you may assign Opal 620 to the C2 channel instead of Opal 570. Do not assign the same fluorophore to more than one channel.

5. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.
6. Remove excess liquid from slides, place in the HybEZ™ or EZ-Batch™ Slide Rack, and add 4–6 drops RNAscope® Multiplex FL v2 HRP blocker to entirely cover each slide.
7. Insert slides into the HybEZ™ Oven for **15 MIN** at **40°C**.

- Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.

Develop HRP-C3 Signal

- Remove excess liquid from slides, place in the HybEZ™ Slide Rack or EZ-Batch™ Slide Rack, and add 4–6 drops RNAscope® Multiplex FL v2 HRP-C3 to entirely cover each slide.
- Insert slides into the HybEZ™ Oven for **15 MIN** at **40°C**.
- Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.
- Remove excess liquid from slides, place in the HybEZ™ or EZ-Batch™ Slide Rack, and add 150–200 µL diluted Opal 620 to each slide, and incubate for **30 MIN** at **40°C**.

Note: You can mix and match channels and fluorophores. Do not assign the same fluorophore to more than one channel.

- Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.
- Remove excess liquid from slides, place in the HybEZ™ or EZ-Batch™ Slide Rack, and add 4–6 drops RNAscope® Multiplex FL v2 HRP blocker to entirely cover each slide.
- Insert slides into the HybEZ™ Oven for **15 MIN** at **40°C**.
- Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.

For the following steps, use reagents from 4-Plex Ancillary Kit [Cat. No. 323120].

Develop HRP-C4 Signal

- Remove excess liquid from slides, place in the HybEZ™ Slide Rack or EZ-Batch™ Slide Rack, and add 4–6 drops RNAscope® Multiplex FL v2 HRP-C4 from 4-Plex Ancillary Kit For Fluorescent Multiplex v2 (Cat. No. 323120) to entirely cover each slide.
- Insert slides into the HybEZ™ Oven for **15 MIN** at **40°C**.
- Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.
- Remove excess liquid from slides, place in the HybEZ™ or EZ-Batch™ Slide Rack, and add 150–200 µL diluted Opal 690 to each slide, and incubate for **30 MIN** at **40°C**.

Note: You can mix and match channels and fluorophores. Do not assign the same fluorophore to more than one channel.

- Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.
- Remove excess liquid from slides, place in the HybEZ™ or EZ-Batch™ Slide Rack, and add 4–6 drops RNAscope® Multiplex FL v2 HRP blocker from 4-Plex Ancillary Kit For Fluorescent Multiplex v2 (Cat. No. 323120) to entirely cover each slide.
- Insert slides into the HybEZ™ Oven for **15 MIN** at **40°C**.
- Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.

Counterstain and Mount the Slides

Note: Do this procedure with no more than five slides at a time.

- Remove excess liquid from slides, and add ~4 drops of DAPI to each slide.
- Incubate for **30 SEC** at **RT**.
- Remove DAPI and *immediately* place 1–2 drops of Prolong Gold antifade mounting medium on the slide (not provided in the kit).
- Carefully place a 24 mm x 50 mm glass coverslip over the tissue section. Avoid trapping air bubbles.
- Dry slides **30 MIN** to **OVERNIGHT** in the dark.
- Store slides in the dark at **2–8°C**.

Note: Image the slides after eight hours or within two weeks.

Imaging Slides

For imaging using multiplexed biomarker imaging systems Vectra® or Mantra™, refer to the guidelines from Perkin Elmer. The following table lists the corresponding filter settings for each fluorophore:

Opal™ fluorophore	Filter setting
Opal 520	FITC
Opal 570	Cy3
Opal 620	Texas Red
Opal 690	Cy5

For Research Use Only. Not For Diagnostic Use.

Advanced Cell Diagnostics, Inc. and/or its affiliate(s) warrant their products as set forth in the ACD General Terms and Conditions of Sale found on the ACD website at <http://www.acdbio.com/store/terms>. Advanced Cell Diagnostics, Inc. reserves the right to change its products and services at any time to incorporate technological developments. This document is subject to change without notice. Although this document has been prepared with every precaution to ensure accuracy, Advanced Cell Diagnostics, Inc. assumes no liability for any errors, omissions, or for any damages resulting from the use of this information.

© 2017 Advanced Cell Diagnostics. All rights reserved. RNAscope® and HybEZ™ are trademarks of Advanced Cell Diagnostics, Inc. All other trademarks belong to their respective owners.

Headquarters

7707 Gateway Blvd Suite 200, Newark, CA 94545 Phone 1-510-576-8800 Toll Free 1-877-576-3636

For support, email support@acdbio.com

www.acdbio.com