

# Preparing Fresh Frozen Tissue for the RNAscope® and BaseScope™ Assays

## Introduction

This Technical Note provides guidelines to prepare OCT embedded fresh frozen tissue for the RNAscope® 2.5 and BaseScope™ Assays. Wear appropriate

protective eyewear, clothing, and gloves. For the latest services and support information, go to: [www.acdbio.com/support](http://www.acdbio.com/support).

## Workflow

### Part 1: Frozen Tissue Processing

1. Harvest animal tissues within **5 MIN** of sacrifice. Cut dissected tissues into pieces  $\leq 5$  mm in thickness, and keep on ice.
2. Select an appropriately sized cryomold to leave at least 5 mm of space between the tissue and the block edge.
3. Place at least 1 lb of dry ice in a styrofoam box.
4. Create a flat surface in the dry ice to stabilize the cryomold. For example, use a metal block and let the block temperature equilibrate to that of dry ice for **10 MIN**.  
Alternatively, pour isopentane into the container until the level of liquid is slightly above the layer of dry ice. Let the bubbling subside before use.
5. Place tissues on a paper towel to remove excess liquid.
6. Add one drop of OCT to the bottom of the mold, and place the tissue in the desired orientation on top. Tissue will be sectioned starting from the bottom of the mold.
7. Carefully add more OCT until the tissue is completely covered and the cryomold is filled. Remove any bubbles that form near the tissue with a forceps.
8. Place the mold on the flat surface you created in the dry ice, or in cooled isopentane.

9. Let the OCT solidify completely. Keep the mold on dry ice for another **5 MIN**. Wrap frozen blocks with labeled foil, and store at **-70°C**. You may also proceed directly to frozen tissue sectioning.

### Part 2: Frozen Tissue Sectioning

1. Set cryostat temperature to **-15°C** to **-20°C**. Temperature is tissue dependent.
2. Remove and discard the microtome blade. Wipe down the knife holder and anti-roll plate with 100% ethanol.
3. Install a new microtome blade.
4. Place a slide box on dry ice, or inside the cryostat.
5. Set the cutting thickness to 10  $\mu\text{m}$  (+/- 5  $\mu\text{m}$ ).
6. Transfer the frozen tissue from the freezer to the cryostat, transporting on dry ice if necessary.
7. Remove the tissue block from the plastic mold, and equilibrate it to the temperature of the cryostat **15-30 MIN**.
8. Mount the tissue block onto the specimen holder using OCT. Allow the OCT to completely solidify.
9. Section the block, and mount the sections onto room temperature SuperFrost® Plus slides. Immediately place each slide into the slide box.

---

**IMPORTANT!** You must use SuperFrost® Slides (Fisher Scientific, Cat # 12-550-15). Do not allow the slides to come to room temperature. Keep cold at all times.

---

## TECHNICAL NOTE: Fresh Frozen Sample Preparation

10. Allow the slides to dry for a minimum of **1 HR** in the slide box before proceeding to the RNAscope® or BaseScope™ assay.

---

**IMPORTANT!** Do not dry slides at room temperature.

---

**NOTE:** You may store the slides at **-70°C** for up to three months. Use the slides within three months to ensure RNA quality.

---

### *Obtaining Support*

For the latest services and support information, go to:

<https://acdbio.com/technical-support/support-overview>.

At the website, you can:

- Access telephone and fax numbers to contact Technical Support and Sales
- Search through FAQs
- Submit a question directly to Technical Support

NOTICE TO PURCHASER: PLEASE REFER TO THE RNASCOPE® OR BASESCOPE™ USER MANUAL FOR LIMITED USE LABEL LICENSE OR DISCLAIMER INFORMATION. Advanced Cell Diagnostics, Inc. reserves the right to change its products and services at any time to incorporate technological developments. This manual is subject to change without notice. Although this manual has been prepared with every precaution to ensure accuracy, Advanced Cell Diagnostics, Inc. assumes no liability for any errors, omissions, or for any damages resulting from the use of this information.

© 2017 Advanced Cell Diagnostics. All rights reserved. RNAscope® and BASESCOPE™ are trademarks of Advanced Cell Diagnostics, Inc. All other trademarks belong to their respective owners.

#### Headquarters

7707 Gateway Blvd Suite 200, Newark, CA 94545 Phone 1-510-576-8800 Toll Free 1-877-576-3636

For support, email [support@acdbio.com](mailto:support@acdbio.com)

[www.acdbio.com](http://www.acdbio.com)

