

BaseScope™ VS Assay

For Ventana DISCOVERY™ ULTRA System

RED

Document Number 323700-USM-ULT

For Research Use Only (RUO), Not for Diagnostic Use

Trademarks

BaseScope™ is a trademark of Advanced Cell Diagnostics, Inc. VENTANA and DISCOVERY are trademarks of Roche. All other trademarks belong to their respective owners.

Disclaimers

Advanced Cell Diagnostics, Inc. reserves the right to change its products and services at any time to incorporate technological developments. This manual is subject to change without notice.

Although this manual has been prepared with every precaution to ensure accuracy, Advanced Cell Diagnostics, Inc. assumes no liability for any errors, omissions, or for any damages resulting from the use of this information.

Copyright

© 2018. Advanced Cell Diagnostics, Inc. All rights reserved.

Contents

Chapter 1. Product Information.....	5
About this guide	5
Product description	5
Background.....	5
Overview	5
Kit contents and storage.....	6
BaseScope™ VS Probes.....	6
RNAscope® VS Control Slides.....	7
BaseScope™ VS Reagents	7
Required materials from Roche Diagnostics.....	8
Equipment and buffers	9
User-supplied materials.....	9
 Chapter 2. Before You Begin	 10
Important procedural guidelines	10
 Chapter 3. Prepare and Pretreat Samples	 11
Prepare FFPE sections	11
Materials required	11
Fix the sample	11
Dehydrate, embed, and cut the sample.....	11
 Chapter 4. Automated BaseScope™ VS Assay	 13
Workflow.....	13
Prepare the materials.....	14
Materials required	14
Prepare the instrument.....	14
Dilute bulk reagents	14
Register new reagents	14
Prepare instrument reagents	14
Create an instrument protocol	15
Print the labels.....	17
Run the BaseScope™ VS Assay	18
Materials required	18
Load the reagents	18
Start the run	18
Prepare detergent	19
Prepare dehydrating reagents.....	19
Complete the run	19
Wash the slides	19
Mount the samples	20
Recommended guidelines.....	20

Chapter 5. Evaluate the Results	21
Scoring guidelines	21
Troubleshooting.....	21
Staining example.....	22
Appendix A. Semi-automated BaseScope™ VS Assay	23
Workflow.....	23
Kit contents and storage.....	24
Prepare the materials.....	24
Materials required	24
Prepare the instrument.....	24
Dilute bulk reagents	24
Register new reagents	25
Prepare instrument reagents	25
Prepare deparaffinization reagents	25
Prepare 1X Target Retrieval	25
Create an instrument protocol	25
Print the labels.....	27
Manually pretreat the samples	27
Materials required	27
Bake the slides.....	28
Deparaffinize FFPE sections	28
Pretreat the slides.....	28
Run the BaseScope™ VS Assay.....	30
Materials required	30
Load the reagents	30
Start the run	30
Prepare detergent	31
Prepare dehydrating reagents.....	31
Complete the run	31
Wash the slides	31
Mount the samples.....	31
Recommended guidelines.....	32
Appendix B. Safety	33
Chemical safety	33
In the U.S.:	33
In the EU:	34
Documentation and Support	35
Obtaining SDSs.....	35
Obtaining support	35
Contact information	35
Limited product warranty.....	35

1

Chapter 1. Product Information



Before using this product, read and understand the information in **Appendix B. Safety** on page 33 in this document.

IMPORTANT! We recommend reading the entire user manual before beginning any protocols.

About this guide

This user manual provides two versions of the BaseScope™ VS Assay:

- **Chapter 4. Automated BaseScope™ VS Assay** starting on page 13.
- **Appendix A. Semi-automated BaseScope™ VS Assay** starting on page 23.

Product description


Background

The BaseScope™ Assay uses a novel and proprietary method of *in situ* hybridization (ISH) to visualize single RNA molecules for splice variants and short targets in samples mounted on slides. BaseScope™ Assays do not require the RNA-free environment used for traditional ISH. The assays are based on ACD's patented signal amplification and background suppression technology. Compared with the RNAscope® 2.5 Assay, the BaseScope™ Assay incorporates an additional signal amplification step, which makes it possible to detect RNA splicing variants, point mutations, small insertions or deletions, and short RNA targets (50–300 nucleotides).

Overview

Figure 1 on page 6 illustrates the BaseScope™ VS Assay procedure. You can complete the procedure on the Ventana DISCOVERY™ ULTRA System in ~12–13 hours. Start with properly prepared samples, pretreat them, and then hybridize RNA-specific probes to target RNA. Amplify the signal using multiple steps, followed by hybridization to alkaline phosphatase (AP)-labeled probes and detection using a red chromogenic substrate. Each single RNA transcript appears as a distinct dot of chromogen precipitate visible when using a common bright-field microscope at 40X magnification.

Figure 1. Procedure overview



Kit contents and storage

The BaseScope™ VS Assay requires the BaseScope™ VS Probes and the BaseScope™ VS Detection Reagent Kit. Probes and Reagent Kits are available separately.

IMPORTANT! BaseScope™ VS Probes must be used with the BaseScope™ VS Detection Reagent Kit. RNAscope® VS probes are incompatible with the BaseScope™ Detection Reagent Kit.

BaseScope™ VS Probes

The BaseScope™ VS Probes consist of the user-specified Target Probe and the Positive and Negative Control Probes. Visit <https://acdbio.com/products> to find a gene-specific Target Probe or appropriate Control Probes. Each probe is sufficient for staining ~30 standard slides. The probes have a shelf life of two years from the manufacturing date when stored as indicated in the following tables:

Target Probes					
<input checked="" type="checkbox"/>	Reagent	Cat. No.	Content	Quantity	Storage
	BaseScope™ VS Target Probe ([species]– [gene])	Various	Probe targeting specific RNA	7 mL x 1 bottle	2–8°C
Control Probes					
<input checked="" type="checkbox"/>	Reagent	Cat. No.	Content	Quantity	Storage
	BaseScope™ VS Positive Control Probe-Human (Hs)-PPIB-3ZZ	701039	Probe targeting common housekeeping gene	7 mL x 1 bottle	2–8°C
	BaseScope™ VS Positive Control Probe-Mouse (Mm)-Ppib-3ZZ	701079	Probe targeting common housekeeping gene	7 mL x 1 bottle	2–8°C
	BaseScope™ VS Positive Control Probe-Human (Hs)-PPIB-1ZZ	701049	Probe targeting common housekeeping gene	7 mL x 1 bottle	2–8°C
	BaseScope™ VS Positive Control Probe-Mouse (Mm)- Ppib-1ZZ	701089	Probe targeting common housekeeping gene	7 mL x 1 bottle	2–8°C
	BaseScope™ VS Negative Control Probe-DapB-3ZZ	701019	Probe targeting bacterial gene dapB	7 mL x 1 bottle	2–8°C
	BaseScope™ VS Negative Control Probe-DapB-1ZZ	701029	Probe targeting bacterial gene dapB	7 mL x 1 bottle	2–8°C

IMPORTANT! When running the BaseScope™ VS assay, make sure that your control probes contain the same number of ZZ pairs as your target probe. Consult support at support.acd@bio-technne.com.

RNAscope® VS Control Slides

The RNAscope® VS Control Slides (Cat. No. 310045 for Human control slide, HeLa; Cat. No. 310023 for Mouse control slide, 3T3) contain FFPE cell pellets sectioned and mounted on slides. The control slides can be used for assay control with the BaseScope™ VS Positive Control Probes and the BaseScope™ VS Negative Control Probes. The slides have a shelf life of nine months from the manufacturing date when stored at 2–8°C with desiccants.

BaseScope™ VS Reagents

BaseScope™ VS Reagent kits provide enough reagents to stain ~60 standard slides. You will receive two kits when you order the BaseScope™ VS Reagent Kit (Cat. No. 323700). BaseScope™ VS Reagents include:

- BaseScope™ VS Detection Reagents (Cat. No. 323710)
- RNAscope® VS Universal Sample Prep Reagent Kit (Cat. No. 323220)
- RNAscope® VS Accessory Kit (Cat. No. 320630)

The reagents are Ready-To-Use (RTU), and have a shelf life of nine months from the manufacturing date when stored as indicated in the following tables:

BaseScope™ VS Detection Reagents (Cat. No. 323710)				
<input checked="" type="checkbox"/>	Reagent	Cat. No.	Quantity	Storage
	BaseScope™ VS AMP 1	323711	14 mL x 1 bottle	2–8°C
	BaseScope™ VS AMP 2	323712	14 mL x 1 bottle	2–8°C
	BaseScope™ VS AMP 3	323713	14 mL x 1 bottle	2–8°C
	BaseScope™ VS AMP 4	323714	14 mL x 1 bottle	2–8°C
	BaseScope™ VS AMP 5	323715	14 mL x 1 bottle	2–8°C
	BaseScope™ VS AMP 6	323716	14 mL x 1 bottle	2–8°C
	BaseScope™ VS AMP 7	323717	14 mL x 1 bottle	2–8°C
	BaseScope™ VS AMP 8	323718	14 mL x 1 bottle	2–8°C
	RNAscope® VS Protease	322218	14 mL x 1 bottle	2–8°C
RNAscope® VS Universal Sample Prep Reagents (Cat. No. 323741)				
<input checked="" type="checkbox"/>	Reagent	Cat. No.	Quantity	Storage
	RNAscope® VS Universal Target Retrieval v2	323741	14 mL x 2 bottles	Room Temp (15–30°C)
	RNAscope® VS Universal Dewax	323742	14 mL x 1 bottle	Room Temp (15–30°C)
RNAscope® VS Accessory Kit (Cat. No. 320630)				
<input checked="" type="checkbox"/>	Reagent	Cat. No.	Quantity	Storage
	RNAscope® VS Hematoxylin	320631	7 mL x 1 bottle	2–8°C
	RNAscope® VS Bluing Reagent	320632	7 mL x 1 bottle	2–8°C

IMPORTANT! Dewax must be in solution and at room temperature before use on the instrument. If stored cold, place at **37°C** for **15 MIN** before each use regardless of the prior storage condition, since it may precipitate during shipment.

IMPORTANT! BaseScope™ VS and RNAscope® VS Universal assays share some of the same reagents including RNAscope® VS Protease, RNAscope® VS Target Retrieval v2, RNAscope® VS Universal Dewax, RNAscope® VS Hematoxylin, and RNAscope® VS Bluing Reagent. Only these reagents may be interchanged among kits. Do not interchange other reagents.

Required materials from Roche Diagnostics

The BaseScope™ VS Assay requires specific materials and equipment available *only* from Roche Diagnostics (Ventana Medical Systems, Inc.). Catalog Numbers are valid in the U.S. only. For other regions, please check Catalog or ordering numbers with your local lab supplier.

Probe Dispensers (Cat. No. 960-761 to 960-785; for Ordering Code, please contact local Roche representative)		
<input checked="" type="checkbox"/>	Component	Storage
	250 Test Probe #1–20 dispensers — fill dispensers with BaseScope™ VS Probes. Use up to 20 probes at one time.	Room Temp (15–30°C)
mRNA Sample Prep Kit (Cat. No. 760-248; Ordering Code 08127166001)		
<input checked="" type="checkbox"/>	Component	Storage
	mRNA Target Retrieval dispenser-fill dispenser with RNAscope® VS Universal Target Retrieval v2	Room Temp (15–30°C)
	mRNA Dewax dispenser — fill dispenser with RNAscope® VS Universal Dewax	Room Temp (15–30°C)
	mRNA Protease dispenser — fill dispenser with RNAscope® VS Protease	Room Temp (15–30°C)
mRNA RED Probe Amplification Kit (Cat. No. 760-236; Ordering Code 7095341001)		
<input checked="" type="checkbox"/>	Component	Storage
	ACD RED AMP 1 dispenser — fill dispenser with BaseScope™ VS AMP 1	Room Temp (15–30°C)
	ACD RED AMP 2 dispenser — fill dispenser with BaseScope™ VS AMP 2	Room Temp (15–30°C)
	ACD RED AMP 3 dispenser — fill dispenser with BaseScope™ VS AMP 3	Room Temp (15–30°C)
	ACD RED AMP 4 dispenser — fill dispenser with BaseScope™ VS AMP 4	Room Temp (15–30°C)
	ACD RED AMP 5 dispenser — fill dispenser with BaseScope™ VS AMP 5	Room Temp (15–30°C)
	ACD RED AMP 6 dispenser — fill dispenser with BaseScope™ VS AMP 6	Room Temp (15–30°C)
	ACD RED AMP 7 dispenser — fill dispenser with BaseScope™ VS AMP 7	Room Temp (15–30°C)
Ancillary Dispensers (Cat. No.771-758, Ordering Code 05271916001)		
<input checked="" type="checkbox"/>	Component	Storage
	250 Test Option #8 — fill dispenser with BaseScope™ VS AMP 8	Room Temp (15–30°C)
mRNA RED Detection Kit (Cat. No. 760-234; Ordering Code 7099037001)		
<input checked="" type="checkbox"/>	Component	Storage
	mRNA Inhibitor–prefilled	2–8°C
	mRNA Activator dispenser–prefilled	2–8°C
	mRNA Naphthol dispenser–prefilled	2–8°C
	mRNA Fast Red dispenser–prefilled	2–8°C

Generic Dispensers (Cat. No. 771-741; Ordering Code 05271720001, Cat. No. 771-742; Ordering Code 05271738001)

<input checked="" type="checkbox"/>	Component	Storage
	250 Test Counterstain 1 dispenser — fill dispenser with VS Hematoxylin	Room Temp (15–30°C)
	250 Test Counterstain 2 dispenser — fill dispenser with VS Bluing Reagent	Room Temp (15–30°C)

Equipment and buffers

<input checked="" type="checkbox"/>	Component	Cat. No.	Ordering Code
	10X DISCOVERY Wash (RUO)	950-510	07311079001
	ULTRA LCS (Predilute)	650-210	05424534001
	SSC (10X)	950-110	05353947001
	Reaction Buffer (10X)	950-300	05353955001
	DISCOVERY CC1	950-500	06414575001

IMPORTANT! To run the BaseScope™ VS assay successfully, use DISCOVERY Wash (950-510). Do not use DISCOVERY EZ Prep. Place 2X SSC (950-110) in the SSC bulk container instead of Ribowash. You may fill the option bulk container with reaction buffer.

User-supplied materials

IMPORTANT! Do not substitute other materials for the SuperFrost® Plus Slides listed in the following table.

<input checked="" type="checkbox"/>	Description	Supplier	Cat. No.
	SuperFrost® Plus Slides (required)	Fisher Scientific	12-550-15
	100% ethanol (EtOH)	American Master Tech Scientific/MLS*	ALREAGAL
	Xylene	Fisher Scientific/MLS	X3P-1GAL
	10% neutral-buffered formalin (NBF)	MLS	—
	Paraffin wax	MLS	—
	1X PBS	MLS	—
	Microtome	MLS	—
	Drying oven, capable of holding temperature at 60 +/- 1°C	MLS	—
	EcoMount	Biocare	EM897L
	Tissue-Tek® Vertical 24 Slide Rack	American Master Tech Scientific/MLS	LWSRA24
	Tissue-Tek® Staining Dishes	American Master Tech Scientific/MLS	LWT4457EA
	Tissue-Tek® Clearing Agent Dishes, xylene resistant	American Master Tech Scientific/MLS	LWT4456EA
	Cover Glass 24 x 50 mm	Fisher Scientific/MLS	12--545-F
	Distilled water	MLS	—
	Dawn detergent or similar detergent	MLS	—
	Fume hood	MLS	—
	Optional: Glass beaker (1 or 2 L)	MLS	—
	Optional: Hot plate	Fisher Scientific/MLS	11-300-49SHP

* Major Laboratory Supplier in North America. For other regions, please check Catalog Numbers with your local lab supplier.

2

Chapter 2. Before You Begin

Prior to running the BaseScope™ VS Assay on your samples for the first time, we recommend that you:

- Be familiar with the Ventana™ DISCOVERY™ ULTRA system. Refer to the Ventana™ System User Manual.
- Run the assay on FFPE RNAscope® VS Control Slides (Cat. No. 310045 for Human control slide, HeLa; Catalog No. 310023 for Mouse control slide, 3T3) using the BaseScope™ Positive and Negative Control Probes.

Important procedural guidelines

- Start with properly fixed and prepared sections. Refer to **Chapter 3. Prepare and Pretreat Samples** on page 11, **Recommended guidelines** on page 20, and to our sample preparation and pretreatment user guides available at <https://acdbio.com/technical-support/user-manuals>.
- Regularly maintain and clean your automated staining instrument.
- Always run positive and negative control probes on your sample to assess sample RNA quality and optimal permeabilization.
- Do *not* substitute required materials. Assay has been validated with these materials only.
- Follow the protocol exactly for best results.
- Do *not* let your sections dry out during the procedure unless specified in the protocol.
- Use good laboratory practices and follow all necessary safety procedures. Refer to **Appendix B. Safety** on page 33 in this document for more information.

3

Chapter 3. Prepare and Pretreat Samples

Formalin-fixed, paraffin-embedded (FFPE) sample preparation is described in the following protocol.

IMPORTANT! We highly recommend following these guidelines. We cannot guarantee assay results with other preparation methods.

For samples treated differently from the following protocol, you may need to optimize pretreatment conditions. Refer to **Recommended guidelines** on page 20, and to <https://acdbio.com/technical-support/solutions>.

Prepare FFPE sections

Materials required

-
- 10% neutral buffered formalin (NBF)
 - 1X PBS
 - Paraffin wax
 - 100% ethanol (EtOH)
 - Xylene
 - Microtome
 - Water bath
 - SuperFrost® Plus slides
-

Fix the sample

1. Immediately following dissection, fix tissue in 10% NBF for **16–32 HRS** at **ROOM TEMPERATURE (RT)**. Fixation time will vary depending on tissue type and size.

 **CAUTION!** Handle biological specimens appropriately.

IMPORTANT! Fixation for **<16 HRS** or **>32 HRS** will impair the performance of the RNAscope® VS Universal Assay.

Dehydrate, embed, and cut the sample


IMPORTANT! Use fresh reagents.

1. Wash sample with 1X PBS.
2. Dehydrate sample using a standard ethanol series, followed by xylene.
3. Embed sample in paraffin using standard procedures.

Note: Embedded samples may be stored at 15–25°C with desiccation. To better preserve RNA quality over a long period (>1 yr), we recommend storing at 2–8°C with desiccation.

4. Trim paraffin blocks as needed, and cut embedded tissue into 5 +/- 1 µm sections using a microtome.

5. Place paraffin ribbon in a **40–45°C** water bath, and mount sections on **SUPERFROST® PLUS SLIDES**. Place tissue as shown for optimal staining:



IMPORTANT! Do not mount more than one section per slide. Place sections in the center of the slide.

6. Air dry slides **OVERNIGHT** at **RT**. Do NOT bake slides unless they will be used for BaseScope™ within 1 week.

OPTIONAL STOPPING POINT Store sections with desiccants, either at room temperature for up to three months, or at -20°C – 4°C for up to six months.


4

Chapter 4. Automated BaseScope™ VS Assay

IMPORTANT! We strongly recommend you run the RNAscope® VS Control Slides (Cat. No. 310045 or 310023) using the BaseScope™ VS Positive and Negative Control Probes along with your samples in every run.

Appendix A. Semi-automated BaseScope™ VS Assay on page 23 describes an offline boiling procedure for use with Cat. No.322000 (RNAscope Target Retrieval Reagents).

Workflow



Prepare the materials

Materials required

Materials Provided by Advanced Cell Diagnostics	Materials Provided by Ventana™ Medical Systems	Other Materials and Equipment
<ul style="list-style-type: none"> • BaseScope™ VS Target Probe • BaseScope™ VS Positive Control Probe • BaseScope™ VS Negative Control Probe • RNAscope® VS Universal Dewax • RNAscope® VS Target Retrieval • RNAscope® VS Protease • BaseScope™ VS AMP 1 • BaseScope™ VS AMP 2 • BaseScope™ VS AMP 3 • BaseScope™ VS AMP 4 • BaseScope™ VS AMP 5 • BaseScope™ VS AMP 6 • BaseScope™ VS AMP 7 • BaseScope™ VS AMP 8 • RNAscope® VS Hematoxylin • RNAscope® VS Bluing Reagent 	<ul style="list-style-type: none"> • DISCOVERY™ ULTRA — automated slide stainer • DISCOVERY Wash Buffer 10X • ULTRA LCS (Predilute) • SSC Buffer 10X • DISCOVERY CC1 • Reaction Buffer (10X) • Probe dispensers • mRNA Sample Prep Kit • mRNA Red Probe Amplification Kit • mRNA Red Detection Kit • User fillable dispensers • Option 8 dispenser 	<ul style="list-style-type: none"> • Distilled water • Dawn detergent or similar detergent • Fume hood • Xylene • Tissue-Tek® Staining Dish • Tissue-Tek® Clearing Agent Dish, xylene-resistant • Tissue-Tek® Vertical 24 Slide Rack • EcoMount • Cover Glass, 24 mm x 50 mm

Prepare the instrument

If the instrument has not been used for ≥ 1 week, follow the guidelines for instrument maintenance in the *Ventana™ System User Manual*.

Dilute bulk reagents

Please prepare the bulk fluids as per manufacturer's instructions.

Register new reagents

Reagent dispensers come with appropriate barcode labels and registration buttons for dispensing BaseScope® VS Reagents. Refer to the *Ventana™ DISCOVERY ULTRA System User Manual* for details.

To register reagents:

1. Log all ACD reagents and probes into the software as "log user-fillable reagents" or "log user-fillable probes".
2. Use the wand that comes with the instrument to register *new* reagent kits.

Prepare instrument reagents

Refer to the tables on pages 8–9 to determine the proper dispenser for each reagent.

1. For BaseScope™ VS reagents AMP1 to AMP 7, transfer the entire volume of each component into the correspondingly labeled mRNA Red Probe Amplification Kit dispenser.
2. For BaseScope™ VS AMP 8, transfer the entire volume into the user-fillable Option 8 dispenser.

- Transfer the BaseScope™ VS Target Probe, BaseScope™ VS Positive Control Probe, BaseScope™ VS Negative Control Probe, RNAscope® VS Universal Dewax, RNAscope® VS Protease, both bottles of RNAscope® VS Target Retrieval v2, RNAscope® VS Hematoxylin, and RNAscope® VS Bluing Reagent to the correspondingly labeled dispensers.

IMPORTANT! Avoid cross contamination between reagents. Dewax must be warmed to room temperature and be completely in solution before use.


- Follow the dispenser product insert instructions to properly prime and handle the dispensers.
- Store tightly-capped dispensers (except the mRNA Dewax dispenser) at **4°C** when not in use. Store the tightly-capped mRNA Dewax dispenser at **15–30°C**.

IMPORTANT! Do not use expired reagents.


- Empty the waste bottle if needed.

Create an instrument protocol

- Open the VS software and click on the **Protocol** button.
- Click on **Create/Edit Protocols**, go to the Procedure drop down menu and select **mRNA Universal**.
- Main protocol steps appear as shown:



4. After selecting the main protocol steps, drop down menus become available. Select the appropriate protocol steps by clicking on the associated check boxes as shown:




The screenshot shows the 'Protocols' software interface. The 'Procedure' dropdown is set to 'mRNA Universal'. The 'Protocol Quick Find' table lists the following items:

Number	Name
1461	P1 UnivBA 16-16 A7=12m
1462	P2 UnivBA 16-16 A7=12m
1612	P1 UnivBA 24-16 A7=12m
1613	P2 UnivBA 24-16 A7=12m

The 'Protocol Steps for procedure mRNA Universal' are listed below:

- List Only Registered Products
- mRNA Universal Procedure v2]
- [All staining done using this procedure is for Research Use Only]
- Wet Slide Load
- Delay
- [Delay refers to a time delayed start: Select time until run start]
- Baking
 - [RECOMMENDED: Set time to 32 minutes]
 - Warmup Slide to 60 Deg C, and Incubate for [0 Hr 32 Min] (Baking)
- Deparaffinization
 - [Dewax]
- Cell Conditioning
 - [RECOMMENDED: Set to 97 C and 16min for FFPE cell pellets or 24min for normal FFPE tissue]
 - 8 Minutes
 - [Target: Retrieval]
 - Warmup Slide to [97 Deg C] from All Temperatures (Cycle 1)
 - 16 Minutes
 - 24 minutes
 - 32 minutes
- [Only select one mRNA option as multiple selections will yield negative results]
- BaseScope
 - mRNA Duplex
 - mRNA AP Detection
 - mRNA HRP Detection
 - Dual Sequence
 - Counterstain
 - Slide Cleaning



The screenshot shows the 'Protocols' software interface with the 'Procedure' dropdown set to 'mRNA Universal'. The 'Protocol Quick Find' table is the same as in the previous screenshot.

The 'Protocol Steps for 1461 P1 UnivBA 16-16 A7=12m' are listed below:

- BaseScope
 - [BaseScope Utilizes Standard Pretreatment]
 - BaseScope 3rd Pretreatment
 - [RECOMMENDED: Set to 37 C and 16min for normal FFPE samples]
 - Warmup Slide to [37 Deg C], and Incubate for [16 Minutes] (Pretreatment #3 Temp RB)
 - [Specify Probe, hybridization time and Temperature]
 - Apply Two Drops of [PROBE 1] (Probe #1), Apply Coverslip, and Incubate for 4 Minutes
 - BaseScope Probe - Temp = 50 C Recommended
 - Warmup Slide to [50 Deg C], and Incubate for 2 Hours (Hybridization)
 - BaseScope Amp 1 - Temp = 53 C Recommended
 - Warmup Slide to [49 Deg C] and Incubate for 32 Minutes (Detection #1 Temp)
 - BaseScope Amp 2 - Temp = 53 C Recommended
 - Warmup Slide to [49 Deg C] and Incubate for 32 Minutes (Detection #2 Temp)
 - BaseScope Amp 4 - Temp = 53 C Recommended
 - Warmup Slide to [49 Deg C] and Incubate for 32 Minutes (Detection #3 Temp)
 - BaseScope Amp 5 - Temp = 53 C Recommended
 - Warmup Slide to [49 Deg C] and Incubate for 32 Minutes (Detection #4 Temp)
 - BaseScope Amp 7 - Recommended 4min
 - Incubate for [0 Hr 12 Min]
- mRNA Duplex
- mRNA AP Detection
- mRNA HRP Detection
- Dual Sequence
- Counterstain
 - Apply One Drop of [COUNTERSTAIN 1] (Counterstain), Apply Coverslip, and Incubate for [8 Minutes]
 - Post Counterstain
 - Apply One Drop of [COUNTERSTAIN 2] (Post Counterstain), Apply Coverslip, and Incubate for [4 Minutes]
- HE Stain

- Select the appropriate assay conditions from the drop down menus according to the following tables:

Tissue Type	Cell Conditioning Temperature	Time
Brain	97°C	16 MIN
Cell pellet	97°C	16 MIN
Colon	97°C	24 MIN
Heart	97°C	24 MIN
Intestine	97°C	24 MIN
Kidney	97°C	24 MIN
Liver	97°C	24 MIN
Lung	97°C	24 MIN
Prostate	97°C	24 MIN
Spleen	97°C	24 MIN
Tonsil	97°C	24 MIN

Standard Temperatures/Times	
VS Protease	37°C, 16 MIN
Standard probe temperatures	50°C
Standard BaseScope™ AMP 1 temperature	49°C !
Standard BaseScope™ AMP 2 temperature	49°C !
Standard BaseScope™ AMP 4 temperature	49°C !
Standard BaseScope™ AMP 5 temperature	49°C †
Standard BaseScope™ AMP 7 incubation time*	4 MIN

*Calibrate the staining intensity for your instrument using the BaseScope™ AMP 7 incubation time. We suggest titrating signal intensity using incubation times of 4 MIN, 12 MIN, and 24 MIN. If the instrument settings have been previously optimized for the mRNA Universal software, you can use the same time setting that you used for RNAscope® AMP 5.

! Please note that these are the latest recommendations from ACD and may differ from the recommendations prompted by the VS software.

† If a hazy nuclear background is observed please increase the Amp 5 incubation temperature to 53°C.

- Click **Save as**, then select a protocol number from the drop down menu and choose a protocol name for each probe. Click **Save**.
- Click **Close** to go back to the main screen.
- Assign probe number from the list to each probe of interest. For each probe selected, assign a protocol.

Print the labels

- Select the **Print Label** icon from the upper right corner of the home page screen.

2. Select your preferred template or create a new template. To create a new template, refer to the *Ventana™ DISCOVERY ULTRA System User Manual* for details.
3. Click on **Protocol**.
4. Select the protocol you created for the BaseScope™ VS Assay and click on the **Add** button. When the protocols for all of the slides have been assigned, click on **Close/Print**.
5. Fill in the template for each slide. Select **Print** when completed.
6. Proceed to the following procedure **Load the reagents**.

Run the BaseScope™ VS Assay

Materials required

-
- Prepared slides
 - Prepared instrument reagents
 - Prepared detergent
 - Distilled water
 - Prepared dehydrating materials
 - Tissue-Tek® Vertical 24 Slide Rack
 - Tissue-Tek® Staining Dish
 - EcoMount
 - Cover Glass, 24 mm x 50 mm
 - Fume hood
 - Xylene
-

Load the reagents


1. Remove the nozzle caps from the filled dispensers and place each cap on the post located on the back of the dispenser.
2. If needed, remove any air bubbles at the nozzle tips by pushing down on the nozzle until the liquid reaches the tip of the nozzle or forms a small meniscus at the tip of the nozzle

IMPORTANT! Do not dispense any drops as this could compromise your drop inventory.

3. Load dispensers onto the reagent racks.
4. Remove the yellow locking ring from the dispensers in the prefilled mRNA Red Detection Kit. Refer to the instructions provided by Ventana™ Medical Systems.
5. Load the reagent racks onto the reagent carousel.

Start the run

1. Click the **Ready** button.




2. Eject slide drawers.
3. Load each slide onto a heater pad with the label facing upward and inward. Ensure that the slides sit securely on the pads.

IMPORTANT! Prior to loading the slides, ensure heater pads are completely dry. Wipe off any solution using laboratory tissue paper.

4. Close slide drawers.

5. Click the **Running** button. Automated assay will finish in **~12-13 HRS**.



IMPORTANT! Before leaving the instrument unattended, ensure all reagents and slides are successfully registered and the instrument is running.

Prepare detergent

1. Prepare 200 mL of diluted detergent by adding 1-2 drops detergent to 200 mL distilled water in a container with a cap.
2. Mix well by inverting the container four to five times.
3. Add diluted detergent to a Tissue-Tek® Staining Dish.

Note: Store diluted detergent at RT.

Prepare dehydrating reagents

- In a fume hood, fill a clearing agent dish with ~200 mL fresh xylene.

Note: Ensure all containers remain covered when not in use.

Complete the run

1. After the run is complete, remove the mRNA Dewax dispenser, place nozzle cap on the dispenser, and store at room temperature.
2. For the remaining reagents, place nozzle caps back on the dispensers and place racks onto magnet locking tray.

IMPORTANT! Store reagent racks at **4°C** until next use. Store the mRNA Dewax dispenser at room temperature.

Wash the slides

1. Submerge a Tissue-Tek® Slide Rack into the Tissue-Tek® Staining Dish containing 200 mL diluted detergent.
2. Open the instrument slide drawers and unload slides.
3. Decant solution on the slides into the slide drawer, then *immediately* load slides into the Tissue-Tek® Slide Rack submerged in detergent.
4. Rinse oil off the slides by moving the slide rack up and down in the dish 10 times.
5. Replace the detergent with distilled water and rinse slides by moving the slide rack up and down a minimum of 10 times.

- Repeat Step 5 three to five times.

Mount the samples

- Remove the slide rack from the staining dish and dry slides in a **60°C** dry oven for **30 MIN**.

IMPORTANT! The Red substrate is alcohol sensitive. Do not dehydrate the slides in alcohol.

- Cool the slides for **5 MIN** at **RT**.
- Briefly dip one slide into fresh pure xylene and *immediately* place 1–2 drops of EcoMount on the slide before the xylene dries.

IMPORTANT! Use the EcoMount mounting medium only.

- Carefully place 24 mm x 50 mm coverslip over the tissue section. Avoid trapping air bubbles.
- Repeat steps 2 and 3 for each slide.
- Air dry slides for at least **5 MIN**.
- Proceed to **Chapter 5. Evaluate the Results** on page 21.

Recommended guidelines

We highly recommend following the guidelines for Cell Conditioning (Target Retrieval), and Pretreatment #3 (Protease) conditions for:

- Any new or previously untested FFPE tissue types.
- Samples prepared differently than the protocol in **Chapter 3. Prepare and Pretreat Samples** on page 11.

IMPORTANT! VS pretreatment reagents and conditions are shared across BaseScope™ and RNAscope® assays. If pretreatment conditions have been previously established for a given tissue with the RNAscope® Universal assay, the same pretreatment conditions may be applied to the BaseScope™ VS assay.

- Stain four representative slides using the positive and negative control probes according to the following table:

Slide No.	Probe	Target Retrieval	Protease
1	Positive control	16 MIN	16 MIN
2	Negative control	16 MIN	16 MIN
3	Positive control	24 MIN	16 MIN
4	Negative control	24 MIN	16 MIN

- Evaluate staining and tissue morphology as in **Chapter 5. Evaluate the Results** and determine which pretreatment condition yielded the highest positive control signal and the lowest negative control signal. Using 1zz PPIB, the positive control signal should have a staining score of 1 or higher, and the negative control signal should be 0.
- Use the optimized pretreatment conditions to run the assay with the target probe.
- If none of the first pretreatment conditions are satisfactory, proceed to the following table for further optimization:

Slide No.	Probe	Target Retrieval	Protease
1	Positive control	24 MIN	32 MIN/50°C
2	Negative control	24 MIN	32 MIN/50°C
3	Positive control	48 MIN	16 MIN/37°C
4	Negative control	48 MIN	16 MIN/37°C
3	Positive control	48 MIN	32 MIN/50°C
4	Negative control	48 MIN	32 MIN/50°C

5. If no satisfactory pretreatment conditions can be found, please contact technical support at support.acd@bio-techne.com.

5

Chapter 5. Evaluate the Results

Examine tissue sections under a standard bright field microscope at 20–40X magnification:

- Assess tissue and cell morphology.
- Assess positive control signal strength. Positive control signal should be visible as punctate dots within the cell cytoplasm at 20–40X magnification.
- Assess negative control background. One dot to every 10 cells displaying background staining per 40X microscope field is acceptable.
- Evaluate target probe signal using the scoring guidelines in the next section.

Scoring guidelines

The BaseScope™ Assay enables a semi-quantitative scoring guideline utilizing the estimated number of punctate dots present within each cell boundary.

An example of how to develop such a guideline for semi-quantitative assessment of BaseScope™ staining intensity is presented below for a gene with expression level varying between 1 to > 10 copies per cell.

Note: If your gene expression level is higher or lower than this range, you may need to scale the criteria accordingly.

Categorize staining into five grades: **0**, **1+**, **2+**, **3+**, and **4+** according to the following table:

Staining Score	Microscope Objective Scoring*
0	No staining, or less than 1 dot/20 cells (40X magnification)
1	1 dot/cell (visible at 20–40X magnification)
2	2–3 dots/cell. No or very few dot clusters (visible at 20–40X magnification)
3	4–10 dots/cell. Less than 10% positive cells have dot clusters (visible at 20X magnification)
4	>10 dots/cell. More than 10% positive cells have dot clusters (visible at 20X magnification)

* Discount cells with artificially high nuclear background staining.


Troubleshooting

For troubleshooting information, please contact technical support at support.acd@bio-technne.com.


Staining example

If the assay is successful, the staining should look like the following images:

Figure 2. BaseScope™ VS Assay results in HeLa cells



1zz Hs-POLR2A (Positive Control)




1zz DapB (Negative Control)



Appendix A. Semi-automated BaseScope™ VS Assay

Most sample types can be fully automated on the Discovery ULTRA. Manual pretreatment may give a better result in some cases. Use the semi-automated procedure for tissues that do not have a satisfactory result when using the fully automated procedure. See **Chapter 4. Automated BaseScope™ VS Assay** on page 13.

Workflow



Kit contents and storage

For Offline Boiling: RNAscope® Target Retrieval Reagents				
<input checked="" type="checkbox"/>	Cat. No.	Reagent	Quantity	Storage
	322000	RNAscope® Target Retrieval Reagents*	70 mL x 4 bottles	Room Temp (15–30°C)

* Not provided with the kit and needs to be purchased separately.

IMPORTANT! Do not substitute the reagent components of the BaseScope™ VS Reagent Kit with those of other BaseScope™ or RNAscope® Reagent Kits, even those having the same name. The Target Retrieval solution in the BaseScope™ VS Reagent Kit CANNOT be used for offline boiling. Please separately purchase the RNAscope® Target Retrieval Reagents (Cat. No. 322000) to boil samples off the instrument.

Prepare the materials

Materials can be prepared ahead of time or while baking the slides, unless otherwise stated. See **Bake the slides** on page 28.

Materials required

Materials provided by Advanced Cell Diagnostics	Materials Provided by Ventana™ Medical Systems	Other Materials and Equipment
<ul style="list-style-type: none"> • RNAscope® 2.5 VS Target Probe • RNAscope® 2.5 VS Positive Control Probe • RNAscope® 2.5 VS Negative Control Probe • RNAscope® Target Retrieval Reagents • RNAscope® VS Protease • BaseScope™ VS AMP 1 • BaseScope™ VS AMP 2 • BaseScope™ VS AMP 3 • BaseScope™ VS AMP 4 • BaseScope™ VS AMP 5 • BaseScope™ VS AMP 6 • BaseScope™ VS AMP 7 • BaseScope™ VS AMP 8 • RNAscope® VS Hematoxylin • RNAscope® VS Bluing Reagent 	<ul style="list-style-type: none"> • DISCOVERY™ ULTRA — automated slide stainer • DISCOVERY Wash • ULTRA LCS (Predilute) • SSC Buffer 10X • Reaction Buffer • DISCOVERY CC1 • Probe dispensers • mRNA Sample Prep Kit • mRNA RED Probe Amplification Kit • mRNA Red Detection Kit • Option 8 dispenser • User fillable dispensers 	<ul style="list-style-type: none"> • Distilled water • Glass beaker (1 or 2 L) • Hot plate • Dawn detergent or similar detergent • Fume hood • Xylene • 100% ethanol (EtOH) • Tissue-Tek® Staining Dishes • Tissue-Tek® Clearing Agent Dishes, xylene-resistant • Tissue-Tek® Vertical 24 Slide Rack • EcoMount • Cover Glass, 24 mm x 50 mm

Prepare the instrument

If the instrument has not been used for ≥ 1 week, follow the guidelines for instrument maintenance in the *Ventana™ System User Manual*.

Dilute bulk reagents

Please prepare the bulk fluids as per manufacturer's instructions.

Register new reagents

Reagent dispensers come with appropriate barcode labels and registration buttons for dispensing BaseScope™ VS Reagents. Refer to the *Ventana™ DISCOVERY ULTRA System User Manual* for details. To register reagents:

1. Log all ACD reagents and probes into the software as “log user-fillable reagents” or “log user-fillable probes”.
2. Use the wand that comes with the instrument to register *new* reagent kits.

Prepare instrument reagents

Refer to the tables on pages 8–9 to determine the proper dispenser for each reagent.

1. For BaseScope™ VS AMP 1 to AMP 7, transfer the entire volume of each component into the correspondingly labeled mRNA Red Probe Amplification kit dispenser.
2. For BaseScope™ VS AMP 8, transfer the entire volume into the user-fillable Option 8 dispenser.
3. Transfer the BaseScope™ VS Target Probe, BaseScope™ VS Positive Control Probe, BaseScope™ VS Negative Control Probe, VS Protease, VS Hematoxylin, and VS Bluing Reagent to the correspondingly labeled dispensers.
4. Follow the dispenser product insert instructions to properly prime and handle the dispensers.
5. Store tightly-capped dispensers at **4°C** when not in use.
6. Check solution levels: DISCOVERY Wash, 2X SSC, Reaction Buffer, ULTRA LCS (Predilute), and CC1. Refill if they are less than half full. To run a full tray of 30 slides, fully fill the 2X SSC and Reaction Buffer containers.

IMPORTANT! Do not use expired reagents.

7. Empty the waste carboy if needed.

Prepare deparaffinization reagents

- In a fume hood, fill two clearing agent dishes with ~200 mL fresh xylene.
- In a fume hood, fill two staining dishes with ~200 mL fresh 100% EtOH.

Note: Ensure all containers remain covered when not in use.

Prepare 1X Target Retrieval



Prepare 1X Target Retrieval while FFPE slides are baking at 60°C, or the following day if you choose the optional stopping point on page 28. 1X Target Retrieval is used in manual cell conditioning (CC).

1. Prepare 700 mL of fresh 1X Target Retrieval by adding 630 mL distilled water to 1 bottle (70 mL) 10X 1X Target Retrieval solution in the beaker.
2. Mix well and cover the beaker with foil.

IMPORTANT! Do not use RNAscope® VS Universal Target Retrieval v2 for offline boiling.

Create an instrument protocol

1. Open the VS software and click on the **Protocol** button.
2. Click on **Create/Edit Protocols**, go to the Procedure drop down menu, and select **mRNA Universal**.
3. Main protocol steps appear as shown:



IMPORTANT! Do not select Baking, Deparaffinization, or Cell Conditioning.

- After selecting the main protocol steps, drop down menus become available. Select the appropriate protocol steps by clicking on the associated check boxes as shown above.

- Select the appropriate assay conditions from the drop down menus according to the following tables:

Standard Temperatures/Times	
VS Protease	37°C, 16 MIN
Standard probe temperatures	50°C
Standard BaseScope™ AMP 1	49°C !
Standard BaseScope™ AMP 2	49°C !
Standard BaseScope™ AMP 4	49°C !
Standard BaseScope™ AMP 5	49°C !‡
Standard BaseScope™ AMP 7 incubation time*	4 MIN

*Calibrate the staining intensity for your instrument using the BaseScope™ AMP 7 incubation time. We suggest titrating signal intensity using incubation times of 4 MIN, 12 MIN, and 24 MIN. If the instrument settings have been previously optimized for the mRNA Universal software, you can use the same time setting that you used for RNAscope® AMP 5.

! Please note that these are the latest recommendations from ACD and may differ from the recommendations prompted by the VS software.

‡ If a hazy nuclear background is observed please increase the Amp 5 incubation temperature to 53°C.

- Click **Save As**, then select a protocol number from the drop-down menu and choose a protocol name for each probe. Click **Save**.
- Click **Close** to go back to the main screen.
- Assign a probe number from the list to each probe of interest. For each probe selected, assign a protocol.

Print the labels

- Select the **Print Label** icon from the bottom of the home page screen.
- Select your preferred template or create a new template. To create a new template, refer to the *Ventana™ System User Manual* for details.
- Select the protocol you created for the BaseScope™ VS Assay.
- Click on **Protocol** to add and print the label.

Manually pretreat the samples

Materials required

Materials Provided by the Target Retrieval Reagents (Cat. No. 322000)	Other Materials and Equipment
<ul style="list-style-type: none"> • RNAscope® Target Retrieval Reagents 	<ul style="list-style-type: none"> • Drying oven • FFPE slides • Tissue-Tek® Vertical 24 Slide Rack • Distilled water • Fume hood • Xylene • 100% ethanol (EtOH) • Tissue-Tek® Clearing Agent Dish • Tissue-Tek® Staining Dish • Glass beaker (1 or 2 L) • Hot plate

Bake the slides

1. Bake slides in a dry oven for **30-60 MIN** at **60°C**.

OPTIONAL STOPPING POINT Use immediately or store at **RT** with desiccants for ≤ 1 week. Prolonged storage may degrade sample RNA.

IMPORTANT! If you continue, prepare the materials for the following protocols while the slides are baking: **Deparaffinize FFPE sections**, **Pretreat the slides**, and **Run the BaseScope™ VS Assay**.

Deparaffinize FFPE sections

IMPORTANT! If you have not done so already, create a protocol for your instrument and print slide labels during this procedure. See pages 25–27.

1. Place slides in a Tissue-Tek® Slide Rack and submerge in the first xylene-containing clearing agent dish in the fume hood.
2. Incubate the slides in xylene for **5 MIN** at **RT**. Agitate the slides by occasionally lifting the slide rack up and down in the clearing agent dish.
3. Remove the slide rack from the first xylene-containing dish and *immediately* place in the second xylene-containing clearing agent dish in the fume hood.
4. Repeat Step 2.
5. Remove the slide rack from the second xylene-containing dish and *immediately* place in the staining dish containing 100% EtOH.
6. Incubate the slides in 100% EtOH for **1 MIN** at **RT** with agitation.
7. Repeat Step 6 with fresh 100% EtOH.
8. Remove the slides from the rack and place on absorbent paper with the section face-up. Air dry for **5 MIN** at **RT**.
9. While slides are drying, place printed labels on the slides.
10. Insert the slides into a Tissue-Tek® Slide Rack and proceed to condition the slides.

Pretreat the slides

Begin heating 1X Target Retrieval Buffer while FFPE slides are baking at 60°C or during the previous section.

IMPORTANT! Do not boil 1X Target Retrieval more than **30 MIN** before use.

1. Heat 1X Target Retrieval Buffer to **98–104°C**:
 - a. Place the beaker containing 1X Target Retrieval Buffer on the hot plate. Cover the beaker with foil and turn the hot plate on high for **10–15 MIN**.
 - b. Once 1X Target Retrieval Buffer reaches a slow boil (**98–104°C**), turn the hot plate to a lower setting to maintain the correct temperature. Check the temperature with a thermometer.
2. With a pair of forceps *very slowly* submerge the slide rack containing the slides into the boiling 1X Target Retrieval Buffer solution. Cover the beaker with foil and boil the slides for the amount of time specified in the following table:

Tissue Type	Target Retrieval Time
Brain and spinal cord	15 MIN
Breast cancer	15 MIN
Cell lines	10 MIN
Colon	15 MIN
GI tract	15 MIN
Head and neck cancer	15 MIN
Heart	15 MIN
Kidney	15 MIN
Liver	30 MIN
Lung	15 MIN
Lymphoma	10 MIN
Placenta	15 MIN
Prostate	15 MIN
Skin	15 MIN
Stomach	15 MIN
Thymus	10 MIN
Tonsil	10 MIN
Xenograft derived from cell lines	7 MIN
Xenograft derived from primary tumor	15 MIN

3. Immediately transfer the hot slide rack from the 1X Target Retrieval Buffer to a staining dish containing distilled water. Do not let the slides cool in Target Retrieval.
4. Wash slides 3–5 times by moving the Tissue-Tek® Slide Rack up and down in the distilled water.
5. Repeat Step 4 with fresh distilled water.
6. Remove the slides from the rack and place on absorbent paper with the section face-up. Air dry for **5 MIN** at **RT**.
7. While slides are drying, place printed labels on the slides.

IMPORTANT! Labels must be in place prior to the next section.

8. Proceed directly to **Load the reagents** on page 30.

Run the BaseScope™ VS Assay

Materials required

-
- Prepared slides
 - Prepared instrument reagents
 - Prepared detergent
 - Distilled water
 - Prepared dehydrating materials
 - Tissue-Tek® Vertical 24 Slide Rack
 - EcoMount
 - Cover Glass, 24 mm x 50 mm
-

Load the reagents


1. Remove the nozzle caps from the filled dispensers and place each cap on the post located on the back of the dispenser.
2. If needed, remove any air bubbles at the nozzle tips by pushing down on the nozzle until the liquid reaches the tip of the nozzle or forms a small meniscus at the tip of the nozzle

IMPORTANT! Do not dispense any drops as this could compromise your drop inventory.

3. Load dispensers onto the reagent racks.
4. Remove the yellow locking ring from the dispensers in the prefilled mRNA RED Detection Kit. Refer to the instructions provided by Ventana™ Medical Systems.
5. Load the reagent racks onto the reagent carousel.

Start the run


1. Click the **Ready** button.



2. Eject slide drawers.
3. Load each slide onto a heater pad with the label facing upward and inward. Ensure that the slides sit securely on the pads.

IMPORTANT! Prior to loading the slides, ensure heater pads are completely dry. Wipe off any solution using laboratory tissue paper.

4. Close slide drawers.
5. Click the **Running** button. Semi-automated assay will finish in **~10 HRS.**



IMPORTANT! Before leaving the instrument unattended, ensure all reagents and slides are successfully registered and the instrument is running.

Prepare detergent

1. Prepare 200 mL of diluted detergent by adding 1 to 2 drops detergent to 200 mL distilled water in a container with a cap.
2. Mix well by inverting the container 4–5 times.
3. Add diluted detergent to a Tissue-Tek® Staining Dish.

Note: Store diluted detergent at RT.

Prepare dehydrating reagents

IMPORTANT! Do not use deparaffinization solutions for dehydration.

- In a fume hood, fill a clearing agent dish with ~200 mL fresh xylene.

Note: Ensure all containers remain covered when not in use.

Complete the run

1. After the run is complete, place nozzle caps back on the dispensers.
2. Store reagent racks at **4°C** until next use.

Wash the slides

1. Submerge a Tissue-Tek® Slide Rack into the Tissue-Tek® Staining Dish containing 200 mL diluted detergent.
2. Open the instrument slide drawer and unload slides.
3. Decant solution on the slides into the slide drawer, then *immediately* load slides into the Tissue-Tek® Slide Rack submerged in detergent.
4. Rinse oil off the slides by moving the slide rack up and down in the dish 10 times.
5. Replace the detergent with distilled water and rinse slides by moving the slide rack up and down a minimum of 10 times.
6. Repeat Step 5 three to five times.

Mount the samples

1. Remove the slide rack from the staining dish and dry slides in a **60°C** dry oven for **30 MIN**.

IMPORTANT! The RED substrate is alcohol sensitive. Do not dehydrate the slides in alcohol.

2. Cool the slides for **5 MIN** at **RT**.
3. Briefly dip one slide into fresh pure xylene and *immediately* place 1–2 drops of EcoMount on the slide before the xylene dries.

IMPORTANT! Use the EcoMount mounting medium only.

4. Carefully place 24 mm x 50 mm coverslip over the tissue section. Avoid trapping air bubbles.
5. Repeat steps 2 and 3 for each slide.
6. Air dry slides for at least **5 MIN**.

7. Proceed to **Chapter 5. Evaluate the Results**.

Recommended guidelines

We highly recommend following the guidelines for Cell Conditioning (Target Retrieval) and Protease conditions for:

- Any new or previously untested FFPE tissue types.
- Samples prepared differently than the protocol in **Chapter 3. Prepare and Pretreat Samples** on page 11.

1. Stain six representative slides using the positive and negative control probes according to the following table:

Slide No.	Probe	Target Retrieval	Protease
1	Positive control	10 MIN	24 MIN
2	Negative control	10 MIN	24 MIN
3	Positive control	10 MIN	16 MIN
4	Negative control	10 MIN	16 MIN
5	Positive control	15 MIN	16 MIN
6	Negative control	15 MIN	16 MIN

2. Evaluate staining and tissue morphology as in **Chapter 5. Evaluate the Results**, and determine which pretreatment condition yielded the highest positive control signal and lowest negative control signal. Using 1zz PPIB, positive control signal should have a staining score of 1 or higher, and the negative control signal should be 0.
3. Use the optimized pretreatment conditions to run the assay with the target probe.
4. If none of the conditions are satisfactory, contact technical support at support.acd@bio-techne.com.



Appendix B. Safety

Chemical safety



WARNING!

GENERAL CHEMICAL HANDLING. To minimize hazards, ensure laboratory personnel read and practice the general safety guidelines for chemical usage, storage, and waste provided below, and consult the relevant SDS for specific precautions and instructions:

- Read and understand the Safety Data Sheets (SDSs) before you store, handle, or work with any chemicals or hazardous materials. To obtain SDSs, see <https://acdbio.com/technical-support/user-manuals>.
- Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example, safety glasses, gloves, or protective clothing).
- Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with adequate ventilation (for example, fume hood).
- Characterize (by analysis if necessary) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
- Ensure that the waste is stored, transferred, transported, and disposed of according to all local, state/provincial, and/or national regulations.
- **IMPORTANT!** Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.

Biological hazard safety



WARNING!

BIOHAZARD. Biological samples such as tissues, body fluids, infectious agents, and blood of humans and other animals have the potential to transmit infectious diseases. Follow all applicable local, state/provincial, and/or national regulations. Wear appropriate protective equipment, which includes but is not limited to: protective eyewear, face shield, clothing/lab coat, and gloves. All work should be conducted in properly equipped facilities using the appropriate safety equipment (for example, physical containment devices). Individuals should be trained according to applicable regulatory and company/institution requirements before working with potentially infectious materials. Read and follow the applicable guidelines and/or regulatory requirements in the following:

In the U.S.:

- U.S. Department of Health and Human Services guidelines published in Biosafety in Microbiological and Biomedical Laboratories found at: <https://www.cdc.gov/biosafety/>
- Occupational Safety and Health Standards, Bloodborne Pathogens (29 CFR § 1910.1030), found at: https://www.osha.gov/pls/oshaweb/owadisp.show_document?p_id=10051&p_table=STANDARDS
- Your company's/institution's Biosafety Program protocols for working with/handling potentially infectious materials.

- Additional information about biohazard guidelines is available at:
<https://www.cdc.gov/biosafety/>

In the EU:

- Check local guidelines and legislation on biohazard and biosafety precaution and refer to the best practices published in the World Health Organization (WHO) Laboratory Biosafety Manual, third edition, found at:
http://www.who.int/csr/resources/publications/biosafety/WHO_CDS_CSR_LYO_2004_11/en/
- Information about the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH) can be found at: **<https://echa.europa.eu/regulations/reach>**

Documentation and Support

Obtaining SDSs

Safety Data Sheets (SDSs) are available at: <https://acdbio.com/technical-support/user-manuals>. For the SDSs of chemicals not distributed by Advanced Cell Diagnostics, contact the chemical manufacturer.

Obtaining support

For the latest services and support information, go to: <https://acdbio.com/technical-support/support-overview>.

At the website, you can:

- Access telephone and fax numbers to contact Technical Support and Sales facilities.
- Search through frequently asked questions (FAQs).
- Submit a question directly to Technical Support.
- Search for user documents, SDSs, application notes, citations, training videos, and other product support documents.
- Find out information about customer training events.

Contact information

Advanced Cell Diagnostics, Inc.
7707 Gateway Blvd Suite 200
Newark, CA 94560
Toll Free: 1-877-576-3636
Direct: 1-510-576-8800
Fax: 1-510-576-8801
Information: info.acd@bio-techne.com
Orders: orders.acd@bio-techne.com
Support Email: support.acd@bio-techne.com

Limited product warranty

Advanced Cell Diagnostics, Inc. and/or its affiliate(s) warrant their products as set forth in the ACD General Terms and Conditions of Sale found on the ACD website. If you have any questions, please contact Advanced Cell Diagnostics at <https://acdbio.com/about/contact>.

Headquarters

7707 Gateway Blvd Suite 200, Newark, CA 94560 Phone 1-510-576-8800 Toll Free 1-877-576-3636
For support, email support.acd@bio-technie.com
www.acdbio.com

