

# VS Universal RNAscope™ Assay combined with Immunohistochemistry: RNA-Protein Co-Detection on DISCOVERY ULTRA


## Introduction

This Technical Note provides guidelines for performing automated chromogenic co-detection of RNA and protein on the Ventana™ DISCOVERY ULTRA System. The Integrated Co-Detection Workflow (ICW) combines ACD's VS Universal RNAscope AP, VS Universal RNAscope HRP, VS RNAscope Duplex, or BaseScope™ VS detection assays with fully automated immunohistochemistry (IHC). RNA-Protein Co-Detection on the DISCOVERY ULTRA requires primary antibody, an RNAscope or BaseScope reagent kit, RNAscope or BaseScope target and control probes, the VS RNA-Protein Co-Detection Ancillary Kit, and additional materials from Roche Tissue Diagnostics. This workflow is compatible with NexES software versions 12.5.3 and greater from Roche Tissue Diagnostics, using the mRNA Universal procedure v6.03 or above. For every chemical, read the Safety Data Sheet (SDS) and follow handling instructions. Wear appropriate protective eyewear, clothing, and gloves. For the latest service and support information, go to [www.acdbio.com/support](http://www.acdbio.com/support).

## Getting Started

Before you begin RNA-Protein Co-Detection, we recommend establishing a working IHC protocol with your primary antibody, tissue, and chromogen of interest on the DISCOVERY ULTRA, and confirming the RNA quality of your tissue using your RNAscope/BaseScope assay of interest with corresponding control probes. You will also need to create staining protocols for RNA-Protein Co-Detection under the mRNA Universal procedure within the NexES software with the help of your ACD FAS. Ensure you have all necessary materials.

Protease treatment is required for RNA detection but can be disruptive to some protein targets. To determine the optimal workflow for RNA-Protein Co-Detection, we recommend performing a qualification run to determine whether the protein target is sensitive to protease, using Chapter 1: Protein Target Qualification. For protease-tolerant targets, proceed to Chapter 2: Sequential ISH-IHC. For protease-sensitive targets, proceed to Chapter 3: Integrated Co-Detection Workflow.



DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022


For Research Use Only. Not for use in diagnostics procedures.


## Integrated Co-Detection Workflow

## Sequential ISH-IHC Workflow

**General Workflow for RNA-Protein Co-Detection on DISCOVERY ULTRA:**

This **workflow** varies as indicated based on whether your protein target is protease-sensitive or protease-tolerant





## Chromogen Combinations for RNA-Protein Co-Detection on DISCOVERY ULTRA

For available options for ISH – IHC chromogen combinations, see the following table:

| ACD ISH Assay                                 | Reagents for ISH Detection   | IHC Detection System/Reagents  |
|---|--|--|
| VS Universal<br>RNAscope AP                   | RNAscope VS Universal AP Reagent Kit<br>and<br><b>mRNA Red Detection Kit</b>   | DISCOVERY AP-conjugated secondary<br>with <b>DISCOVERY Yellow (AP)</b> IHC chromogen<br>or<br>DISCOVERY HRP-conjugated secondary<br>with <b>DISCOVERY Teal HRP</b><br>or<br>DISCOVERY HRP-conjugated secondary<br>with <b>DISCOVERY Green HRP</b>  |
| VS Universal<br>RNAscope HRP<br>(chromogenic) | RNAscope VS Universal HRP Reagent Kit<br>and<br><b>mRNA Green Detection Kit</b><br>with DISCOVERY Inhibitor<br>or<br><b>mRNA Teal Detection Kit</b><br>with DISCOVERY Inhibitor<br>or<br><b>mRNA Purple Detection Kit</b><br>with DISCOVERY Inhibitor<br>or<br><b>mRNA DAB Detection Kit</b> | DISCOVERY AP-conjugated secondary<br>with <b>DISCOVERY Yellow (AP)</b> IHC Chromogen<br>or<br>DISCOVERY AP-conjugated secondary<br>with <b>ChromoMap Red (AP)</b><br>or<br>DISCOVERY HRP-conjugated secondary<br>with <b>DISCOVERY Teal HRP</b><br>or<br>DISCOVERY HRP-conjugated secondary<br>with <b>DISCOVERY Green HRP</b><br>or<br>DISCOVERY HRP-conjugated secondary<br>with <b>DISCOVERY Purple HRP</b> |
| VS Universal<br>RNAscope HRP<br>(fluorescent) | RNAscope VS Universal HRP Reagent Kit<br>and<br>DISCOVERY Fluorescent Detection Kit<br>with DISCOVERY Inhibitor  | DISCOVERY HRP-conjugated secondary<br>with DISCOVERY Fluorescent Detection Kit   |
| VS Universal<br>RNAscope Duplex               | RNAscope VS Duplex Reagent Kit<br>and<br><b>mRNA Red Detection Kit</b><br>and<br><b>mRNA Green Detection Kit</b><br>or<br><b>mRNA Teal Detection Kit</b><br>or<br><b>mRNA DAB Detection Kit</b>  | DISCOVERY AP-conjugated secondary<br>with <b>DISCOVERY Yellow (AP)</b> IHC Chromogen<br>or<br>DISCOVERY HRP-conjugated secondary<br>with <b>DISCOVERY Teal HRP</b><br>or<br>DISCOVERY HRP-conjugated secondary<br>with <b>DISCOVERY Green HRP</b>  |
| VS BaseScope (AP)                             | BaseScope VS Universal AP Reagent Kit<br>and<br><b>mRNA Red Detection Kit</b>  | DISCOVERY AP-conjugated secondary<br>with <b>DISCOVERY Yellow (AP)</b> IHC chromogen<br>or<br>DISCOVERY HRP-conjugated secondary<br>with <b>DISCOVERY Teal HRP</b><br>or<br>DISCOVERY HRP-conjugated secondary<br>with <b>DISCOVERY Green HRP</b>  |

## Materials for RNA-Protein Co-Detection on DISCOVERY ULTRA

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.

## ACD VS Chromogenic ISH Detection Kits

### RNAscope VS Reagent Kits

RNAscope 2.5 VS Probes and BaseScope VS Probes are available separately. Roche RNA-Protein Co-Detection can be performed with the RNAscope VS Universal AP Reagent Kit (Cat. No. 323250), RNAscope VS Universal HRP Reagent Kit (Cat. No. 323200), RNAscope VS Duplex Reagent Kit (Cat. No. 323300), and BaseScope VS Reagent Kit (Cat. No. 323700). Each RNAscope VS Reagent Kit includes the following Ready-To-Use (RTU) reagents to stain ~60 standard slides.

- RNAscope VS Universal Sample Prep Reagent Kit v2 (Cat. No. 323740)
- RNAscope VS Accessory Kit (Cat. No. 320630)
- Assay-Specific Detection Reagents

RNAscope VS Universal Sample Prep Reagent Kit v2 (Cat. No. 323740) and RNAscope VS Accessory Kit (Cat. No. 320630) are stored as indicated in the following tables:

| RNAscope VS Sample Prep Reagent Kit v2 (Cat. No. 323740) |   |                  |                     |
|--|---|------------------|---------------------|
| <input checked="" type="checkbox"/>                      | Reagent                                   | Quantity         | Storage             |
|  | RNAscope VS Universal Target Retrieval v2 | 10 mL x 2 bottle | Room Temp (15–30°C) |
|  | RNAscope VS Universal Dewax               | 14 mL x 1 bottle | Room Temp (15–30°C) |
| RNAscope VS Accessory Kit (Cat. No. 320630)              |   |                  |                     |
| <input checked="" type="checkbox"/>                      | Reagent                                   | Quantity         | Storage             |
|  | RNAscope VS Hematoxylin - RTU             | 7 mL x 1 bottle  | 2–8°C               |
|  | RNAscope VS Bluing Reagent - RTU          | 7 mL x 1 bottle  | 2–8°C               |

RTU Assay-specific Detection Reagents within each Reagent Kit are stored as indicated in the following tables:

| RNAscope VS Universal AP Detection Reagents (Cat. No. 323260) |                                |                  |         |
|---|--------------------------------|------------------|---------|
| <input checked="" type="checkbox"/>                           | Reagent                        | Quantity         | Storage |
|   | RNAscope VS Universal AP AMP 1 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 2 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 3 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 4 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 5 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 6 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 7 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Protease           | 14 mL x 1 bottle | 2–8°C   |

| RNAscope VS Universal HRP Detection Reagents (Cat. No. 323210) |                                 |                  |         |
|--|---------------------------------|------------------|---------|
| <input checked="" type="checkbox"/>                            | Reagent                         | Quantity         | Storage |
|  | RNAscope VS Universal HRP AMP 1 | 14 mL x 1 bottle | 2–8°C   |
|  | RNAscope VS Universal HRP AMP 2 | 14 mL x 1 bottle | 2–8°C   |
|  | RNAscope VS Universal HRP AMP 3 | 14 mL x 1 bottle | 2–8°C   |

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.

| RNAscope VS Universal HRP Detection Reagents (Cat. No. 323210) |                                 |                  |         |
|--|---------------------------------|------------------|---------|
| <input checked="" type="checkbox"/>                            | Reagent                         | Quantity         | Storage |
|  | RNAscope VS Universal HRP AMP 4 | 14 mL x 1 bottle | 2–8°C   |
|  | RNAscope VS Universal HRP AMP 5 | 14 mL x 1 bottle | 2–8°C   |
|  | RNAscope VS Universal HRP AMP 6 | 14 mL x 1 bottle | 2–8°C   |
|  | RNAscope VS Universal HRP AMP 7 | 14 mL x 1 bottle | 2–8°C   |
|  | RNAscope VS Protease            | 14 mL x 1 bottle | 2–8°C   |

| RNAscope VS Duplex Detection Reagents (Cat. No. 323310) |                             |                   |         |
|---|-----------------------------|-------------------|---------|
| <input checked="" type="checkbox"/>                     | Reagent                     | Quantity          | Storage |
|   | RNAscope VS Duplex AMP 1    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 2    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 3    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 4    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 5    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 6    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 7    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 8    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 9    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP Wash | 14 mL x 2 bottles | 2–8°C   |
|   | RNAscope VS Protease        | 14 mL x 1 bottle  | 2–8°C   |

| BaseScope VS Detection Reagents (Cat. No. 323710) |                                       |                  |         |
|---|---------------------------------------|------------------|---------|
| <input checked="" type="checkbox"/>               | Reagent                               | Quantity         | Storage |
|   | BaseScope VS AMP 1                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 2                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 3                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 4                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 5                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 6                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 7                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 8                    | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope 2.5 VS Pretreat 3 - Protease | 14 mL x 1 bottle | 2–8°C   |

## ACD Reagents for RNA-Protein Co-Detection on DISCOVERY ULTRA

The VS RNA-Protein Co-Detection Ancillary Kit includes materials for combined RNA and protein detection. Materials needed from this kit for ISH-IHC depend on the protease sensitivity of the protein target and if the staining is performed in an integrated or sequential manner.

**Co-Detection Antibody Diluent** helps preserve RNA prior to RNAscope detection and is recommended for use with the integrated co-detection workflow (ICW).

**VS Co-Detection Protease** helps counteract the impact of post-primary fixation on RNA accessibility and is recommended for use with the integrated co-detection workflow (ICW) on tissue samples.

**VS Co-Detection Inhibitor** is a chemical enzyme inhibitor recommended for use following the RNAscope/BaseScope assay prior to IHC detection. Use this reagent with both sequential and integrated co-detection workflows.

| VS RNA-Protein Co-Detection Ancillary Kit - Cat No. 323760 |  |                        |                         |         |
|--|--|------------------------|-------------------------|---------|
| <input checked="" type="checkbox"/>                        | Reagent  | Source / Ordering Info | Quantity                | Storage |
|  | Co-Detection Antibody Diluent                      | ACD/Cat No. 323160     | 120 mL x 1 bottle       | 2–8°C   |
|  | VS RNA-Protein Co-Detection Protease and Inhibitor | ACD/Cat No. 323190     | See the following table |         |

| VS RNA-Protein Co-Detection Protease and Inhibitor - Cat No. 323190 |                           |                  |         |
|---|---------------------------|------------------|---------|
| <input checked="" type="checkbox"/>                                 | Reagent                   | Quantity         | Storage |
|   | VS Co-Detection Protease  | 14 mL x 1 bottle | 2–8°C   |
|   | VS Co-Detection Inhibitor | 14 mL x 1 bottle | 2–8°C   |

## Required Materials from Roche Tissue Diagnostics

The Integrated Co-Detection Workflow (ICW) requires specific materials and equipment available *only* from Roche Tissue Diagnostics (Ventana Medical Systems, Inc.). Catalog numbers are valid in the United States only. For other regions, please check catalog or ordering numbers with your local lab supplier.

| Roche Materials Required for RNA-Protein Co-Detection |                                  |                    |                                    |   |  |
|---|----------------------------------|--------------------|------------------------------------|---|--|
| <input checked="" type="checkbox"/>                   | Component                        | Cat. No.           | Ordering Code                      | Fill with:  | Application:   |
|   | Probe Dispensers                 | 960-761 to 960-780 | Contact local Roche representative | RNAscope VS Probes  | Integrated Co-Detection Workflow<br>Sequential Co-Detection Workflow                                 |
|   | mRNA Sample Prep Kit             | 760-248            | 08127166001                        | RNAscope VS Sample Prep Reagent Kit v2 and VS Protease reagent from Detection Kit | Protein Target Qualification<br>Integrated Co-Detection Workflow<br>Sequential Co-Detection Workflow |
|   | Antibody Dispensers              | 770-001 to 770-099 | Contact local Roche representative | Primary Antibody Concentrate diluted in Co-Detection Antibody Diluent             | Protein Target Qualification<br>Integrated Co-Detection Workflow<br>Sequential Co-Detection Workflow |
|   | User-fillable Enzyme 1 Dispenser | 771-721            | 05271614001                        | VS Co-Detection Protease  | Integrated Co-Detection Workflow   |

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.

| Roche Materials Required for RNA-Protein Co-Detection |   |                           |               |                               |  |
|---|---|---------------------------|---------------|-------------------------------|--|
| <input checked="" type="checkbox"/>                   | Component                               | Cat. No.                  | Ordering Code | Fill with:                    | Application:   |
|   | User-fillable Fixative 1 Dispenser      | 771-731                   | 05271614001   | 10% Neutral Buffered Formalin | Integrated Co-Detection Workflow   |
|   | User-fillable Pretreatment 1 Dispenser  | 960-901                   | 5280095001    | VS Co-Detection Inhibitor     | Integrated Co-Detection Workflow<br>Sequential Co-Detection Workflow                                 |
|   | Counterstain 1 dispenser                | 771-741                   | 05271720001   | VS Hematoxylin                | Protein Target Qualification<br>Integrated Co-Detection Workflow<br>Sequential Co-Detection Workflow |
|   | Counterstain 2 dispenser                | 771-742                   | 05271738001   | VS Bluing Reagent             | Protein Target Qualification<br>Integrated Co-Detection Workflow<br>Sequential Co-Detection Workflow |
|   | DISCOVERY Antibody Block                | 760-4204                  | 05268869001   | Pre-filled                    | Protein Target Qualification<br>Integrated Co-Detection Workflow<br>Sequential Co-Detection Workflow |
|   | Additional ISH Assay-Specific Materials | (See the following table) |               |                               | Integrated Co-Detection Workflow<br>Sequential Co-Detection Workflow                                 |
|   | Additional IHC Detection Materials      | (See the following table) |               |                               | Protein Target Qualification<br>Integrated Co-Detection Workflow<br>Sequential Co-Detection Workflow |

| Additional Roche Materials Required for VS Universal AP ISH |                                  |          |               |   |
|---|----------------------------------|----------|---------------|---|
| <input checked="" type="checkbox"/>                         | Component                        | Cat. No. | Ordering Code | Fill with:  |
|   | mRNA RED Probe Amplification Kit | 760-236  | 07095341001   | RNAscope VS Universal AP Detection Reagents AMP 1-7 |
|   | mRNA RED Detection Kit           | 760-234  | 07099037001   | Pre-filled  |

| Additional Roche Materials Required for VS Universal HRP ISH |                                |          |               |  |
|--|--------------------------------|----------|---------------|--|
| <input checked="" type="checkbox"/>                          | Component                      | Cat. No. | Ordering Code | Fill with:   |
|  | mRNA Probe Amplification Kit   | 760-222  | 06614337001   | RNAscope VS Universal HRP Detection Reagents AMP 1-7 |
|  | DISCOVERY Inhibitor*           | 760-4840 | 07017944001   | Pre-filled   |
|  | mRNA Purple HRP Detection Kit† | 760-255  | 08127166001   | Pre-filled   |
|  | mRNA Green HRP Detection Kit†  | 760-278  | 08952612001   | Pre-filled   |
|  | mRNA Teal HRP Detection Kit†   | 760-256  | 08352941001   | Pre-filled   |
|  | mRNA DAB Detection Kit†        | 760-224  | 06614353001   | Pre-filled   |
|  | DISCOVERY DCC Kit†             | 760-240  | 07988192001   | Pre-filled   |
|  | DISCOVERY FAM Kit†             | 760-243  | 07988150001   | Pre-filled   |
|  | DISCOVERY FITC Kit†            | 760-232  | 07259212001   | Pre-filled   |
|  | DISCOVERY Rhodamine Kit†       | 760-233  | 07259883001   | Pre-filled   |
|  | DISCOVERY Rhodamine 6G Kit†    | 760-244  | 07988168001   | Pre-filled   |

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.



### Additional Roche Materials Required for VS Universal HRP ISH

|  |                        |         |             |            |
|--|------------------------|---------|-------------|------------|
|  | DISCOVERY Red 610 Kit† | 760-245 | 07988176001 | Pre-filled |
|  | DISCOVERY Cy5 Kit†     | 760-238 | 07551215001 | Pre-filled |

\* DISCOVERY Inhibitor is not required if using mRNA DAB Detection Kit

† Choose one mRNA Detection Kit for RNAscope HRP detection

### Additional Roche Materials Required for VS Duplex ISH

| <input checked="" type="checkbox"/> | Component                     | Cat. No. | Ordering Code | Fill with:  |
|-------------------------------------|-------------------------------|----------|---------------|---|
|                                     | mRNA Duplex Amp Kit           | 760-249  | 08127174001   | RNAscope VS Duplex Detection Reagents AMP 1-9, AMP Wash |
|                                     | mRNA Link (pre-filled)        | 760-6014 | 08127115001   | Pre-filled  |
|                                     | mRNA RED Detection Kit        | 760-234  | 07099037001   | Pre-filled  |
|                                     | mRNA Green HRP Detection Kit* | 760-278  | 08952612001   | Pre-filled  |
|                                     | mRNA Teal HRP Detection Kit*  | 760-256  | 08352941001   | Pre-filled  |
|                                     | mRNA DAB Detection Kit*       | 760-224  | 06614353001   | Pre-filled  |

\* Choose one mRNA Detection Kit for VS Duplex Channel 1 detection

### Additional Roche Materials Required for VS BaseScope ISH

| <input checked="" type="checkbox"/> | Component                        | Cat. No. | Ordering Code | Fill with:                              |
|-------------------------------------|----------------------------------|----------|---------------|---|
|                                     | mRNA RED Probe Amplification Kit | 760-236  | 07095341001   | BaseScope VS Detection Reagents AMP 1-7 |
|                                     | Option 8 dispenser               | 771-758  | 05271916001   | BaseScope VS AMP 8 reagent              |
|                                     | mRNA RED Detection Kit           | 760-234  | 07099037001   | Pre-filled                              |

### Reagent Options for Ventana IHC AP Detection

| <input checked="" type="checkbox"/> | Component                            | Cat. No. | Ordering Code | Storage |
|-------------------------------------|--------------------------------------|----------|---------------|---------|
|                                     | DISCOVERY UltraMap anti-Ms Alk Phos* | 760-4312 | 05269687001   | 2-8°C   |
|                                     | DISCOVERY UltraMap anti-Rb Alk Phos* | 760-4314 | 05269709001   | 2-8°C   |
|                                     | DISCOVERY Yellow Kit †               | 760-239  | 07698445001   | 2-8°C   |
|                                     | DISCOVERY ChromoMap Red Kit*         | 760-160  | 05266653001   | 2-8°C   |

| Reagent Options for Ventana IHC HRP Detection |                                  |          |               |         |
|---|----------------------------------|----------|---------------|---------|
| <input checked="" type="checkbox"/>           | Product                          | Cat. No. | Ordering Code | Storage |
|   | DISCOVERY UltraMap anti-Ms HRP*  | 760-4313 | 05269695001   | 2–8°C   |
|   | DISCOVERY UltraMap anti-Rb HRP*  | 760-4315 | 05269717001   | 2–8°C   |
|   | DISCOVERY UltraMap anti-Rat HRP* | 760-4456 | 05891884001   | 2–8°C   |
|   | DISCOVERY UltraMap anti-Gt HRP*  | 760-4648 | 06607241001   | 2–8°C   |
|   | DISCOVERY Purple Kit †           | 760-229  | 07053983001   | 2–8°C   |
|   | DISCOVERY Green HRP Kit †        | 760-278  | 07053983001   | 2–8°C   |
|   | DISCOVERY Teal HRP Kit †         | 760-247  | 08254338001   | 2–8°C   |
|   | DISCOVERY ChromoMap DAB Kit †    | 760-159  | 05266645001   | 2–8°C   |
|   | DISCOVERY DCC Kit †              | 760-240  | 07988192001   | 2–8°C   |
|   | DISCOVERY FAM Kit †              | 760-243  | 07988150001   | 2–8°C   |
|   | DISCOVERY FITC Kit†              | 760-232  | 07259212001   | 2–8°C   |
|   | DISCOVERY Rhodamine Kit †        | 760-233  | 07259883001   | 2–8°C   |
|   | DISCOVERY Rhodamine 6G Kit †     | 760-244  | 07988168001   | 2–8°C   |
|   | DISCOVERY Red 610 Kit†           | 760-245  | 07988176001   | 2–8°C   |
|   | DISCOVERY Cy5 Kit †              | 760-238  | 07551215001   | 2–8°C   |

\* Choose one secondary detection antibody depending on the primary antibody species and desired IHC chromogen

† Choose one IHC chromogen

### Instrument buffers

| <input checked="" type="checkbox"/> | Component                | Cat. No. | Ordering Code | Storage             |
|-------------------------------------|--------------------------|----------|---------------|---------------------|
|                                     | 10X DISCOVERY Wash (RUO) | 950-510  | 7311079001    | Room Temp (15–30°C) |
|                                     | ULTRA LCS (predilute)    | 650-210  | 5424534001    | Room Temp (15–30°C) |
|                                     | SSC Buffer (10X)         | 950-110  | 5353947001    | Room Temp (15–30°C) |
|                                     | Reaction Buffer (10X)    | 760-107  | 5266262001    | Room Temp (15–30°C) |
|                                     | DISCOVERY CC1            | 950-500  | 6414575001    | Room Temp (15–30°C) |

### Additional Required Materials

#### Materials for Co-Detection

| <input checked="" type="checkbox"/> | Reagent                      | Source / Ordering Info | Quantity  | Usage   |
|-------------------------------------|------------------------------|------------------------|-----------|---|
|                                     | Primary Antibody (RTU)       | User                   | As needed | <ul style="list-style-type: none"> <li>Protein Target Qualification (<i>optional</i>)</li> <li>Sequential ISH-IHC Workflow (<i>optional</i>)</li> <li>Integrated Co-Detection Workflow (<i>optional</i>)</li> </ul> |
|                                     | Primary Antibody Concentrate | User                   | As needed | <ul style="list-style-type: none"> <li>Protein Target Qualification (<i>optional</i>)</li> <li>Sequential ISH-IHC Workflow (<i>optional</i>)</li> </ul>   |

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.

| <input checked="" type="checkbox"/> | Reagent                       | Source / Ordering Info | Quantity | Usage  |
|-------------------------------------|-------------------------------|------------------------|----------|--|
|                                     |                               |                        |          | <ul style="list-style-type: none"> <li>Integrated Co-Detection Workflow (<i>optional</i>)</li> </ul> |
|                                     | 10% Neutral Buffered Formalin | User                   | 10 mL    | <ul style="list-style-type: none"> <li>Integrated Co-Detection Workflow</li> </ul>                   |

### Other materials

| <input checked="" type="checkbox"/> | Description   | Supplier                            | Cat. No.     |
|-------------------------------------|---|-------------------------------------|--------------|
|                                     |   |                                     |              |
|                                     | SuperFrost® Plus Slides (required)                            | Fisher Scientific                   | 12-550-15    |
|                                     | ProLong® Gold Antifade Reagent                                | Life Technologies                   | P36930       |
|                                     | EcoMount (if using Red detection)                             | Biocare                             | EM897L       |
|                                     | Tissue-Tek® Vertical 24 Slide Rack                            | American Master Tech Scientific/MLS | LWSRA24      |
|                                     | Tissue-Tek Staining Dishes                                    | American Master Tech Scientific/MLS | LWT4457EA    |
|                                     | Cover Glass 24 x 50 mm  | Fisher Scientific/MLS               | 12--545-F    |
|                                     | Distilled water   | MLS                                 | —            |
|                                     | Mild liquid dishwashing detergent (Dawn detergent or similar) | MLS                                 | —            |
|                                     | Drying oven, capable of holding temperature at 60 +/- 1°C     | MLS                                 | —            |
|                                     | Fume hood   | MLS                                 | —            |
|                                     | 100% ethanol (EtOH)   | MLS                                 | —            |
|                                     | Xylene  | MLS                                 | —            |
|                                     | Tissue-Tek Clearing Agent Dishes, xylene-resistant            | American Master Tech Scientific/MLS | LWT4456EA    |
|                                     | Optional: Glass beaker (1 or 2 L)                             | MLS                                 | —            |
|                                     | Optional: Hot plate   | Fisher Scientific/MLS               | 11-300-49SHP |

## Chapter 1: Protein Target Qualification

### Materials

The following tables provide general information on materials needed for protein target qualification. The full material list varies depending on your choice of primary antibody, IHC protein detection reagents, and IHC chromogen. Please refer to [Materials for RNA-Protein Co-Detection on DISCOVERY ULTRA](#).

#### ACD reagents for protein target qualification

Protein target qualification requires the following Ready-To-Use (RTU) reagents from ACD:

- RNAscope VS Protease (found within RNAscope or BaseScope Reagent Kit)
- RNAscope VS Universal Sample Prep Reagent Kit v2 (Cat. No. 323740)
- RNAscope VS Accessory Kit (Cat. No. 320630)

RNAscope VS Universal Sample Prep Reagent Kit v2 (Cat. No. 323740) and RNAscope VS Accessory Kit (Cat. No. 320630) are stored as indicated in the following tables:

| RNAscope VS Sample Prep Reagent Kit v2 (Cat. No. 323740) |   |                  |                     |
|--|---|------------------|---------------------|
| <input checked="" type="checkbox"/>                      | Reagent                                   | Quantity         | Storage             |
|  | RNAscope VS Universal Target Retrieval v2 | 10 mL x 2 bottle | Room Temp (15–30°C) |
|  | RNAscope VS Universal Dewax               | 14 mL x 1 bottle | Room Temp (15–30°C) |

RNAscope VS Protease reagent is a component of each ACD Detection Kit and is stored as indicated in the following tables:

| RNAscope Protease Component         |   |          |                                       |                  |         |
|-------------------------------------|---|----------|---------------------------------------|------------------|---------|
| <input checked="" type="checkbox"/> | ACD Detection Kit                               | Part No. | Reagent                               | Quantity         | Storage |
|                                     | RNAscope VS Universal AP Detection Reagents     | 323260   | RNAscope VS Protease                  | 14 mL x 1 bottle | 2–8°C   |
|                                     | RNAscope VS Universal HRP Detection Reagents    | 323610   | RNAscope VS Protease                  | 14 mL x 1 bottle | 2–8°C   |
|                                     | RNAscope VS Universal Duplex Detection Reagents | 323310   | RNAscope VS Protease                  | 14 mL x 1 bottle | 2–8°C   |
|                                     | BaseScope VS Detection Reagents                 | 323710   | RNAscope 2.5 VS Pretreat 3 - Protease | 14 mL x 1 bottle | 2–8°C   |

#### Optional ACD reagents for protein target qualification

The following additional reagents from ACD can be used for protein target qualification on Roche DISCOVERY ULTRA:

- Co-Detection Antibody Diluent from VS RNA-Protein Co-Detection Ancillary Kit (Cat No. 323760)
- RNAscope VS Accessory Kit (Cat. No. 320630)

| VS RNA-Protein Co-Detection Ancillary Kit - Cat No. 323760 |                               |                        |                   |         |
|--|-------------------------------|------------------------|-------------------|---------|
| <input checked="" type="checkbox"/>                        | Reagent                       | Source / Ordering Info | Quantity          | Storage |
|  | Co-Detection Antibody Diluent | ACD/Cat No. 323160     | 120 mL x 1 bottle | 2–8°C   |

| RNAscope VS Accessory Kit (Cat. No. 320630) |                            |                 |         |
|---|----------------------------|-----------------|---------|
| <input checked="" type="checkbox"/>         | Reagent                    | Quantity        | Storage |
|   | RNAscope VS Hematoxylin    | 7 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Bluing Reagent | 7 mL x 1 bottle | 2–8°C   |

### Roche materials for protein target qualification

Protein target qualification can be performed with either Roche RTU primary antibody or your choice of primary antibody concentrate diluted in the ACD Co-Detection Antibody Diluent. For a list of available Roche RTU primary antibodies and ordering information, please contact your local Roche representative.

The following additional materials from Roche are used for protein target qualification on DISCOVERY ULTRA. Catalog numbers are valid in the United States only. For other regions, please check catalog or ordering numbers with your local lab supplier.

### Reagent dispensers

| User-filled Dispensers Required for Protein Target Qualification |  |                    |                    |                                    |   |
|--|--|--------------------|--------------------|------------------------------------|---|
| <input checked="" type="checkbox"/>                              | Component  | Source             | Cat. No.           | Ordering Code                      | Fill with:  |
|  | mRNA Sample Prep Kit                                   | Roche              | 760-248            | 08127166001                        | RNAscope VS Sample Prep Reagent Kit v2 and VS Protease reagent from Detection Kit |
|  | User-fillable Antibody Dispensers ( <i>Optional*</i> ) | 770-001 to 770-099 | 770-001 to 770-099 | Contact local Roche representative | Primary antibody concentrate diluted in Co-Detection Antibody Diluent             |
|  | Counterstain 1 dispenser                               | 771-741            | 771-741            | 05271720001                        | VS Hematoxylin  |
|  | Counterstain 2 dispenser                               | 771-742            | 771-742            | 05271738001                        | VS Bluing Reagent   |

\*Protein target qualification can be performed with either Roche RTU primary antibody or your choice of primary antibody.

| Reagent Options for Ventana IHC AP Detection |                                      |          |               |         |
|--|--------------------------------------|----------|---------------|---------|
| <input checked="" type="checkbox"/>          | Component                            | Cat. No. | Ordering Code | Storage |
|  | DISCOVERY UltraMap anti-Ms Alk Phos* | 760-4312 | 05269687001   | 2–8°C   |
|  | DISCOVERY UltraMap anti-Rb Alk Phos* | 760-4314 | 05269709001   | 2–8°C   |
|  | DISCOVERY Yellow Kit*                | 760-239  | 07698445001   | 2–8°C   |
|  | DISCOVERY Red Kit                    | 760-228  | 07425333001   | 2–8°C   |
|  | DISCOVERY ChromoMap Red Kit*         | 760-160  | 05266653001   | 2–8°C   |

| Reagent Options for Ventana IHC HRP Detection |                                  |          |               |         |
|---|----------------------------------|----------|---------------|---------|
| <input checked="" type="checkbox"/>           | Product                          | Cat. No. | Ordering Code | Storage |
|   | DISCOVERY UltraMap anti-Ms HRP*  | 760-4313 | 05269695001   | 2–8°C   |
|   | DISCOVERY UltraMap anti-Rb HRP*  | 760-4315 | 05269717001   | 2–8°C   |
|   | DISCOVERY UltraMap anti-Rat HRP* | 760-4456 | 05891884001   | 2–8°C   |
|   | DISCOVERY UltraMap anti-Gt HRP*  | 760-4648 | 06607241001   | 2–8°C   |
|   | DISCOVERY Purple Kit †           | 760-229  | 07053983001   | 2–8°C   |
|   | DISCOVERY Green HRP Kit †        | 760-278  | 07053983001   | 2–8°C   |
|   | DISCOVERY Teal HRP Kit †         | 760-247  | 08254338001   | 2–8°C   |
|   | DISCOVERY ChromoMap DAB Kit †    | 760-159  | 05266645001   | 2–8°C   |
|   | DISCOVERY DCC Kit †              | 760-240  | 07988192001   | 2–8°C   |
|   | DISCOVERY FAM Kit †              | 760-243  | 07988150001   | 2–8°C   |
|   | DISCOVERY FITC Kit†              | 760-232  | 07259212001   | 2–8°C   |
|   | DISCOVERY Rhodamine Kit †        | 760-233  | 07259883001   | 2–8°C   |
|   | DISCOVERY Rhodamine 6G Kit †     | 760-244  | 07988168001   | 2–8°C   |
|   | DISCOVERY Red 610 Kit†           | 760-245  | 07988176001   | 2–8°C   |
|   | DISCOVERY Cy5 Kit †              | 760-238  | 07551215001   | 2–8°C   |

\* Choose one secondary detection antibody depending on the primary antibody species and desired IHC chromogen

† Choose one IHC chromogen

### Instrument buffers

| <input checked="" type="checkbox"/> | Component                | Cat. No. | Ordering Code | Storage             |
|-------------------------------------|--------------------------|----------|---------------|---------------------|
|                                     | 10X DISCOVERY Wash (RUO) | 950-510  | 7311079001    | Room Temp (15–30°C) |
|                                     | ULTRA LCS (Predilute)    | 650-210  | 5424534001    | Room Temp (15–30°C) |
|                                     | SSC Buffer (10X)         | 950-110  | 5353947001    | Room Temp (15–30°C) |
|                                     | Reaction Buffer (10X)    | 760-107  | 5266262001    | Room Temp (15–30°C) |
|                                     | DISCOVERY CC1            | 950-500  | 6414575001    | Room Temp (15–30°C) |

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.

## User-Sourced Materials

| <input checked="" type="checkbox"/> | Description   | Supplier                            | Cat. No.     |
|-------------------------------------|---|-------------------------------------|--------------|
|                                     | Primary Antibody (RTU)  | User                                | Various      |
|                                     | Primary Antibody Concentrate                                  | User                                | Various      |
|                                     | SuperFrost Plus Slides (required)                             | Fisher Scientific                   | 12-550-15    |
|                                     | ProLong Gold Antifade Reagent                                 | Life Technologies                   | P36930       |
|                                     | EcoMount (if using Red detection)                             | Biocare                             | EM897L       |
|                                     | Tissue-Tek Vertical 24 Slide Rack                             | American Master Tech Scientific/MLS | LWSRA24      |
|                                     | Tissue-Tek Staining Dishes                                    | American Master Tech Scientific/MLS | LWT4457EA    |
|                                     | Cover Glass 24 x 50 mm  | Fisher Scientific/MLS               | 12--545-F    |
|                                     | Distilled water   | MLS                                 | —            |
|                                     | Mild liquid dishwashing detergent (Dawn detergent or similar) | MLS                                 | —            |
|                                     | Drying oven, capable of holding temperature at 60 +/- 1°C     | MLS                                 | —            |
|                                     | Fume hood   | MLS                                 | —            |
|                                     | 100% ethanol (EtOH)   | MLS                                 | —            |
|                                     | Xylene  | MLS                                 | —            |
|                                     | Tissue-Tek Clearing Agent Dishes, xylene-resistant            | American Master Tech Scientific/MLS | LWT4456EA    |
|                                     | Optional: Glass beaker (1 or 2 L)                             | MLS                                 | —            |
|                                     | Optional: Hot plate   | Fisher Scientific/MLS               | 11-300-49SHP |

## Assay Procedure

### Prepare the DISCOVERY ULTRA

#### Prepare the instrument

If the instrument has not been used for > 1 week, follow guidelines for instrument maintenance from Roche Tissue Diagnostics. Before use, empty the waste carboys if needed.

#### Dilute instrument bulk reagents

1. Prepare the instrument bulk fluids according to the manufacturer's instructions.
2. Fill bulk solution containers for 1X DISCOVERY Wash, ULTRA LCS (predilute), and CC1 (predilute) to be at least half full. Fully fill bulk solution containers for 2X SSC and 1X Reaction Buffer.

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.

---

**IMPORTANT!** Do not use expired reagents.

---

### *Register new reagents*

Reagent dispensers come with appropriate barcode labels and registration buttons for dispensing RNAscope VS Universal Reagents. Refer to the *Ventana DISCOVERY ULTRA System User Manual* for details. To register reagents:

- Log all ACD reagents and probes into the software as **log user-fillable reagents** and **log user-fillable probes**, respectively.
- Use the reagent registration wand that comes with the instrument to register new reagent kits from Roche Tissue Diagnostics

### *Prepare user-fillable reagents for protein target qualification*

Refer to the table on page 13 to determine the proper dispenser for each user-fillable reagent.

---

**IMPORTANT!** Avoid cross contamination between reagents. Dewax must be warmed to room temperature and be completely in solution before use.

---


1. If working with primary antibody concentrate, prepare the primary antibody:
  - a. Dilute the primary antibody in Co-Detection Antibody Diluent using previously established conditions for IHC staining on the DISCOVERY ULTRA or antibody manufacturer's recommendations.
  - b. Transfer the diluted primary antibody to an Antibody dispenser.
2. Fill the mRNA Sample Prep Kit:
  - a. Transfer the contents of both bottles of VS Target Retrieval v2 from the RNAscope VS Sample Prep Reagent Kit v2 to the mRNA Target Retrieval Dispenser.
  - b. Transfer the VS Dewax reagent from the RNAscope VS Sample Prep Reagent Kit v2 to the mRNA Dewax Dispenser.

**Note:** Leave the mRNA Protease dispenser empty. This dispenser is not used but must be on-instrument for the run to start successfully.
3. Fill the user-fillable dispensers:
  - a. Transfer the VS Protease from the Detection Kit to the Enzyme 2 dispenser
  - b. If using VS Accessory reagents for counterstaining, transfer the VS Hematoxylin and VS Bluing to the Counterstain 1 and Counterstain 2 dispensers.
4. Follow the dispenser product insert instructions to properly prime and handle the dispensers.
5. Store tightly capped dispensers (except the mRNA Dewax dispenser) at 4°C when not in use.
6. Store tightly capped mRNA Dewax dispenser at room temperature when not in use.


### *Create instrument protocols for protein target qualification*

1. Open the NexES software and click the **Protocol** button.
2. Click **Create/Edit Protocols**, go to the Procedure drop-down menu and select **mRNA Universal**. Main protocol selections appear as shown:






3. To create a protocol to test the sensitivity of your protein target to RNAscope VS Protease, select the appropriate pretreatment conditions as shown in the following screenshot. After the main step selections, drop-down menus become available for further selection. For Cell Conditioning, we recommend the same conditions previously used on your tissue for ISH-only staining.



4. Select Dual Sequence to view and select IHC options.

- Select **DS Pretreatment** followed by **Protease-ISH** to enable use of the mRNA Protease dispenser outside of the mRNA assay selections. Select **37°C** and **16 MIN** as your protease conditions.



- For HRP detection, select **DS Inhibitor** followed by **Neutralize**. You do not need to select these options for AP detection.
- Select **DS Antibody**, followed by the desired primary antibody incubation time, temperature, and barcode according to your established IHC protocol.
- Choose a secondary antibody selection corresponding to your desired secondary detection method.  
*For example, if using an anti-rabbit HRP multimer, select **DS Multimer HRP** and choose **UMap anti-Rb HRP**. Select the secondary incubation time according to your established IHC protocol.*
- Select the appropriate detection kit and available incubation conditions according to your established IHC protocol. For a list of compatible IHC chromogen detection kits, see the table on page 10.
- Select your preferred counterstain/post-counterstain reagents and assay time.
- Click **Save As**, then select a unique protocol number from the drop-down menu and choose a protocol name. Click **Active**, add relevant comments in the available field, and click **Save**.
- To create a control protocol without the RNAscope VS Protease application, de-select **DS Pretreatment**. Click **Save As**, then select a unique protocol number from the drop-down menu and choose a protocol name. Click **Active**, add relevant comments in the available field, and click **Save**.


### Print the labels

- Select the **Print Label** icon from the upper right corner of the home screen.
- Select your preferred template or create a new template. To create a new template, refer to the *Ventana DISCOVERY ULTRA System User Manual* for details.
- Click **Protocol**.
- Select the protocol with the RNAscope VS Protease application and the control protocol without the RNAscope VS Protease application created in the previous section. Click the **Add** button after selecting each protocol.
- When the protocols for all slides have been assigned, click **Close/Print**.
- Fill in the template for each slide. Click **Print** when completed

### Load the reagents

- Remove the nozzle caps from the filled dispensers and place each cap on the post located on the back of the dispenser.
- Prime the user-fillable dispensers. For guidance, refer to the instructions provided by Roche Tissue Diagnostics.

3. If needed, remove any air bubbles at the nozzle tip by pushing down on the nozzle until the liquid reaches the tip of the nozzle or forms a small meniscus at the tip of the nozzle.
4. Remove the yellow locking ring from the dispensers in all the prefilled dispensers. Refer to the instructions provided by Roche Tissue Diagnostics.
5. Load the dispensers onto the reagent racks.
6. Load the reagent racks onto the reagent carousel.
7. Select the **Ready** button.




8. Open the slide drawers.
9. Load each slide onto a heater pad with the label facing upward and inward. Ensure that the slides sit securely on the pads.

---

**IMPORTANT!** Prior to loading the slides, ensure heater pads are completely dry. Wipe off any liquid using laboratory tissue paper.

---

10. Close the slide drawers.
11. Select the **Running** button.



12. The assay duration varies from **15 – 20 HRS** based on assay selections.

---

**IMPORTANT!** Before leaving the instrument unattended, ensure all reagents and slides are successfully registered and the instrument is running.

---

### *Complete the run*

1. After the run is complete, remove the Dewax reagent, place nozzle cap on the dispenser, and store at room temperature.
2. For the remaining reagents, place nozzle caps back on the dispensers and place racks onto the magnet locking tray.

---

**IMPORTANT!** Store reagent racks at 4°C until next use. Store the Dewax dispenser at room temperature.

---

### *Wash and dry the slides*

1. Prepare 200 mL of diluted detergent by adding 1–2 drops detergent to 200 mL distilled water in a container with a cap.
2. Mix well by inverting the container 4–5 times.
3. Add diluted detergent to a Tissue-Tek Staining Dish.

**Note:** Store diluted detergent at **RT**.

1. Submerge a Tissue-Tek Slide Rack into the Tissue-Tek Staining Dish containing 200 mL diluted detergent.
2. Open the instrument slide drawers and unload slides.

3. Decant solution on the slides into the slide drawer, then *immediately* load slides into the Tissue-Tek Slide Rack submerged in detergent.
4. Rinse oil off the slides by moving the slide rack up and down in the dish 10 times.
5. Replace the detergent with distilled water and rinse slides by moving the slide rack up and down a minimum of **10** times.
6. Repeat Step 8 three to five times.
7. Transfer the slides into a Tissue-Tek Staining Dish containing **200 mL** distilled water.
8. Place slides in a drying oven at **60°C** for at least **30 MIN**.

### *Mount the samples*

1. In a fume hood, fill two clearing agent dishes with ~200 mL fresh xylene.
2. Once slides are dry, move the Tissue-Tek Slide rack into the staining dish containing xylene for **1 MIN** with occasional agitation.
3. Move the Tissue-Tek Slide rack into the staining dish containing xylene for **1 MIN** with occasional agitation.
4. Lay each slide flat with the sections facing up in the fume hood then add 1–2 drops of EcoMount or other chromogen-compatible xylene-based mounting medium. Carefully place a 24 mm x 50 mm coverslip over the section and avoid trapping air bubbles.
5. Air dry slides for at least **15 MIN** before evaluation.

---

**IMPORTANT!** mRNA Teal, mRNA Green, DISCOVERY Teal HRP and DISCOVERY Green HRP chromogens are light sensitive and may fade over time. For best results, protect stored slides from the light and image within one week of staining. .

---

### *Interpret the results of protein target qualification:*

1. Compare the staining patterns with and without RNAscope VS Protease application.
2. If the staining is comparable, proceed to **Chapter 2: Sequential RNA-Protein Co-Detection**.
3. If protein staining intensity, cell positivity rate, or staining distribution is impacted by RNAscope VS Protease treatment, proceed to **Chapter 3: Integrated Co-Detection Workflow**.

## Chapter 2: Sequential RNA-Protein Co-Detection

In this workflow, tissue samples undergo RNAscope pretreatment including RNAscope VS Protease digestion, followed by ISH detection, and then undergo IHC detection and counterstaining. First complete **Chapter 1: Protein Target Qualification** to confirm your protein epitope(s) of interest are protease-tolerant prior to proceeding with this workflow.

### Materials

#### ACD reagents for Sequential RNA-Protein Co-Detection

Sequential RNA-Protein Co-Detection on the Roche DISCOVERY ULTRA requires the following reagents from ACD:

- RNAscope VS Universal Sample Prep Reagent Kit v2 (Cat. No. 323740)
- ACD Assay-specific Detection Reagent Kit:
  - RNAscope VS Universal AP Detection Reagents (Cat. No. 323260)
  - RNAscope VS Universal HRP Detection Reagents (Cat. No. 323210)
  - RNAscope VS Duplex Detection Reagents (Cat. No. 323310)
  - BaseScope VS Detection Reagents (Cat. No. 323710)
- VS Co-Detection Inhibitor from VS RNA-Protein Co-Detection Ancillary Kit (Cat. No. 323760)

RNAscope VS Universal Sample Prep Reagent Kit v2 (Cat. No. 323740) and RNAscope VS Accessory Kit (Cat. No. 320630) are RTU and stored as indicated in the following tables:

| RNAscope VS Sample Prep Reagent Kit v2 (Cat. No. 323740) |   |                  |                     |
|--|---|------------------|---------------------|
| <input checked="" type="checkbox"/>                      | Reagent                                   | Quantity         | Storage             |
|  | RNAscope VS Universal Target Retrieval v2 | 10 mL x 2 bottle | Room Temp (15–30°C) |
|  | RNAscope VS Universal Dewax               | 14 mL x 1 bottle | Room Temp (15–30°C) |

Assay-specific Detection Reagents within each Reagent Kit are stored as indicated in the following tables:

| RNAscope VS Universal AP Detection Reagents (Cat. No. 323260) |                                |                  |         |
|---|--------------------------------|------------------|---------|
| <input checked="" type="checkbox"/>                           | Reagent                        | Quantity         | Storage |
|   | RNAscope VS Universal AP AMP 1 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 2 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 3 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 4 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 5 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 6 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 7 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Protease           | 14 mL x 1 bottle | 2–8°C   |

| RNAscope VS Universal HRP Detection Reagents (Cat. No. 323210) |                                 |                  |         |
|--|---------------------------------|------------------|---------|
| <input checked="" type="checkbox"/>                            | Reagent                         | Quantity         | Storage |
|  | RNAscope VS Universal HRP AMP 1 | 14 mL x 1 bottle | 2–8°C   |
|  | RNAscope VS Universal HRP AMP 2 | 14 mL x 1 bottle | 2–8°C   |
|  | RNAscope VS Universal HRP AMP 3 | 14 mL x 1 bottle | 2–8°C   |
|  | RNAscope VS Universal HRP AMP 4 | 14 mL x 1 bottle | 2–8°C   |
|  | RNAscope VS Universal HRP AMP 5 | 14 mL x 1 bottle | 2–8°C   |
|  | RNAscope VS Universal HRP AMP 6 | 14 mL x 1 bottle | 2–8°C   |
|  | RNAscope VS Universal HRP AMP 7 | 14 mL x 1 bottle | 2–8°C   |
|  | RNAscope VS Protease            | 14 mL x 1 bottle | 2–8°C   |

| RNAscope VS Duplex Detection Reagents (Cat. No. 323310) |                             |                   |         |
|---|-----------------------------|-------------------|---------|
| <input checked="" type="checkbox"/>                     | Reagent                     | Quantity          | Storage |
|   | RNAscope VS Duplex AMP 1    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 2    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 3    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 4    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 5    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 6    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 7    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 8    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 9    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP Wash | 14 mL x 2 bottles | 2–8°C   |
|   | RNAscope VS Protease        | 14 mL x 1 bottle  | 2–8°C   |

| BaseScope VS Detection Reagents (Cat. No. 323710) |                                       |                  |         |
|---|---------------------------------------|------------------|---------|
| <input checked="" type="checkbox"/>               | Reagent                               | Quantity         | Storage |
|   | BaseScope VS AMP 1                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 2                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 3                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 4                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 5                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 6                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 7                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 8                    | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope 2.5 VS Pretreat 3 - Protease | 14 mL x 1 bottle | 2–8°C   |

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.

| VS RNA-Protein Co-Detection Ancillary Kit - Cat No. 323760 |  |                        |                   |         |
|--|--|------------------------|-------------------|---------|
| <input checked="" type="checkbox"/>                        | Reagent  | Source / Ordering Info | Quantity          | Storage |
|  | Co-Detection Antibody Diluent*                     | ACD/Cat No. 323160     | 120 mL x 1 bottle | 2–8°C   |
|  | VS RNA-Protein Co-Detection Protease and Inhibitor | ACD/Cat No. 323190     | See the following |         |
|  | VS Co-Detection Protease†                          |                        | 14 mL x 1 bottle  | 2–8°C   |
|  | VS Co-Detection Inhibitor‡                         |                        | 14 mL x 1 bottle  | 2–8°C   |

\*Use of Co-Detection Antibody Diluent is optional for Sequential RNA-Protein Co-Detection

†VS Co-Detection Protease is not needed for Sequential RNA-Protein Co-Detection

‡ VS Co-Detection Inhibitor is a chemical enzyme inhibitor recommended for use between any RNAscope/BaseScope assay and IHC detection. We recommend this reagent for both sequential and integrated co-detection workflows.

### Optional ACD reagents for Sequential RNA-Protein Co-Detection

The following components from ACD are optional for Sequential RNA-Protein Co-Detection on the Roche DISCOVERY ULTRA:

- Co-Detection Antibody Diluent from VS RNA-Protein Co-Detection Ancillary Kit (Cat No. 323760)
- RNAscope VS Accessory Kit (Cat. No. 320630)

| VS RNA-Protein Co-Detection Ancillary Kit - Cat No. 323760 |                               |                        |                   |         |
|--|-------------------------------|------------------------|-------------------|---------|
| <input checked="" type="checkbox"/>                        | Reagent                       | Source / Ordering Info | Quantity          | Storage |
|  | Co-Detection Antibody Diluent | ACD/Cat No. 323160     | 120 mL x 1 bottle | 2–8°C   |

| RNAscope VS Accessory Kit (Cat. No. 320630) |                            |                 |         |
|---|----------------------------|-----------------|---------|
| <input checked="" type="checkbox"/>         | Reagent                    | Quantity        | Storage |
|   | RNAscope VS Hematoxylin    | 7 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Bluing Reagent | 7 mL x 1 bottle | 2–8°C   |

### Roche materials for Sequential RNA-Protein Co-Detection:

Sequential RNA-Protein Co-Detection can be performed with either Roche RTU primary antibody or with your choice of primary antibody concentrate diluted in the ACD Co-Detection Antibody Diluent. For a list of available Roche RTU Primary Antibodies and ordering information, please contact your local Roche representative.

The following additional materials from Roche can be used for Sequential RNA-Protein Co-Detection on DISCOVERY ULTRA. Catalog numbers are valid in the United States only. For other regions, please check catalog or ordering numbers with your local lab supplier.



### Reagent dispensers

| Roche Materials for Sequential RNA-Protein Co-Detection |  |                            |                                    |  |
|---|--|----------------------------|------------------------------------|--|
| <input checked="" type="checkbox"/>                     | Component                                | Cat. No.                   | Ordering Code                      | Fill with:   |
|   | Probe Dispensers                         | 960-761 to 960-780         | Contact local Roche representative | RNAscope VS Probes   |
|   | mRNA Sample Prep Kit                     | 760-248                    | 08127166001                        | RNAscope VS Sample Prep Reagent Kit v2 and VS Protease reagent from Detection Kit  |
|   | Antibody Dispensers ( <i>Optional*</i> ) | 770-001 to 770-099         | Contact local Roche representative | User-sourced Primary Antibody Concentrate diluted in Co-Detection Antibody Diluent |
|   | Pretreatment 1 Dispenser                 | 960-901                    | 05280095001                        | VS Co-Detection Inhibitor  |
|   | Counterstain 1 dispenser                 | 771-741                    | 05271720001                        | VS Hematoxylin   |
|   | Counterstain 2 dispenser                 | 771-742                    | 05271738001                        | VS Bluing Reagent  |
|   | Roche RTU Counterstain Reagents          | Various                    | Contact local Roche representative | Pre-filled   |
|   | DISCOVERY Antibody Block                 | 760-4204                   | 05268869001                        | Pre-filled   |
|   | Additional ISH Assay-Specific Materials  | (See the following tables) |                                    |  |
|   | Additional IHC Detection Materials       | (See the following tables) |                                    |  |

\*Sequential RNA-Protein Co-Detection can be performed with either Roche RTU primary antibody or your choice of primary antibody.

### Additional ISH assay-specific materials

| Additional Roche Materials Required for VS Universal AP ISH |                                  |          |               |   |
|---|----------------------------------|----------|---------------|---|
| <input checked="" type="checkbox"/>                         | Component                        | Cat. No. | Ordering Code | Fill with:  |
|   | mRNA RED Probe Amplification Kit | 760-236  | 07095341001   | RNAscope VS Universal AP Detection Reagents AMP 1-7 |
|   | mRNA RED Detection Kit           | 760-234  | 07099037001   | Pre-filled  |

| Additional Roche Materials Required for VS Universal HRP ISH |                                |          |               |  |
|--|--------------------------------|----------|---------------|--|
| <input checked="" type="checkbox"/>                          | Component                      | Cat. No. | Ordering Code | Fill with:   |
|  | mRNA Probe Amplification Kit   | 760-222  | 06614337001   | RNAscope VS Universal HRP Detection Reagents AMP 1-7 |
|  | DISCOVERY Inhibitor*           | 760-4840 | 07017944001   | Pre-filled   |
|  | mRNA Purple HRP Detection Kit* | 760-255  | 08127166001   | Pre-filled   |
|  | mRNA Green HRP Detection Kit†  | 760-278  | 08952612001   | Pre-filled   |
|  | mRNA Teal HRP Detection Kit†   | 760-256  | 08352941001   | Pre-filled   |
|  | mRNA DAB Detection Kit†        | 760-224  | 06614353001   | Pre-filled   |
|  | DISCOVERY DCC Kit†             | 760-240  | 07988192001   | Pre-filled   |
|  | DISCOVERY FAM Kit†             | 760-243  | 07988150001   | Pre-filled   |
|  | DISCOVERY FITC Kit†            | 760-232  | 07259212001   | Pre-filled   |
|  | DISCOVERY Rhodamine Kit†       | 760-233  | 07259883001   | Pre-filled   |
|  | DISCOVERY Rhodamine 6G Kit†    | 760-244  | 07988168001   | Pre-filled   |

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.



| Additional Roche Materials Required for VS Universal HRP ISH |                         |         |             |            |
|--|-------------------------|---------|-------------|------------|
|  | DISCOVERY Red 610 Kitt† | 760-245 | 07988176001 | Pre-filled |
|  | DISCOVERY Cy5 Kitt†     | 760-238 | 07551215001 | Pre-filled |

\* DISCOVERY Inhibitor is not required if using mRNA DAB Detection Kit

† Choose one mRNA Detection Kit for RNAscope HRP detection

| Additional Roche Materials Required for VS Duplex ISH |                               |          |               |   |
|---|-------------------------------|----------|---------------|---|
| <input checked="" type="checkbox"/>                   | Component                     | Cat. No. | Ordering Code | Fill with:  |
|   | mRNA Duplex Amp Kit           | 760-249  | 08127174001   | RNAscope VS Duplex Detection Reagents AMP 1–9, AMP Wash |
|   | mRNA Link (Pre-filled)        | 760-6014 | 08127115001   | Pre-filled  |
|   | mRNA RED Detection Kit        | 760-234  | 07099037001   | Pre-filled  |
|   | mRNA Green HRP Detection Kit* | 760-278  | 08952612001   | Pre-filled  |
|   | mRNA Teal HRP Detection Kit*  | 760-256  | 08352941001   | Pre-filled  |
|   | mRNA DAB Detection Kit*       | 760-224  | 06614353001   | Pre-filled  |

\* Choose one mRNA Detection Kit for VS Duplex Channel 1 detection

| Additional Roche Materials Required for VS BaseScope ISH |                                  |          |               |   |
|--|----------------------------------|----------|---------------|---|
| <input checked="" type="checkbox"/>                      | Component                        | Cat. No. | Ordering Code | Storage                                 |
|  | mRNA RED Probe Amplification Kit | 760-236  | 07095341001   | BaseScope VS Detection Reagents AMP 1–7 |
|  | Option 8 dispenser               | 771-758  | 05271916001   | BaseScope VS AMP 8 reagent              |
|  | mRNA RED Detection Kit           | 760-234  | 07099037001   | Pre-filled                              |

### Additional IHC assay-specific materials

| Reagent Options for Ventana IHC AP Detection |                                      |          |               |         |
|--|--------------------------------------|----------|---------------|---------|
| <input checked="" type="checkbox"/>          | Component                            | Cat. No. | Ordering Code | Storage |
|  | DISCOVERY UltraMap anti-Ms Alk Phos* | 760-4312 | 05269687001   | 2–8°C   |
|  | DISCOVERY UltraMap anti-Rb Alk Phos* | 760-4314 | 05269709001   | 2–8°C   |
|  | DISCOVERY Yellow Kit*                | 760-239  | 07698445001   | 2–8°C   |
|  | DISCOVERY Red Kit                    | 760-228  | 07425333001   | 2–8°C   |
|  | DISCOVERY ChromoMap Red Kit*         | 760-160  | 05266653001   | 2–8°C   |

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.

| Reagent Options for Ventana IHC HRP Detection |                                  |          |               |         |
|---|----------------------------------|----------|---------------|---------|
| <input checked="" type="checkbox"/>           | Product                          | Cat. No. | Ordering Code | Storage |
|   | DISCOVERY UltraMap anti-Ms HRP*  | 760-4313 | 05269695001   | 2–8°C   |
|   | DISCOVERY UltraMap anti-Rb HRP*  | 760-4315 | 05269717001   | 2–8°C   |
|   | DISCOVERY UltraMap anti-Rat HRP* | 760-4456 | 05891884001   | 2–8°C   |
|   | DISCOVERY UltraMap anti-Gt HRP*  | 760-4648 | 06607241001   | 2–8°C   |
|   | DISCOVERY Purple Kit*            | 760-229  | 07053983001   | 2–8°C   |
|   | DISCOVERY Green HRP Kit †        | 760-278  | 07053983001   | 2–8°C   |
|   | DISCOVERY Teal HRP Kit †         | 760-247  | 08254338001   | 2–8°C   |
|   | DISCOVERY ChromoMap DAB Kit †    | 760-159  | 05266645001   | 2–8°C   |
|   | DISCOVERY DCC Kit †              | 760-240  | 07988192001   | 2–8°C   |
|   | DISCOVERY FAM Kit †              | 760-243  | 07988150001   | 2–8°C   |
|   | DISCOVERY FITC Kit†              | 760-232  | 07259212001   | 2–8°C   |
|   | DISCOVERY Rhodamine Kit †        | 760-233  | 07259883001   | 2–8°C   |
|   | DISCOVERY Rhodamine 6G Kit †     | 760-244  | 07988168001   | 2–8°C   |
|   | DISCOVERY Red 610 Kit†           | 760-245  | 07988176001   | 2–8°C   |
|   | DISCOVERY Cy5 Kit †              | 760-238  | 07551215001   | 2–8°C   |

\* Choose one secondary detection antibody depending on the primary antibody species and desired IHC chromogen

† Choose one IHC chromogen

### Instrument buffers

| <input checked="" type="checkbox"/> | Component                | Cat. No. | Ordering Code | Storage             |
|-------------------------------------|--------------------------|----------|---------------|---------------------|
|                                     | 10X DISCOVERY Wash (RUO) | 950-510  | 7311079001    | Room Temp (15–30°C) |
|                                     | ULTRA LCS (Predilute)    | 650-210  | 5424534001    | Room Temp (15–30°C) |
|                                     | SSC Buffer (10X)         | 950-110  | 5353947001    | Room Temp (15–30°C) |
|                                     | Reaction Buffer (10X)    | 760-107  | 5266262001    | Room Temp (15–30°C) |
|                                     | DISCOVERY CC1            | 950-500  | 6414575001    | Room Temp (15–30°C) |

### User-Sourced Materials

| <input checked="" type="checkbox"/> | Description                       | Supplier          | Cat. No.  |
|-------------------------------------|-----------------------------------|-------------------|-----------|
|                                     | Primary Antibody (RTU)            | User              | Various   |
|                                     | Primary Antibody Concentrate      | User              | Various   |
|                                     | SuperFrost Plus Slides (required) | Fisher Scientific | 12-550-15 |
|                                     | ProLong Gold Antifade Reagent     | Life Technologies | P36930    |
|                                     | EcoMount (if using Red detection) | Biocare           | EM897L    |

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.

| <input checked="" type="checkbox"/> | Description   | Supplier                               | Cat. No.     |
|-------------------------------------|---|--|--------------|
|                                     | Tissue-Tek Vertical 24 Slide Rack                             | American Master<br>Tech Scientific/MLS | LWSRA24      |
|                                     | Tissue-Tek Staining Dishes                                    | American Master<br>Tech Scientific/MLS | LWT4457EA    |
|                                     | Cover Glass 24 x 50 mm  | Fisher<br>Scientific/MLS               | 12--545-F    |
|                                     | Distilled water   | MLS                                    | —            |
|                                     | Mild liquid dishwashing detergent (Dawn detergent or similar) | MLS                                    | —            |
|                                     | Drying oven, capable of holding temperature at 60 +/- 1°C     | MLS                                    | —            |
|                                     | Fume hood   | MLS                                    | —            |
|                                     | 100% ethanol (EtOH)   | MLS                                    | —            |
|                                     | Xylene  | MLS                                    | —            |
|                                     | Tissue-Tek Clearing Agent Dishes, xylene-resistant            | American Master<br>Tech Scientific/MLS | LWT4456EA    |
|                                     | Optional: Glass beaker (1 or 2 L)                             | MLS                                    | —            |
|                                     | Optional: Hot plate   | Fisher<br>Scientific/MLS               | 11-300-49SHP |

## Assay Procedure

### Prepare the DISCOVERY ULTRA

#### Prepare the instrument

If the instrument has not been used for > 1 week, follow guidelines for instrument maintenance from Roche Tissue Diagnostics. Before use, empty the waste carboys if needed.

#### Dilute instrument bulk reagents

1. Prepare the instrument bulk fluids according to the manufacturer's instructions.
2. Fill bulk solution containers for 1X DISCOVERY Wash, ULTRA LCS (predilute), and CC1 (predilute) to be at least half full. Fully fill bulk solution containers for 2X SSC and 1X Reaction Buffer.

---

**IMPORTANT!** Do not use expired reagents.

---

#### Register new reagents

Reagent dispensers come with appropriate barcode labels and registration buttons for dispensing RNAscope VS Universal Reagents. Refer to the *Ventana DISCOVERY ULTRA System User Manual* for details. To register reagents:

- Log all ACD reagents and probes into the software as **log user-fillable reagents** and **log user-**

**fillable probes**, respectively.

- Use the reagent registration wand that comes with the instrument to register new reagent kits from Roche Tissue Diagnostics

### *Prepare user-fillable reagents for RNA-Protein Co-Detection*

Refer to the table on pages 6–7 to determine the proper dispenser for each reagent.

---


**IMPORTANT!** Avoid cross contamination between reagents. Dewax must be warmed to room temperature and be completely in solution before use.

---


1. Transfer the entire volume of each AMP component of the Detection Kit to the corresponding labeled dispenser from the appropriate mRNA Amplification kit (see pages 6–7 for details)
2. Fill the mRNA Sample Prep Kit:
  - a. Transfer the contents of both bottles of VS Target Retrieval v2 from the RNAscope VS Sample Prep Reagent Kit v2 to the mRNA Target Retrieval Dispenser.
  - b. Transfer the VS Dewax reagent from the RNAscope VS Sample Prep Reagent Kit v2 to the mRNA Dewax Dispenser.
  - c. Transfer the VS Protease from the Detection Kit to the mRNA Protease Dispenser
3. Fill the User Fillable Dispensers:
  - a. Transfer the RNAscope 2.5 or BaseScope VS Target Probe and control probes to the corresponding probe dispensers.
  - b. Transfer the VS Co-Detection Inhibitor from the VA RNA-Protein Co-Detection Ancillary kit to the Pretreatment 1 dispenser.
  - c. Transfer the VS Hematoxylin and VS Bluing to the Counterstain 1 and Counterstain 2 dispensers.
4. Follow the dispenser product insert instructions to properly prime and handle the dispensers.
5. Store tightly capped dispensers (except the mRNA Dewax dispenser) at 4°C when not in use.
6. Store mRNA Dewax dispenser at room temperature when not in use.

### *Create an instrument protocol*

1. Open the NexES software and click the **Protocol** button.
2. Click **Create/Edit Protocols**, go to the Procedure drop-down menu and select **mRNA Universal**. Main protocol selections appear as shown:



3. Select the appropriate pretreatment conditions as shown in the following screenshot. After the main step selections, drop-down menus become available for further selection. For Cell Conditioning, we recommend the same conditions previously used on your tissue for ISH-only staining.



4. Select one mRNA Detection assay of interest.
  - For BaseScope VS, select **BaseScope**.
  - For RNAscope VS Duplex, select **mRNA Duplex**.
  - For RNAscope VS Universal AP, select **mRNA AP Detection**.
  - For RNAscope VS Universal HRP, select **mRNA HRP Detection**.

Select only one mRNA option per protocol. Choosing multiple selections yields negative results.

[ Only select one mRNA option as multiple selections will yield negative results ]

☐ BaseScope

☐ mRNA Duplex

☐ mRNA AP Detection

☐ mRNA HRP Detection

☐ Dual Sequence

☐ Counterstain

☐ Slide Cleaning

5. Select the appropriate Detection Inhibitor:

- For **BaseScope** or **mRNA AP Detection**, mRNA Red Inhibitor reagent is applied automatically.
- For **mRNA Duplex**, select the **Inhibitor** check box to display the available selections. Choose the appropriate inhibitor selection based on the chromogen combination in use for Duplex ISH detection, **Red-DAB Inhibitor**, **Red-Green Inhibitor**, or **Red-Teal Inhibitor**.

[ Only select one mRNA option as multiple selections will yield negative results ]

☐ BaseScope

☒ mRNA Duplex

[ For mRNA Duplex AP / HRP Detection ]

☒ Inhibitor

☐ Red-DAB Inhibitor

☐ Red-Green Inhibitor

☐ Red-Teal Inhibitor

☐ mRNA Duplex 3rd Pretreatment

☐ mRNA Duplex ICW Protease

[ Target Probe Cocktail ]

[ RECOMMENDED: Set to 43 C for probe incubation ]

☒ Apply Two Drops of [ Probe ] ( Probe #1 ), Apply Coverslip, and Incubate for 4 Minutes

☒ Warmup Slide to [ Low Temperature ], and Incubate for 2 Hours ( Hybridization )

- For **mRNA HRP Detection**, select mRNA HRP Detection Inhibitor if using the mRNA DAB Detection Kit. Select DISCOVERY Inhibitor if using an alternative HRP-based chromogen such as the mRNA Purple Detection Kit, mRNA Green Detection Kit, or mRNA Teal Detection Kit.

[ Only select one mRNA option as multiple selections will yield negative results ]

☐ BaseScope

☐ mRNA Duplex

☐ mRNA AP Detection

☒ mRNA HRP Detection

[ Select an Inhibitor NOTE: DISCOVERY Inhibitor to be used with fluorescence ]

☒ mRNA HRP Detection Inhibitor

[ Inhibitor for mRNA DAB Detection will be applied ]

☐ DISCOVERY Inhibitor

☐ HRP detection 3rd Pretreatment



6. Under the mRNA assay of interest, select the corresponding 3<sup>rd</sup> **Pretreatment** check box. Do not select the detection ICW Protease option. Recommended incubation is 37°C for 16 MIN.

[ Only select one mRNA option as multiple selections will yield negative results ]

☐ BaseScope

☐ mRNA Duplex

☒ mRNA AP Detection

[ Inhibitor for mRNA Red Detection will be applied ]

☒ AP detection 3rd Pretreatment

[ RECOMMENDED: Set to 37 C and 16min for normal FFPE samples ]

[ Protease ]

✓ Warmup Slide to [ 37 Deg C ], and Incubate for [ 0 Hr 16 Min ] ( Pretreatment #3 Temp RB )

☐ AP Detection ICW Protease

[ Target Probe: RECOMMENDED: Set to 43 C for probe incubation ]

✗ Apply Two Drops of [ Probe ] ( Probe #1 ), Apply Coverslip, and Incubate for 4 Minutes

✗ Warmup Slide to [ Low Temperature ], and Incubate for 2 Hours ( Hybridization )

7. Select the ISH staining conditions:
  - If running **BaseScope**, select the following ISH staining conditions:

[ Only select one mRNA option as multiple selections will yield negative results ]

☒ BaseScope

[ BaseScope Utilizes Standard Pretreatment ]

☒ BaseScope 3rd Pretreatment

[ RECOMMENDED: Set to 37 C and 16min for normal FFPE samples ]

✓ Warmup Slide to [ 37 Deg C ], and Incubate for [ 16 Minutes ] ( Pretreatment #3 Temp RB )

☐ BaseScope ICW Protease

[ Specify Probe, hybridization time and temperature ]

✓ Apply Two Drops of [ PROBE 1 ] ( Probe #1 ), Apply Coverslip, and Incubate for 4 Minutes

[ BaseScope Probe - Temp = 50 C Recommended ]

✓ Warmup Slide to [ 50 Deg C ], and Incubate for 2 Hours ( Hybridization )

[ BaseScope Amp 1 - Temp = 49 C Recommended ]

✓ Warmup Slide to [ 49 Deg C ], and Incubate for 32 Minutes ( Detection #1 Temp )

[ BaseScope Amp 2 - Temp = 49 C Recommended ]

✓ Warmup Slide to [ 49 Deg C ], and Incubate for 32 Minutes ( Detection #2 Temp )

[ BaseScope Amp 4 - Temp = 49 C Recommended ]

✓ Warmup Slide to [ 49 Deg C ], and Incubate for 32 Minutes ( Detection #3 Temp )

[ BaseScope Amp 5 - Temp = 39 C Recommended ]

✓ Warmup Slide to [ 39 Deg C ], and Incubate for 32 Minutes ( Detection #4 Temp )

[ BaseScope Amp 7 - Recommended 4min ]

✓ Incubate for [ 0 Hr 4 Min ]

☐ mRNA Duplex

☐ mRNA AP Detection

☐ mRNA HRP Detection

☐ Dual Sequence

| Standard Temperatures/Times for BaseScope VS |              |
|--|--------------|
| VS Protease temperature and time             | 37°C, 16 MIN |
| Standard BaseScope probe temperature         | 50°C         |
| Standard BaseScope AMP 1 temperature         | 49°C         |
| Standard BaseScope AMP 2 temperature         | 49°C         |
| Standard BaseScope AMP 4 temperature         | 49°C         |
| Standard BaseScope AMP 5 temperature         | 39°C         |
| Standard BaseScope AMP 7 incubation time*    | 4 MIN        |
| mRNA Red Chromogen                           | Default†     |

\* BaseScope Amp 7 incubation time is determined by instrument calibration. Use the instrument setting previously optimized for the mRNA Universal software. For assistance, consult your local ACD FAS for more information.

†For all BaseScope assays, mRNA Red Detection is applied automatically. No chromogen selections are needed.

- If running **mRNA Duplex**, select the following ISH staining conditions. For the purpose of demonstration, selections for the Red-Green chromogen combination are shown:

☐ BaseScope

☒ mRNA Duplex

[ For mRNA Duplex AP / HRP Detection ]

☒ Inhibitor

☐ Red-DAB Inhibitor

☒ Red-Green Inhibitor

☐ Red-Teal Inhibitor

☒ mRNA Duplex 3rd Pretreatment

[ RECOMMENDED: Set to 37 C and 16min for normal FFPE samples ]

[ Protease ]

✓ Warmup Slide to [ 37 Deg C ], and Incubate for [ 0 Hr 16 Min ] ( Pretreatment #3 Temp RB )

☐ mRNA Duplex ICW Protease

[ Target Probe Cocktail ]

[ RECOMMENDED: Set to 43 C for probe incubation ]

✓ Apply Two Drops of [ PROBE 1 ] ( Probe #1 ), Apply Coverslip, and Incubate for 4 Minutes

✓ Warmup Slide to [ 43 Deg C ], and Incubate for 2 Hours ( Hybridization )

[ Amp 1 Duplex ]

[ RECOMMENDED: Temp = 39 C for most samples ]

✓ Warmup Slide to [ 39 Deg C ], and Incubate for 32 Minutes ( Hybridization #1 )

[ Amp 2 Duplex ]

[ RECOMMENDED: Temp = 39 C for most samples ]

✓ Warmup Slide to [ 39 Deg C ], and Incubate for 32 Minutes ( Hybridization #3 )

[ Amp 5 incubation time: RECOMMENDED is 4 min ]

✓ Incubate for [ 0 Hr 4 Min ] ( Hybridization #5 )

[ Amp 8 incubation time: RECOMMENDED: 4 min ]

✓ Incubate for [ 0 Hr 4 Min ] ( Hybridization #6 )

☐ Red-DAB AP-HRP Duplex

☐ Red-Teal AP-HRP Duplex

☒ Red-Green AP-HRP Duplex

[ Recommended Settings are 20 min for both Green H2O2 and Green Activator ]

✓ Apply One Drop of mRNA Green H2O2, and Incubate for [ 20 Minutes ]

✓ Apply One Drop of mRNA Green Act, and Incubate for [ 20 Minutes ]

☐ mRNA AP Detection



| Standard Temperatures/Times for mRNA Duplex Detection  |  |
|--|--|
| VS Protease temperature and time                       | 37°C, 16 MIN   |
| Suggested probe temperature                            | 43°C   |
| Suggested RNAscope VS Duplex AMP 1 & AMP 2 temperature | 39°C   |
| RNAscope VS Duplex AMP 5 incubation time               | 4 MIN*   |
| RNAscope VS Duplex AMP 8 incubation time               | 4 MIN*   |
| Chromogen settings for Red-DAB AP-HRP Duplex†          | Preset, no incubation selections   |
| Chromogen settings for Red-Green AP-HRP Duplex†        | mRNA Green H <sub>2</sub> O <sub>2</sub> — 20 MIN<br>mRNA Green Act — 20 MIN |
| Chromogen settings for Red-Teal AP-HRP Duplex†         | mRNA Teal H <sub>2</sub> O <sub>2</sub> — 20 MIN<br>mRNA Teal Act — 20 MIN   |

\* VS Duplex Amp 5 and Amp8 incubation times are determined by instrument calibration. Use the instrument setting previously optimized for the mRNA Universal software. For assistance, consult your local ACD FAS for more information.

†Choose one detection combination for mRNA Duplex. All times/temperatures are preset for Red-DAB Detection.

- If choosing **mRNA AP Detection**, select the following ISH staining conditions:

[ Only select one mRNA option as multiple selections will yield negative results ]

☐ BaseScope

☐ mRNA Duplex

☒ mRNA AP Detection

[ Inhibitor for mRNA Red Detection will be applied ]

☒ AP detection 3rd Pretreatment

[ RECOMMENDED: Set to 37 C and 16min for normal FFPE samples ]

[ Protease ]

✓ Warmup Slide to [ 37 Deg C ], and Incubate for [ 0 Hr 16 Min ] ( Pretreatment #3 Temp RB )

☐ AP Detection ICW Protease

[ Target Probe: RECOMMENDED: Set to 43 C for probe incubation ]

✓ Apply Two Drops of [ PROBE 1 ] ( Probe #1 ), Apply Coverslip, and Incubate for 4 Minutes

✓ Warmup Slide to [ 43 Deg C ], and Incubate for 2 Hours ( Hybridization )

[ Amp 1 AP ]

[ RECOMMENDED: Temp = 39 C for most samples ]

✓ Warmup Slide to [ 39 Deg C ], and Incubate for 32 Minutes ( Hybridization #2 )

[ Amp 2 AP ]

[ RECOMMENDED: Temp = 39 C for most samples ]

✓ Warmup Slide to [ 39 Deg C ], and Incubate for 32 Minutes ( Hybridization #4 )

[ Amp 5 AP incubation time: RECOMMENDED: 4 min ]

✓ Incubate for [ 0 Hr 4 Min ] ( Hybridization #5 )

☐ mRNA HRP Detection

☐ Dual Sequence

| Standard Temperatures/Times for mRNA AP Detection |              |
|---|--------------|
| VS Protease temperature and time                  | 37°C, 16 MIN |
| Suggested probe temperature                       | 43°C         |
| Suggested RNAscope AP AMP 1 & AMP 2 temperature   | 39°C         |
| RNAscope AP AMP 5 incubation time                 | 4 MIN*       |
| Chromogen settings for mRNA Red Detection         | Default†     |

\* VS Universal AP Amp 5 incubation time is determined by instrument calibration. Use the instrument setting previously optimized for the mRNA Universal software. For assistance, consult your local ACD FAS for more information.

†For all mRNA AP, mRNA Red Detection is applied automatically. No chromogen selections are needed.

- If choosing **mRNA HRP Detection**, select the following ISH staining conditions:

[ Only select one mRNA option as multiple selections will yield negative results ]

☐ BaseScope

☐ mRNA Duplex

☐ mRNA AP Detection

☒ mRNA HRP Detection

[ Select an Inhibitor NOTE: DISCOVERY Inhibitor to be used with fluorescence ]

☒ mRNA HRP Detection Inhibitor

[ Inhibitor for mRNA DAB Detection will be applied ]

☐ DISCOVERY Inhibitor

☒ HRP detection 3rd Pretreatment

[ RECOMMENDED: Set to 37 C and 16min for normal FFPE samples ]

[ Protease ]

✓ Warmup Slide to [ 37 Deg C ], and Incubate for [ 0 Hr 16 Min ] ( Pretreatment #3 Temp RB )

☐ HRP Detection ICW Protease

[ Target Probe ]

[ Target Probe: RECOMMENDED: Set to 43 C for probe incubation ]

✓ Apply Two Drops of [ PROBE 1 ] ( Probe #1 ), Apply Coverslip, and Incubate for 4 Minutes

✓ Warmup Slide to [ 43 Deg C ], and Incubate for 2 Hours ( Hybridization )

[ Amp 1 HRP ]

[ RECOMMENDED: Temp = 39 C for most samples ]

✓ Warmup Slide to [ 39 Deg C ], and Incubate for 32 Minutes ( Hybridization #5 )

[ Amp 2 HRP ]

[ RECOMMENDED: Temp = 39 C for most samples ]

✓ Warmup Slide to [ 39 Deg C ], and Incubate for 32 Minutes ( Hybridization #6 )

[ Amp 5 HRP incubation time: RECOMMENDED: 4 min ]

✓ Incubate for [ 0 Hr 4 Min ] ( Hybridization #6 )

[ Default detection is mRNA DAB unless a fluor or chromogen is selected ]

☐ Rhodamine

☐ FITC

☐ Rhodamine 6G

☐ Cy5

☐ FAM

☐ Red 610

☐ DCC

☐ mRNA Purple

☐ mRNA Teal

☐ mRNA Green


☐ Dual Sequence

| Standard Temperatures/Times for mRNA HRP Detection |  |
|--|--|
| VS Protease temperature and time                   | 37°C, 16 MIN   |
| Suggested probe temperature                        | 43°C   |
| Suggested RNAscope HRP AMP 1 & AMP 2 temperature   | 39°C   |
| RNAscope HRP AMP 5 incubation time                 | 4 MIN*   |
| Chromogen settings for mRNA DAB†                   | Default*   |
| Chromogen settings for mRNA Purple                 | mRNA Purple H <sub>2</sub> O <sub>2</sub> — 40 MIN                           |
| Chromogen settings for mRNA Teal                   | mRNA Teal H <sub>2</sub> O <sub>2</sub> — 20 MIN<br>mRNA Teal Act — 20 MIN   |
| Chromogen settings for mRNA Green                  | mRNA Green H <sub>2</sub> O <sub>2</sub> — 20 MIN<br>mRNA Green Act — 20 MIN |

\* VS Universal HRP Amp 5 incubation time is determined by instrument calibration. Use the instrument setting previously optimized for the mRNA Universal software. For assistance, consult your local ACD FAS for more information.

†For mRNA HRP, mRNA DAB Detection is applied by default if no other chromogens are selected.

8. After completing ISH assay selections, select **Dual Sequence** to enable selections for secondary antibody detection.



☒ **Dual Sequence**

☐ **Antibody Denaturation**

☐ **DS Pretreatment**

☐ **DS Option**

☐ **DS A/B Block**

☐ **DS Inhibitor**

**[ DISC Inhibitor will not be applied ]**

☐ **DS ISH**

☐ **DS Antibody**

☐ **Disable heat for 2nd Fixative**

☐ **DS Post-Antibody Fixative**

☐ **DS Antibody Amplification**

☐ **DS Linking Antibody**

☐ **Disable heat for 2nd Detection**

☐ **DS Multimer HRP**

☐ **DS Multimer AP**

☐ **DS Proximity Detection**

☐ **DS Enzyme conjugate**

☐ **DS DISCOVERY Amplification**

☐ **DS DAB**

☐ **DS Silver**

☐ **DS DISCOVERY Purple**

☐ **DS DISCOVERY Yellow HRP**

☐ **DS DISCOVERY Teal HRP**

☐ **DS DISCOVERY Green HRP**

☐ **DS DISCOVERY Red HRP**

☐ **DS DISCOVERY Blue HRP**

☐ **DS Red**

☐ **DS DISCOVERY Red**

☐ **DS DISCOVERY Yellow**

☐ **DS Blue**

☐ **DS Rhodamine**

☐ **DS FITC**

☐ **DS Cy5**

☐ **DS DCC**

☐ **DS FAM**

☐ **DS Red 610**

☐ **DS Rhodamine 6G**

☐ **DS Open Detection Kit**

☐ **Triple Stain**

☐ **Counterstain**

9. To apply VS Co-Detection Inhibitor from the VS RNA-Protein Co-Detection Ancillary Kit, select **DS Pretreatment**, **DS 2<sup>nd</sup> Pretreatment**, **Use DW for DS 2<sup>ns</sup> Pretreatment**. Select the Pretreatment 1 barcode to correspond to the VS Co-Detection Inhibitor reagent. Recommended incubation is **40°C for 32 MIN** prior to protein detection.



☒ **Dual Sequence**

☐ Antibody Denaturation

☒ **DS Pretreatment**

☐ DS 1st Pretreatment

☐ Protease-ISH

☐ DS Cell Conditioning

☒ **DS 2nd Pretreatment**

☐ Use RB for DS 2nd Pretreatment

☒ **Use DW for DS 2nd Pretreatment**

☐ DS Disable heat PT2-EZ

[ RECOMMENDED: If using temperatures above 42°C, set incubation for less than 1 hour ]

✓ Warmup Slide to [ **40 Deg C** ], and Incubate for **4 Minutes** ( **DS Pretreatment #2 Temp DW** )

✓ Apply One Drop of [ **PRETREATMENT 1** ] ( **DS Pretreatment #2** ), and Incubate for [ **0 Hr 32 Min** ]

☐ DS Enzyme

☐ DS Option

☐ DS A/B Block

☐ DS Inhibitor

10. To apply Primary antibody, select **DS Antibody** and your desired incubation time and temperature, using previously optimized IHC-only conditions.

☒ **Dual Sequence**

☐ Antibody Denaturation

☐ DS Pretreatment

☐ DS Option

☐ DS A/B Block

☐ DS Inhibitor

[ DISC Inhibitor will not be applied ]

☐ DS ISH

☒ **DS Antibody**

☐ DS Antibody Manual Application

☐ DS Antibody Blocking

☐ Disable heat for DS Antibody

☐ DS High Temp Ab incubation

✓ Warmup Slide to [ **37 Deg C** ] from **Very Low Temperatures** ( **DS Primary Antibody** )

☐ DS Extended Ab incubation

✓ Apply One Drop of [ **ANTIBODY 1** ] ( **DS Antibody** ), and Incubate for [ **32 Minutes** ]

☐ Disable heat for 2nd Fixative

11. Select the settings for secondary detection:

**Note:** We recommend using Roche's Antibody Block reagent directly before applying the secondary antibody to prevent non-specific cross-detection of RNAscope. RNA-Protein Co-Detection performed without Antibody Block could result in a hue shift of the RNA dots and interfere with interpretation of RNA-protein colocalization.

- For **HRP-based secondary detection**, choose a Roche Multimer HRP reagent corresponding to the primary host species. To enable secondary incubation, select **DS Multimer HRP**. To apply Antibody Block, select **DS Multimer HRP Blocker**, then select Antibody Block reagent. Select your Multimer HRP reagent of choice and desired incubation time, using previously optimized IHC-only conditions.

☐ Disable heat for 2nd Detection

☒ **DS Multimer HRP**

☒ **DS Multimer HRP Blocker**

[ Select Multimer blocker and Multimer species ]

[ Note: Recommended Multimer HRP Reagent incubation time for ICW is 32min ]

✓ Apply One Drop of [ [Antibody Block](#) ] ( **DS Multimer HRP Blocking** ), No Coverslip and Incubate for **32 Minutes**

✓ Apply One Drop of [ [UMap anti-Ms HRP](#) ] ( **DS Multimer HRP** ), and Incubate for [ [32 Minutes](#) ]

☐ **DS Multimer AP**

☐ **DS Proximity Detection**

☐ **DS Enzyme conjugate**

☐ **DS DISCOVERY Amplification**

☐ **DS DAB**

- For AP-based secondary detection, choose a Roche Multimer AP reagent corresponding to the primary host species. To enable secondary incubation, select **DS Multimer AP**. To apply Antibody Block, select **DS Multimer AP Blocker**, then select Antibody Block reagent. Select your Multimer AP reagent of choice and desired incubation time, using previously optimized IHC-only conditions.

☐ Disable heat for 2nd Detection

☐ **DS Multimer HRP**

☒ **DS Multimer AP**

☒ **DS Multimer AP Blocker**

[ Select Multimer blocker and Multimer species ]

[ Note: Recommended Multimer AP Reagent incubation time for ICW is 32min ]

✓ Apply One Drop of [ [Antibody Block](#) ] ( **DS Multimer AP Blocking #1** ), No Coverslip and Incubate for **32 Minutes**

✓ Apply One Drop of [ [UMap anti-Ms AP](#) ] ( **DS Multimer AP** ), and Incubate for [ [32 Minutes](#) ]

[ Select Multimer blocker and Multimer species ]

☐ **DS Proximity Detection**

☐ **DS Enzyme conjugate**

☐ **DS DISCOVERY Amplification**

☐ **DS DAB**

- Select chromogens for IHC detection using the following recommendations:

|  |                                       |
|--|---------------------------------------|
| <input type="checkbox"/> <b>DS Enzyme conjugate</b>        | <b>HRP Chromogens</b>                 |
| <input type="checkbox"/> <b>DS DISCOVERY Amplification</b> |                                       |
| <input checked="" type="checkbox"/> <b>DS DAB</b>          |                                       |
| <input type="checkbox"/> <b>DS Silver</b>                  |                                       |
| <input type="checkbox"/> <b>DS DISCOVERY Purple</b>        |                                       |
| <input type="checkbox"/> <b>DS DISCOVERY Yellow HRP</b>    | <b>AP Chromogens</b>                  |
| <input type="checkbox"/> <b>DS DISCOVERY Teal HRP</b>      |                                       |
| <input type="checkbox"/> <b>DS DISCOVERY Green HRP</b>     |                                       |
| <input type="checkbox"/> <b>DS DISCOVERY Red HRP</b>       |                                       |
| <input type="checkbox"/> <b>DS DISCOVERY Blue HRP</b>      |                                       |
| <input type="checkbox"/> <b>DS Red</b>                     | <b>DISCOVERY HRP Fluorescent Kits</b> |
| <input type="checkbox"/> <b>DS DISCOVERY Red</b>           |                                       |
| <input type="checkbox"/> <b>DS DISCOVERY Yellow</b>        |                                       |
| <input type="checkbox"/> <b>DS Blue</b>                    |                                       |
| <input type="checkbox"/> <b>DS Rhodamine</b>               |                                       |
| <input type="checkbox"/> <b>DS FITC</b>                    |                                       |
| <input type="checkbox"/> <b>DS Cy5</b>                     |                                       |
| <input type="checkbox"/> <b>DS DCC</b>                     |                                       |
| <input type="checkbox"/> <b>DS FAM</b>                     |                                       |
| <input type="checkbox"/> <b>DS Red 610</b>                 |                                       |
| <input type="checkbox"/> <b>DS Rhodamine 6G</b>            |                                       |
| <input type="checkbox"/> <b>DS Open Detection Kit</b>      |                                       |


| IHC Chromogen Selection and Settings |                            |   |
|--------------------------------------|----------------------------|---|
| IHC Enzyme                           | IHC Chromogen              | Recommended Chromogen Settings*   |
| AP                                   | DISCOVERY Yellow           | 44 MIN – 2 HRS  |
|                                      | DISCOVERY Red              | 12 MIN†   |
|                                      | DS Red / ChromoMap Red‡    | Default   |
| HRP<br>(chromogenic)                 | DS DAB / ChromoMap DAB     | Default   |
|                                      | DISCOVERY Purple           | 40 MIN  |
|                                      | DISCOVERY Teal HRP         | DISCOVERY Teal H <sub>2</sub> O <sub>2</sub> – 16–32 MIN<br>DISCOVERY Teal Act – 16 MIN   |
|                                      | DISCOVERY Green HRP        | DISCOVERY Green H <sub>2</sub> O <sub>2</sub> – 16–32 MIN<br>DISCOVERY Green Act – 16 MIN |
| HRP<br>(fluorescent)                 | DISCOVERY DCC Kit          | 32 MIN  |
|                                      | DISCOVERY FAM Kit          | 20 MIN  |
|                                      | DISCOVERY FITC Kit         | 20 MIN  |
|                                      | DISCOVERY Rhodamine Kit    | 32 MIN  |
|                                      | DISCOVERY Rhodamine 6G Kit | 32 MIN  |
|                                      | DISCOVERY Red 610 Kit      | 32 MIN  |
|                                      | DISCOVERY Cy5 Kit          | 40 MIN  |

\*We recommend using the same IHC chromogen settings for Sequential ISH-IHC as previously optimized for IHC alone.

†Extending DISCOVERY Red incubation can result in a dot-like background

‡For stronger AP-based Red IHC detection, ChromoMap Red is recommended (select **DS Red**).

13. Select your preferred Counterstain and Post-Counterstain settings. A light counterstain is recommended for best visualization of multiplex chromogenic staining.
14. At the top of the Protocol Steps window, click **Save As**, then select a unique protocol number from the drop-down menu and choose a protocol name. Click **Active**, add relevant comments in the available field, and click **Save**.



15. Make a new protocol for each probe/antibody/chromogen combination. Click **Save**.

### Print the labels

1. Select the **Print Label** icon from the upper right corner of the home screen.
2. Select your preferred template or create a new template. To create a new template, refer to the *Ventana DISCOVERY ULTRA System User Manual* for details.
3. Click **Protocol**.
4. Select the protocols you created in the section above. Click the **Add** button. When the protocols for all slides have been assigned, click **Close/Print**.
5. Fill in the template for each slide. Click **Print** when completed

### Load the reagents


1. Remove the nozzle caps from the filled dispensers and place each cap on the post located on the back of the dispenser.
2. Prime the user-fillable dispensers. For guidance, refer to the instructions provided by Roche Tissue Diagnostics.
3. If needed, remove any air bubbles at the nozzle tip by pushing down on the nozzle until the liquid reaches the tip of the nozzle or forms a small meniscus at the tip of the nozzle.
4. Remove the yellow locking ring from the dispensers in all the prefilled dispensers. Refer to the instructions provided by Roche Tissue Diagnostics.
5. Load the dispensers onto the reagent racks.
6. Load the reagent racks onto the reagent carousel.
7. Select the **Ready** button.

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.






8. Open the slide drawers.
9. Load each slide onto a heater pad with the label facing upward and inward. Ensure that the slides sit securely on the pads.

---

**IMPORTANT!** Prior to loading the slides, ensure heater pads are completely dry. Wipe off any liquid using laboratory tissue paper.

---

10. Close the slide drawers.
11. Select the **Running** button.



12. The assay duration varies based on assay selections, approximately **15 – 20 HRS**.

---

**IMPORTANT!** Before leaving the instrument unattended, ensure all reagents and slides are successfully registered and the instrument is running.

---

### *Complete the run*

1. After the run is complete, remove the Dewax reagent, place nozzle cap on the dispenser, and store at room temperature.
2. For the remaining reagents, place nozzle caps back on the dispensers and place racks onto magnet locking tray.

---

**IMPORTANT!** Store reagent racks at **4°C** until next use. Store the Dewax dispenser at room temperature.

---

### *Wash and dry the slides*

1. Prepare 200 mL of diluted detergent by adding 1–2 drops detergent to 200 mL distilled water in a container with a cap.
2. Mix well by inverting the container 4–5 times.
3. Add diluted detergent to a Tissue-Tek Staining Dish.

**Note:** Store diluted detergent at **RT**.

4. Submerge a Tissue-Tek Slide Rack into the Tissue-Tek Staining Dish containing 200 mL diluted detergent.
5. Open the instrument slide drawers and unload slides.
6. Decant solution on the slides into the slide drawer, then *immediately* load slides into the Tissue-Tek Slide Rack submerged in detergent.
7. Rinse oil off the slides by moving the slide rack up and down in the dish 10 times.
8. Replace the detergent with distilled water and rinse slides by moving the slide rack up and down a minimum of **10** times.
9. Repeat Step 8 three to five times.
10. Transfer the slides into a Tissue-Tek Staining Dish containing **200 mL** distilled water.
11. Place slides in a drying oven at **60°C** for at least **30 MIN**.

*Mount the samples*

1. In a fume hood, fill two clearing agent dishes with ~200 mL fresh xylene.
2. Once slides are dry, move the Tissue-Tek Slide rack into the staining dish containing xylene for **1 MIN** with occasional agitation.
3. Move the Tissue-Tek Slide rack into the staining dish containing xylene for **1 MIN** with occasional agitation.
4. Lay each slide flat with the sections facing up in the fume hood then add 1–2 drops of EcoMount or other chromogen-compatible xylene-based mounting medium. Carefully place a 24 mm x 50 mm coverslip over the section and avoid trapping air bubbles.
5. Air dry slides for at least **15 MIN** before evaluation.

---

**IMPORTANT!** mRNA Teal, mRNA Green, DISCOVERY Teal HRP and DISCOVERY Green HRP chromogens are light sensitive and can fade over time. For best results, protect stored slides from the light and image within one week of staining.

---

### Chapter 3: Integrated RNA-Protein Co-Detection Workflow (ICW)

This workflow is designed to enable broader compatibility of RNA co-detection with protease-tolerant protein epitopes. Tissue samples undergo deparaffinization and RNAscope Target Retrieval, then primary antibody incubation and post-primary fixation for stabilization prior to protease digestion. After protease digestion, use ACD's probes and ISH detection reagents to detect RNA, followed by secondary antibody detection and counterstaining.

## Materials

### ACD reagents for Integrated RNA-Protein Co-Detection

Integrated RNA-Protein Co-Detection Workflow requires the use of the following reagents from ACD for use on Roche DISCOVERY ULTRA:

- RNAscope VS Universal Sample Prep Reagent Kit v2 (Cat. No. 323740)
- ACD Assay-specific Detection Reagent Kit:
  - RNAscope VS Universal AP Detection Reagents (Cat. No. 323260)
  - RNAscope VS Universal HRP Detection Reagents (Cat. No. 323210)
  - RNAscope VS Duplex Detection Reagents (Cat. No. 323310)
  - BaseScope VS Detection Reagents (Cat. No. 323710)
- VS RNA-Protein Co-Detection Ancillary Reagent Kit (Cat. No. 323760)

RNAscope VS Universal Sample Prep Reagent Kit v2 (Cat. No. 323740) and RNAscope VS Accessory Kit (Cat. No. 320630) are RTU and stored as indicated in the following tables:

| RNAscope VS Sample Prep Reagent Kit v2 (Cat. No. 323740) |   |                  |                     |
|--|---|------------------|---------------------|
| <input checked="" type="checkbox"/>                      | Reagent                                   | Quantity         | Storage             |
|  | RNAscope VS Universal Target Retrieval v2 | 10 mL x 2 bottle | Room Temp (15–30°C) |
|  | RNAscope VS Universal Dewax               | 14 mL x 1 bottle | Room Temp (15–30°C) |

Assay-specific Detection Reagents within each Reagent Kit are stored as indicated in the following tables:

| RNAscope VS Universal AP Detection Reagents (Cat. No. 323260) |                                |                  |         |
|---|--------------------------------|------------------|---------|
| <input checked="" type="checkbox"/>                           | Reagent                        | Quantity         | Storage |
|   | RNAscope VS Universal AP AMP 1 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 2 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 3 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 4 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 5 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 6 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Universal AP AMP 7 | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope VS Protease           | 14 mL x 1 bottle | 2–8°C   |

### RNAscope VS Universal HRP Detection Reagents (Cat. No. 323210)

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.

| <input checked="" type="checkbox"/> | Reagent                         | Quantity         | Storage |
|-------------------------------------|---------------------------------|------------------|---------|
|                                     | RNAscope VS Universal HRP AMP 1 | 14 mL x 1 bottle | 2–8°C   |
|                                     | RNAscope VS Universal HRP AMP 2 | 14 mL x 1 bottle | 2–8°C   |
|                                     | RNAscope VS Universal HRP AMP 3 | 14 mL x 1 bottle | 2–8°C   |
|                                     | RNAscope VS Universal HRP AMP 4 | 14 mL x 1 bottle | 2–8°C   |
|                                     | RNAscope VS Universal HRP AMP 5 | 14 mL x 1 bottle | 2–8°C   |
|                                     | RNAscope VS Universal HRP AMP 6 | 14 mL x 1 bottle | 2–8°C   |
|                                     | RNAscope VS Universal HRP AMP 7 | 14 mL x 1 bottle | 2–8°C   |
|                                     | RNAscope VS Protease            | 14 mL x 1 bottle | 2–8°C   |

| RNAscope VS Duplex Detection Reagents (Cat. No. 323310) |                             |                   |         |
|---|-----------------------------|-------------------|---------|
| <input checked="" type="checkbox"/>                     | Reagent                     | Quantity          | Storage |
|   | RNAscope VS Duplex AMP 1    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 2    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 3    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 4    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 5    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 6    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 7    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 8    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP 9    | 14 mL x 1 bottle  | 2–8°C   |
|   | RNAscope VS Duplex AMP Wash | 14 mL x 2 bottles | 2–8°C   |
|   | RNAscope VS Protease        | 14 mL x 1 bottle  | 2–8°C   |

| BaseScope VS Detection Reagents (Cat. No. 323710) |                                       |                  |         |
|---|---------------------------------------|------------------|---------|
| <input checked="" type="checkbox"/>               | Reagent                               | Quantity         | Storage |
|   | BaseScope VS AMP 1                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 2                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 3                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 4                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 5                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 6                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 7                    | 14 mL x 1 bottle | 2–8°C   |
|   | BaseScope VS AMP 8                    | 14 mL x 1 bottle | 2–8°C   |
|   | RNAscope 2.5 VS Pretreat 3 - Protease | 14 mL x 1 bottle | 2–8°C   |

## VS RNA-Protein Co-Detection Ancillary Kit - Cat No. 323760

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.

| <input checked="" type="checkbox"/> | Reagent  | Source / Ordering Info | Quantity          | Storage |
|-------------------------------------|--|------------------------|-------------------|---------|
|                                     | Co-Detection Antibody Diluent*                     | ACD/Cat No. 323160     | 120 mL x 1 bottle | 2–8°C   |
|                                     | VS RNA-Protein Co-Detection Protease and Inhibitor | ACD/Cat No. 323190     | See the following |         |
|                                     | VS Co-Detection Protease                           |                        | 14 mL x 1 bottle  | 2–8°C   |
|                                     | VS Co-Detection Inhibitor                          |                        | 14 mL x 1 bottle  | 2–8°C   |

\*Use of Co-Detection Antibody Diluent is recommended for best preservation of RNA and maintenance of ISH sensitivity when performing Integrated RNA-Protein Co-Detection

### Optional ACD reagents for Integrated RNA-Protein Co-Detection:

The following additional components from ACD are optional for the Integrated RNA-Protein Co-Detection Workflow on Roche DISCOVERY ULTRA:

- RNAscope VS Accessory Kit (Cat. No. 320630)

| RNAscope VS Accessory Kit (Cat. No. 320630) |                            |                 |         |  |
|---|----------------------------|-----------------|---------|--|
| <input checked="" type="checkbox"/>         | Reagent                    | Quantity        | Storage |  |
|   | RNAscope VS Hematoxylin    | 7 mL x 1 bottle | 2–8°C   |  |
|   | RNAscope VS Bluing Reagent | 7 mL x 1 bottle | 2–8°C   |  |

### Roche materials for Integrated RNA-Protein Co-Detection:

Integrated RNA-Protein Co-Detection may be performed with either Roche RTU primary antibody or your choice of primary antibody concentrate diluted in the ACD Co-Detection Antibody Diluent. For a list of available Roche RTU Primary Antibodies and ordering information, please contact your local Roche representative.

The following additional materials from Roche can be used for Integrated RNA-Protein Co-Detection on DISCOVERY ULTRA. Catalog numbers are valid in the United States only. For other regions, please check catalog or ordering numbers with your local lab supplier.

### Reagent dispensers

| Roche Materials for Integrated RNA-Protein Co-Detection |                                 |                    |                                    |  |
|---|---------------------------------|--------------------|------------------------------------|--|
| <input checked="" type="checkbox"/>                     | Component                       | Cat. No.           | Ordering Code                      | Fill with:   |
|   | Probe Dispensers                | 960-761 to 960-780 | Contact local Roche representative | RNAscope VS Probes   |
|   | mRNA Sample Prep Kit            | 760-248            | 08127166001                        | RNAscope VS Sample Prep Reagent Kit v2 and VS Protease reagent from Detection Kit  |
|   | Antibody Dispensers (Optional*) | 770-001 to 770-099 | Contact local Roche representative | User-sourced Primary Antibody Concentrate diluted in Co-Detection Antibody Diluent |
|   | Enzyme 1 Dispenser              | 771-721            | 05271517001                        | VS Co-Detection Protease   |
|   | Pretreatment 1 Dispenser        | 960-901            | 05280095001                        | VS Co-Detection Inhibitor  |
|   | Counterstain 1 dispenser        | 771-741            | 05271720001                        | VS Hematoxylin   |
|   | Counterstain 2 dispenser        | 771-742            | 05271738001                        | VS Bluing Reagent  |
|   | Roche RTU Counterstain          | Various            | Contact local Roche                | Pre-filled   |

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.

| Roche Materials for Integrated RNA-Protein Co-Detection |   |                          |                |            |
|---|---|--------------------------|----------------|------------|
| <input checked="" type="checkbox"/>                     | Component                               | Cat. No.                 | Ordering Code  | Fill with: |
|   | Reagents                                |                          | representative |            |
|   | DISCOVERY Antibody Block                | 760-4204                 | 05268869001    | Pre-filled |
|   | Additional ISH Assay-Specific Materials | see the following tables |                |            |
|   | Additional IHC Detection Materials      | see the following tables |                |            |

\*Protein target qualification can be performed with either Roche RTU primary antibody or your choice of primary antibody.

### Additional ISH assay-specific materials

| Additional Roche Materials Required for VS Universal AP ISH |                                  |          |               |   |
|---|----------------------------------|----------|---------------|---|
| <input checked="" type="checkbox"/>                         | Component                        | Cat. No. | Ordering Code | Fill with:  |
|   | mRNA RED Probe Amplification Kit | 760-236  | 07095341001   | RNAscope VS Universal AP Detection Reagents AMP 1-7 |
|   | mRNA RED Detection Kit           | 760-234  | 07099037001   | Pre-filled  |

| Additional Roche Materials Required for VS Universal HRP ISH |                                |          |               |  |
|--|--------------------------------|----------|---------------|--|
| <input checked="" type="checkbox"/>                          | Component                      | Cat. No. | Ordering Code | Fill with:   |
|  | mRNA Probe Amplification Kit   | 760-222  | 06614337001   | RNAscope VS Universal HRP Detection Reagents AMP 1-7 |
|  | DISCOVERY Inhibitor*           | 760-4840 | 07017944001   | Pre-filled   |
|  | mRNA Purple HRP Detection Kit* | 760-255  | 08127166001   | Pre-filled   |
|  | mRNA Green HRP Detection Kit†  | 760-278  | 08952612001   | Pre-filled   |
|  | mRNA Teal HRP Detection Kit†   | 760-256  | 08352941001   | Pre-filled   |
|  | mRNA DAB Detection Kit†        | 760-224  | 06614353001   | Pre-filled   |
|  | DISCOVERY DCC Kit†             | 760-240  | 07988192001   | Pre-filled   |
|  | DISCOVERY FAM Kit†             | 760-243  | 07988150001   | Pre-filled   |
|  | DISCOVERY FITC Kit†            | 760-232  | 07259212001   | Pre-filled   |
|  | DISCOVERY Rhodamine Kit†       | 760-233  | 07259883001   | Pre-filled   |
|  | DISCOVERY Rhodamine 6G Kit†    | 760-244  | 07988168001   | Pre-filled   |
|  | DISCOVERY Red 610 Kit†         | 760-245  | 07988176001   | Pre-filled   |
|  | DISCOVERY Cy5 Kit†             | 760-238  | 07551215001   | Pre-filled   |

\* DISCOVERY Inhibitor is not required if using mRNA DAB Detection Kit

† Choose one mRNA Detection Kit for RNAscope HRP detection

| Additional Roche Materials Required for VS Duplex ISH |                        |          |               |   |
|---|------------------------|----------|---------------|---|
| <input checked="" type="checkbox"/>                   | Component              | Cat. No. | Ordering Code | Fill with:  |
|   | mRNA Duplex Amp Kit    | 760-249  | 08127174001   | RNAscope VS Duplex Detection Reagents AMP 1-9, AMP Wash |
|   | mRNA Link (Pre-filled) | 760-6014 | 08127115001   | Pre-filled  |
|   | mRNA RED Detection Kit | 760-234  | 07099037001   | Pre-filled  |

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.

| Additional Roche Materials Required for VS Duplex ISH |                               |          |               |            |
|---|-------------------------------|----------|---------------|------------|
| <input checked="" type="checkbox"/>                   | Component                     | Cat. No. | Ordering Code | Fill with: |
|   | mRNA Green HRP Detection Kit* | 760-278  | 08952612001   | Pre-filled |
|   | mRNA Teal HRP Detection Kit*  | 760-256  | 08352941001   | Pre-filled |
|   | mRNA DAB Detection Kit*       | 760-224  | 06614353001   | Pre-filled |

\* Choose one mRNA Detection Kit for VS Duplex Channel 1 detection

| Additional Roche Materials Required for VS BaseScope ISH |                                  |          |               |   |
|--|----------------------------------|----------|---------------|---|
| <input checked="" type="checkbox"/>                      | Component                        | Cat. No. | Ordering Code | Storage                                 |
|  | mRNA RED Probe Amplification Kit | 760-236  | 07095341001   | BaseScope VS Detection Reagents AMP 1-7 |
|  | Option 8 dispenser               | 771-758  | 05271916001   | BaseScope VS AMP 8 reagent              |
|  | mRNA RED Detection Kit           | 760-234  | 07099037001   | Pre-filled                              |

### *Additional IHC assay-specific materials*

| Reagent Options for Ventana IHC AP Detection |                                      |          |               |         |
|--|--------------------------------------|----------|---------------|---------|
| <input checked="" type="checkbox"/>          | Component                            | Cat. No. | Ordering Code | Storage |
|  | DISCOVERY UltraMap anti-Ms Alk Phos* | 760-4312 | 05269687001   | 2-8°C   |
|  | DISCOVERY UltraMap anti-Rb Alk Phos* | 760-4314 | 05269709001   | 2-8°C   |
|  | DISCOVERY Yellow Kitt                | 760-239  | 07698445001   | 2-8°C   |
|  | DISCOVERY Red Kit                    | 760-228  | 07425333001   | 2-8°C   |
|  | DISCOVERY ChromoMap Red Kit*         | 760-160  | 05266653001   | 2-8°C   |

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.



| Reagent Options for Ventana IHC HRP Detection |                                  |          |               |         |
|---|----------------------------------|----------|---------------|---------|
| <input checked="" type="checkbox"/>           | Product                          | Cat. No. | Ordering Code | Storage |
|   | DISCOVERY UltraMap anti-Ms HRP*  | 760-4313 | 05269695001   | 2–8°C   |
|   | DISCOVERY UltraMap anti-Rb HRP*  | 760-4315 | 05269717001   | 2–8°C   |
|   | DISCOVERY UltraMap anti-Rat HRP* | 760-4456 | 05891884001   | 2–8°C   |
|   | DISCOVERY UltraMap anti-Gt HRP*  | 760-4648 | 06607241001   | 2–8°C   |
|   | DISCOVERY Purple Kit*            | 760-229  | 07053983001   | 2–8°C   |
|   | DISCOVERY Green HRP Kit †        | 760-278  | 07053983001   | 2–8°C   |
|   | DISCOVERY Teal HRP Kit †         | 760-247  | 08254338001   | 2–8°C   |
|   | DISCOVERY ChromoMap DAB Kit †    | 760-159  | 05266645001   | 2–8°C   |
|   | DISCOVERY DCC Kit †              | 760-240  | 07988192001   | 2–8°C   |
|   | DISCOVERY FAM Kit †              | 760-243  | 07988150001   | 2–8°C   |
|   | DISCOVERY FITC Kit†              | 760-232  | 07259212001   | 2–8°C   |
|   | DISCOVERY Rhodamine Kit †        | 760-233  | 07259883001   | 2–8°C   |
|   | DISCOVERY Rhodamine 6G Kit †     | 760-244  | 07988168001   | 2–8°C   |
|   | DISCOVERY Red 610 Kit†           | 760-245  | 07988176001   | 2–8°C   |
|   | DISCOVERY Cy5 Kit †              | 760-238  | 07551215001   | 2–8°C   |

\* Choose one secondary detection antibody depending on the primary antibody species and desired IHC chromogen

† Choose one IHC chromogen

### Instrument buffers

| <input checked="" type="checkbox"/> | Component                | Cat. No. | Ordering Code | Storage             |
|-------------------------------------|--------------------------|----------|---------------|---------------------|
|                                     | 10X DISCOVERY Wash (RUO) | 950-510  | 7311079001    | Room Temp (15–30°C) |
|                                     | ULTRA LCS (Predilute)    | 650-210  | 5424534001    | Room Temp (15–30°C) |
|                                     | SSC Buffer (10X)         | 950-110  | 5353947001    | Room Temp (15–30°C) |
|                                     | Reaction Buffer (10X)    | 760-107  | 5266262001    | Room Temp (15–30°C) |
|                                     | DISCOVERY CC1            | 950-500  | 6414575001    | Room Temp (15–30°C) |

### User-Sourced Materials

| <input checked="" type="checkbox"/> | Description                       | Supplier          | Cat. No.  |
|-------------------------------------|-----------------------------------|-------------------|-----------|
|                                     | 10% Neutral Buffered Formalin     | SigmaAldrich      | HT5011    |
|                                     | Primary Antibody (RTU)            | User              | Various   |
|                                     | Primary Antibody Concentrate      | User              | Various   |
|                                     | SuperFrost Plus Slides (required) | Fisher Scientific | 12-550-15 |

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.

| <input checked="" type="checkbox"/> | Description   | Supplier                            | Cat. No.     |
|-------------------------------------|---|-------------------------------------|--------------|
|                                     | ProLong Gold Antifade Reagent                                 | Life Technologies                   | P36930       |
|                                     | EcoMount mounting medium                                      | Biocare                             | EM897L       |
|                                     | Tissue-Tek Vertical 24 Slide Rack                             | American Master Tech Scientific/MLS | LWSRA24      |
|                                     | Tissue-Tek Staining Dishes                                    | American Master Tech Scientific/MLS | LWT4457EA    |
|                                     | Cover Glass 24 x 50 mm  | Fisher Scientific/MLS               | 12-545-F     |
|                                     | Distilled water   | MLS                                 | —            |
|                                     | Mild liquid dishwashing detergent (Dawn detergent or similar) | MLS                                 | —            |
|                                     | Drying oven, capable of holding temperature at 60 +/- 1°C     | MLS                                 | —            |
|                                     | Fume hood   | MLS                                 | —            |
|                                     | 100% ethanol (EtOH)   | MLS                                 | —            |
|                                     | Xylene  | MLS                                 | —            |
|                                     | Tissue-Tek Clearing Agent Dishes, xylene-resistant            | American Master Tech Scientific/MLS | LWT4456EA    |
|                                     | Optional: Glass beaker (1 or 2 L)                             | MLS                                 | —            |
|                                     | Optional: Hot plate   | Fisher Scientific/MLS               | 11-300-49SHP |

## Assay Procedure

### Prepare the DISCOVERY ULTRA

#### Prepare the instrument

If the instrument has not been used for > 1 week, follow guidelines for instrument maintenance from Roche Tissue Diagnostics. Before use, empty the waste carboys if needed.

#### Dilute instrument bulk reagents

1. Prepare the instrument bulk fluids according to the manufacturer's instructions.
2. Fill bulk solution containers for 1X DISCOVERY Wash, ULTRA LCS (Predilute), and CC1 (Predilute) to be at least half full. Fully fill bulk solution containers for 2X SSC and 1X Reaction Buffer.

---

**IMPORTANT!** Do not use expired reagents.

---

#### Register new reagents

Reagent dispensers come with appropriate barcode labels and registration buttons for dispensing RNAscope VS Universal Reagents. Refer to the *Ventana DISCOVERY ULTRA System User Manual* for

details. To register reagents:

- Log all ACD reagents and probes into the software as **log user-fillable reagents** and **log user-fillable probes**, respectively.
- Use the reagent registration wand that comes with the instrument to register new reagent kits from Roche Tissue Diagnostics

### *Prepare user-fillable reagents for RNA-Protein Co-Detection*

Refer to the table on pages 6–7 to determine the proper dispenser for each reagent.

---

**IMPORTANT!** Avoid cross contamination between reagents. Dewax must be warmed to room temperature and be completely in solution before use.

---

1. Transfer the entire volume of each AMP component of the Detection Kit to the corresponding labeled dispenser from the appropriate mRNA Amplification kit (see pages 6-7 for details)
2. Fill the mRNA Sample Prep Kit:
  - a. Transfer the contents of both bottles of VS Target Retrieval v2 from the RNAscope VS Sample Prep Reagent Kit v2 to the mRNA Target Retrieval Dispenser.
  - b. Transfer the VS Dewax reagent from the RNAscope VS Sample Prep Reagent Kit v2 to the mRNA Dewax Dispenser.
  - c. Transfer the VS Protease from the Detection Kit to the mRNA Protease Dispenser
3. Fill the User Fillable Dispensers:
  - a. Transfer the RNAscope 2.5 or BaseScope VS Target Probe and control probes to the correspondingly Probe dispensers.
  - b. Transfer the VS Co-Detection Protease to the Enzyme 1 dispenser.
  - c. Transfer the VS Co-Detection Inhibitor to the Pretreatment 1 dispenser.
  - d. Transfer the VS Hematoxylin and VS Bluing to the Counterstain 1 and Counterstain 2 dispensers.
  - e. In a fume hood, carefully transfer 10% Neutral Buffered Formalin (NBF) to the Fixative 1 dispenser. For best results, use freshly aliquoted 10% NBF for each run.

---


**IMPORTANT!** Keep the Fixative 1 dispenser filled with 10% NBF capped when not on the instrument.

---


4. Follow the dispenser product insert instructions to properly prime and handle the dispensers.
5. Store tightly capped dispensers (except the Dewax dispenser) at 4°C when not in use.
6. Store tightly capped mRNA Dewax dispenser at room temperature when not in use.

### *Create an instrument protocol*

1. Open the NexES software and click the **Protocol** button.
2. Click **Create/Edit Protocols**, go to the Procedure drop-down menu and select **mRNA Universal**. Main protocol selections appear as shown:



3. Select the appropriate pretreatment conditions as shown in the following screenshot. After the main step selections, drop-down menus become available for further selection. For Cell Conditioning, we recommend starting with the same conditions previously used on your tissue for ISH-only staining.



4. To enable Co-Detection, select **Integrated Co-Detection Workflow**. Select Antibody barcode for primary antibody and Fixative 1 barcode for 10% NBF. Recommended incubation times for ICW Primary Antibody and post-primary fixation incubation is **36 MIN** and **32 MIN**, respectively.

☒ **Integrated Co-Detection Workflow**

☐ **Antibody Manual Application**

[ Note: Recommended Primary Antibody incubation time for ICW is 36min ]

✓ Apply One Drop of [ **ANTIBODY 1** ] ( **Antibody** ), Apply Coverslip, and Incubate for [ **36 Minutes** ]

[ Note: 10% NBF uses Fixative barcode ]

[ Note: Recommended time for Post-Primary Fixation for 32min ]

✓ Apply Three Drops of [ **FIXATIVE 1** ] ( **Post-Antibody Fixative** ), Apply Coverslip, and Incubate for [ **32 Minutes** ]

[ Only select one mRNA option as multiple selections will yield negative results ]

**Note:** To preserve RNA quality during primary antibody incubation, we recommend using primary antibody concentrate diluted in Co-Detection Antibody Diluent and that the primary antibody incubation occurs at room temperature. You may need to titrate the primary antibody concentration to strengthen protein detection. If titration by manual antibody application is desired instead of by automated dispenser, select **Antibody Manual Application**.

☒ **Integrated Co-Detection Workflow**

☒ **Antibody Manual Application**

✓ Hand Apply ( **Primary Antibody** ), and Incubate for [ **32 Minutes** ]

[ Note: 10% NBF uses Fixative barcode ]

[ Note: Recommended time for Post-Primary Fixation for 32min ]

✓ Apply Three Drops of [ **FIXATIVE 1** ] ( **Post-Antibody Fixative** ), Apply Coverslip, and Incubate for [ **32 Minutes** ]

[ Only select one mRNA option as multiple selections will yield negative results ]

5. Select one mRNA Detection assay of interest.
  - For BaseScope VS, select **BaseScope**.
  - For RNAscope VS Duplex, select **mRNA Duplex**.
  - For RNAscope VS Universal AP, select **mRNA AP Detection**.
  - For RNAscope VS Universal HRP, select **mRNA HRP Detection**.

Select only one mRNA option per protocol. Choosing multiple selections yields negative results.

[ Only select one mRNA option as multiple selections will yield negative results ]

☐ **BaseScope**

☐ **mRNA Duplex**

☐ **mRNA AP Detection**

☐ **mRNA HRP Detection**

☐ **Dual Sequence**

☐ **Counterstain**

☐ **Slide Cleaning**

6. Select the appropriate Detection Inhibitor:
  - For **BaseScope** or **mRNA AP Detection**, mRNA Red Inhibitor reagent is applied automatically.
  - For **mRNA Duplex**, select the **Inhibitor** check box to display the available selections. Choose the appropriate inhibitor selection based on the chromogen combination in use for Duplex ISH detection, **Red-DAB Inhibitor**, **Red-Green Inhibitor**, or **Red-Teal Inhibitor**.



[ Only select one mRNA option as multiple selections will yield negative results ]

☐ BaseScope

☒ mRNA Duplex

[ For mRNA Duplex AP / HRP Detection ]

☒ Inhibitor

☐ Red-DAB Inhibitor

☐ Red-Green Inhibitor

☐ Red-Teal Inhibitor

☐ mRNA Duplex 3rd Pretreatment

☐ mRNA Duplex ICW Protease

[ Target Probe Cocktail ]

[ RECOMMENDED: Set to 43 C for probe incubation ]

☒ Apply Two Drops of [ Probe ] ( Probe #1 ), Apply Coverslip, and Incubate for 4 Minutes

☒ Warmup Slide to [ Low Temperature ], and Incubate for 2 Hours ( Hybridization )

- For **mRNA HRP Detection**, select mRNA HRP Detection Inhibitor if using the mRNA DAB Detection Kit. Select DISCOVERY Inhibitor if using an alternative HRP-based chromogen such as the mRNA Purple Detection Kit, mRNA Green Detection Kit, or mRNA Teal Detection Kit.

[ Only select one mRNA option as multiple selections will yield negative results ]

☐ BaseScope

☐ mRNA Duplex

☐ mRNA AP Detection

☒ mRNA HRP Detection

[ Select an Inhibitor NOTE: DISCOVERY Inhibitor to be used with fluorescence ]

☒ mRNA HRP Detection Inhibitor

[ Inhibitor for mRNA DAB Detection will be applied ]

☐ DISCOVERY Inhibitor

☐ HRP detection 3rd Pretreatment

- Under the mRNA assay of interest, select the corresponding **ICW Protease** check box. Do not select the detection 3rd Pretreatment option.
- Under **ICW Protease** select the **Enzyme 1** barcode to correspond to the open dispenser filled with VS Co-Detection Protease reagent. Recommended incubation is 37°C for 16 MIN.

**IMPORTANT!** Incubating VS Co-Detection Protease above 37°C is not recommended and may negatively impact protein detection.

[ Only select one mRNA option as multiple selections will yield negative results ]

☐ BaseScope

☐ mRNA Duplex

☒ mRNA AP Detection

[ Inhibitor for mRNA Red Detection will be applied ]

☐ AP detection 3rd Pretreatment

☒ AP Detection ICW Protease

[ RECOMMENDED: Set to 37 C and 16min for normal FFPE samples ]

[ ICW Protease ]

☒ Apply Two Drops of [ ENZYME 1 ] ( Enzyme ), Apply Coverslip, and Incubate for 0 Hr 4 Min

☒ Warmup Slide to [ 37 Deg C ], and Incubate for [ 0 Hr 16 Min ] ( Pretreatment #3 Temp RB )



**IMPORTANT!** For the ICW protocol to start, the complete mRNA Sample Prep dispenser set including the mRNA Protease dispenser must be present on the reagent carousel even if the 3<sup>rd</sup> **Pretreatment** selection is not checked. We do not recommend filling the mRNA Protease dispenser with VS Co-Detection Protease reagent as the instrument cannot distinguish it from RNAscope VS Protease reagent and this could lead to misapplication of the incorrect protease reagent when using other protocols. Instead, use the **ICW Protease** selection and Enzyme 1 dispenser as described above.

9. Select the ISH staining conditions:
  - If running **BaseScope**, select the following ISH staining conditions:

[ Only select one mRNA option as multiple selections will yield negative results ]

☒ **BaseScope**

[ BaseScope Utilizes Standard Pretreatment ]

☐ BaseScope 3rd Pretreatment

☒ **BaseScope ICW Protease**

[ RECOMMENDED: Set to 37 C and 16min for normal FFPE samples ]

[ ICW Protease ]

✓ Apply Two Drops of [ **ENZYME 1** ] ( **Enzyme** ), Apply Coverslip, and Incubate for **0 Hr 4 Min**

✓ Warmup Slide to [ **37 Deg C** ], and Incubate for [ **0 Hr 16 Min** ] ( **Pretreatment #3 Temp RB** )

[ Specify Probe, hybridization time and temperature ]

✓ Apply Two Drops of [ **PROBE 1** ] ( **Probe #1** ), Apply Coverslip, and Incubate for **4 Minutes**

[ BaseScope Probe - Temp = 50 C Recommended ]

✓ Warmup Slide to [ **50 Deg C** ], and Incubate for **2 Hours** ( **Hybridization** )

[ BaseScope Amp 1 - Temp = 49 C Recommended ]

✓ Warmup Slide to [ **49 Deg C** ], and Incubate for **32 Minutes** ( **Detection #1 Temp** )

[ BaseScope Amp 2 - Temp = 49 C Recommended ]

✓ Warmup Slide to [ **49 Deg C** ], and Incubate for **32 Minutes** ( **Detection #2 Temp** )

[ BaseScope Amp 4 - Temp = 49 C Recommended ]

✓ Warmup Slide to [ **49 Deg C** ], and Incubate for **32 Minutes** ( **Detection #3 Temp** )

[ BaseScope Amp 5 - Temp = 39 C Recommended ]

✓ Warmup Slide to [ **39 Deg C** ], and Incubate for **32 Minutes** ( **Detection #4 Temp** )

[ BaseScope Amp 7 - Recommended 4min ]

✓ Incubate for [ **0 Hr 4 Min** ]

☐ mRNA Duplex

☐ mRNA AP Detection

☐ mRNA HRP Detection

☐ Dual Sequence

| Standard Temperatures/Times for BaseScope VS |              |
|--|--------------|
| ICW Protease temperature and time            | 37°C, 16 MIN |
| Standard BaseScope probe temperature         | 50°C         |
| Standard BaseScope AMP 1 temperature         | 49°C         |
| Standard BaseScope AMP 2 temperature         | 49°C         |
| Standard BaseScope AMP 4 temperature         | 49°C         |
| Standard BaseScope AMP 5 temperature         | 39°C         |
| Standard BaseScope AMP 7 incubation time*    | 4 MIN        |
| mRNA Red Chromogen                           | Default†     |

\* BaseScope Amp 7 incubation time is determined by instrument calibration. Use the instrument setting previously optimized for the mRNA Universal software. For assistance, consult your local ACD FAS for more information.

†For all BaseScope assays, mRNA Red Detection is applied automatically. No chromogen selections are needed.

- If running **mRNA Duplex**, select the following ISH staining conditions. In this example, selections are shown for a Red-Green chromogen combination:

☒ **Inhibitor**

☐ Red-DAB Inhibitor

☒ Red-Green Inhibitor

☐ Red-Teal Inhibitor

☐ mRNA Duplex 3rd Pretreatment

☒ mRNA Duplex ICW Protease

[ RECOMMENDED: Set to 37 C and 16min for normal FFPE samples ]

[ ICW Protease ]

✓ Apply Two Drops of [ [ENZYME 1](#) ] ( **Enzyme** ), Apply Coverslip, and Incubate for **0 Hr 4 Min**

✓ Warmup Slide to [ [37 Deg C](#) ], and Incubate for [ **0 Hr 16 Min** ] ( **Pretreatment #3 Temp RB** )

[ Target Probe Cocktail ]

[ RECOMMENDED: Set to 43 C for probe incubation ]

✓ Apply Two Drops of [ [PROBE 1](#) ] ( **Probe #1** ), Apply Coverslip, and Incubate for **4 Minutes**

✓ Warmup Slide to [ [43 Deg C](#) ], and Incubate for **2 Hours** ( **Hybridization** )

[ Amp 1 Duplex ]

[ RECOMMENDED: Temp = 39 C for most samples ]

✓ Warmup Slide to [ [39 Deg C](#) ], and Incubate for **32 Minutes** ( **Hybridization #1** )

[ Amp 2 Duplex ]

[ RECOMMENDED: Temp = 39 C for most samples ]

✓ Warmup Slide to [ [39 Deg C](#) ], and Incubate for **32 Minutes** ( **Hybridization #3** )

[ Amp 5 incubation time: RECOMMENDED is 4 min ]

✓ Incubate for [ [0 Hr 4 Min](#) ] ( **Hybridization #5** )

[ Amp 8 incubation time: RECOMMENDED: 4 min ]

✓ Incubate for [ [0 Hr 4 Min](#) ] ( **Hybridization #6** )

☐ Red-DAB AP-HRP Duplex

☐ Red-Teal AP-HRP Duplex

☒ Red-Green AP-HRP Duplex

[ Recommended Settings are 20 min for both Green H202 and Green Activator ]

✓ Apply One Drop of mRNA Green H202, and Incubate for [ [20 Minutes](#) ]

| Standard Temperatures/Times for mRNA Duplex Detection  |  |
|--|--|
| ICW Protease temperature and time                      | 37°C, 16 MIN   |
| Suggested probe temperature                            | 43°C   |
| Suggested RNAscope VS Duplex AMP 1 & AMP 2 temperature | 39°C   |
| RNAscope VS Duplex AMP 5 incubation time               | 4 MIN*   |
| RNAscope VS Duplex AMP 8 incubation time               | 4 MIN*   |
| Chromogen settings for Red-DAB AP-HRP Duplex†          | Preset, no incubation selections   |
| Chromogen settings for Red-Green AP-HRP Duplex†        | mRNA Green H <sub>2</sub> O <sub>2</sub> — 20 MIN<br>mRNA Green Act — 20 MIN |
| Chromogen settings for Red-Teal AP-HRP Duplex†         | mRNA Teal H <sub>2</sub> O <sub>2</sub> — 20 MIN<br>mRNA Teal Act — 20 MIN   |

\* VS Duplex Amp 5 and Amp8 incubation times are determined by instrument calibration. Use the instrument setting previously optimized for the mRNA Universal software. For assistance, consult your local ACD FAS for more information.

†Choose one detection combination for mRNA Duplex. All times/temperatures are preset for Red-DAB Detection.

- If choosing **mRNA AP Detection**, select the following ISH staining conditions:

[ Only select one mRNA option as multiple selections will yield negative results ]

☐ BaseScope

☐ mRNA Duplex

☒ mRNA AP Detection

[ Inhibitor for mRNA Red Detection will be applied ]

☐ AP detection 3rd Pretreatment

☒ AP Detection ICW Protease

[ RECOMMENDED: Set to 37 C and 16min for normal FFPE samples ]

[ ICW Protease ]

✓ Apply Two Drops of [ ENZYME 1 ] ( Enzyme ), Apply Coverslip, and Incubate for 0 Hr 4 Min

✓ Warmup Slide to [ 37 Deg C ], and Incubate for [ 0 Hr 16 Min ] ( Pretreatment #3 Temp RB )

[ Target Probe: RECOMMENDED: Set to 43 C for probe incubation ]

✓ Apply Two Drops of [ PROBE 1 ] ( Probe #1 ), Apply Coverslip, and Incubate for 4 Minutes

✓ Warmup Slide to [ 43 Deg C ], and Incubate for 2 Hours ( Hybridization )

[ Amp 1 AP ]

[ RECOMMENDED: Temp = 39 C for most samples ]

✓ Warmup Slide to [ 39 Deg C ], and Incubate for 32 Minutes ( Hybridization #2 )

[ Amp 2 AP ]

[ RECOMMENDED: Temp = 39 C for most samples ]

✓ Warmup Slide to [ 39 Deg C ], and Incubate for 32 Minutes ( Hybridization #4 )

[ Amp 5 AP incubation time: RECOMMENDED: 4 min ]

✓ Incubate for [ 0 Hr 4 Min ] ( Hybridization #5 )

☐ mRNA HRP Detection

☐ Dual Sequence

| Standard Temperatures/Times for mRNA AP Detection |              |
|---|--------------|
| ICW Protease temperature and time                 | 37°C, 16 MIN |
| Suggested probe temperature                       | 43°C         |
| Suggested RNAscope AP AMP 1 & AMP 2 temperature   | 39°C         |
| RNAscope AP AMP 5 incubation time                 | 4 MIN*       |
| Chromogen settings for mRNA Red Detection         | Default†     |

\* VS Universal AP Amp 5 incubation time is determined by instrument calibration. Use the instrument setting previously optimized for the mRNA Universal software. For assistance, consult your local ACD FAS for more information.

†For all mRNA AP, mRNA Red Detection is applied automatically. No chromogen selections are needed.

- If choosing **mRNA HRP Detection**, select the following ISH staining conditions:



[ Only select one mRNA option as multiple selections will yield negative results ]

☐ BaseScope

☐ mRNA Duplex

☐ mRNA AP Detection

☒ mRNA HRP Detection

[ Select an Inhibitor NOTE: DISCOVERY Inhibitor to be used with fluorescence ]

☒ mRNA HRP Detection Inhibitor

[ Inhibitor for mRNA DAB Detection will be applied ]

☐ DISCOVERY Inhibitor

☐ HRP detection 3rd Pretreatment

☒ HRP Detection ICW Protease

[ RECOMMENDED: Set to 37 C and 16min for normal FFPE samples ]

[ ICW Protease ]

✓ Apply Two Drops of [ ENZYME 1 ] ( Enzyme ), Apply Coverslip, and Incubate for 0 Hr 4 Min

✓ Warmup Slide to [ 37 Deg C ], and Incubate for [ 0 Hr 32 Min ] ( Pretreatment #3 Temp RB )

[ Target Probe ]

[ Target Probe: RECOMMENDED: Set to 43 C for probe incubation ]

✓ Apply Two Drops of [ PROBE 1 ] ( Probe #1 ), Apply Coverslip, and Incubate for 4 Minutes

✓ Warmup Slide to [ 43 Deg C ], and Incubate for 2 Hours ( Hybridization )

[ Amp 1 HRP ]

[ RECOMMENDED: Temp = 39 C for most samples ]

✓ Warmup Slide to [ 39 Deg C ], and Incubate for 32 Minutes ( Hybridization #5 )

[ Amp 2 HRP ]

[ RECOMMENDED: Temp = 39 C for most samples ]

✓ Warmup Slide to [ 39 Deg C ], and Incubate for 32 Minutes ( Hybridization #6 )

[ Amp 5 HRP incubation time: RECOMMENDED: 4 min ]

✓ Incubate for [ 0 Hr 4 Min ] ( Hybridization #6 )

[ Default detection is mRNA DAB unless a fluor or chromogen is selected ]

☐ Rhodamine

☐ FITC

☐ Rhodamine 6G

☐ Cy5

☐ FAM

☐ Red 610

☐ DCC

☐ mRNA Purple

☐ mRNA Teal

☐ mRNA Green


☐ Dual Sequence

| Standard Temperatures/Times for mRNA HRP Detection |  |
|--|--|
| ICW Protease temperature and time                  | 37°C, 16 MIN   |
| Suggested probe temperature                        | 43°C   |
| Suggested RNAscope HRP AMP 1 & AMP 2 temperature   | 39°C   |
| RNAscope HRP AMP 5 incubation time                 | 4 MIN*   |
| Chromogen settings for mRNA DAB†                   | Default*   |
| Chromogen settings for mRNA Purple                 | mRNA Purple H <sub>2</sub> O <sub>2</sub> — 40 MIN                           |
| Chromogen settings for mRNA Teal                   | mRNA Teal H <sub>2</sub> O <sub>2</sub> — 20 MIN<br>mRNA Teal Act — 20 MIN   |
| Chromogen settings for mRNA Green                  | mRNA Green H <sub>2</sub> O <sub>2</sub> — 20 MIN<br>mRNA Green Act — 20 MIN |

\* VS Universal HRP Amp 5 incubation time is determined by instrument calibration. Use the instrument setting previously optimized for the mRNA Universal software. For assistance, consult your local ACD FAS for more information.

†For mRNA HRP, mRNA DAB Detection is applied by default if no other chromogens are selected.


10. After completing ISH assay selections, select **Dual Sequence** to enable selections for secondary antibody detection.



- ☒ **Dual Sequence**
- ☐ Antibody Denaturation
- ☐ DS Pretreatment
- ☐ DS Option
- ☐ DS A/B Block
- ☐ DS Inhibitor
- [ DISC Inhibitor will not be applied ]
- ☐ DS ISH
- ☐ DS Antibody
- ☐ Disable heat for 2nd Fixative
- ☐ DS Post-Antibody Fixative
- ☐ DS Antibody Amplification
- ☐ DS Linking Antibody
- ☐ Disable heat for 2nd Detection
- ☐ DS Multimer HRP
- ☐ DS Multimer AP
- ☐ DS Proximity Detection
- ☐ DS Enzyme conjugate
- ☐ DS DISCOVERY Amplification
- ☐ DS DAB
- ☐ DS Silver
- ☐ DS DISCOVERY Purple
- ☐ DS DISCOVERY Yellow HRP
- ☐ DS DISCOVERY Teal HRP
- ☐ DS DISCOVERY Green HRP
- ☐ DS DISCOVERY Red HRP
- ☐ DS DISCOVERY Blue HRP
- ☐ DS Red
- ☐ DS DISCOVERY Red
- ☐ DS DISCOVERY Yellow
- ☐ DS Blue
- ☐ DS Rhodamine
- ☐ DS FITC
- ☐ DS Cy5
- ☐ DS DCC
- ☐ DS FAM
- ☐ DS Red 610
- ☐ DS Rhodamine 6G
- ☐ DS Open Detection Kit
- ☐ Triple Stain
- ☐ Counterstain




To apply VS Co-Detection Inhibitor, select **DS Pretreatment**, **DS 2nd Pretreatment**, and **Use DW for DS 2nd Pretreatment**. Select the Pretreatment 1 barcode to correspond to the VS Co-Detection Inhibitor reagent. Recommended incubation is **40°C for 32 MIN** prior to protein detection.



11. Select the settings for secondary detection:

**Note:** We recommend using Roche's Antibody Block reagent directly prior to secondary application, to prevent non-specific cross-detection of RNAscope. RNA-Protein Co-Detection can be performed without Antibody Block, but this could result in a hue shift of the RNA dots and interfere with interpretation of RNA-protein colocalization.

- For **HRP-based secondary detection**, choose a Roche Multimer HRP reagent corresponding to the primary host species. To enable secondary incubation, select **DS Multimer HRP**. To apply Antibody Block, select **DS Multimer HRP Blocker**, then select Antibody Block reagent. Select your Multimer HRP reagent of choice and desired incubation time, using previously optimized IHC-only conditions.



- For **AP-based secondary detection**, choose a Roche Multimer AP reagent corresponding to the primary host species. To enable secondary incubation, select **DS Multimer AP**. To apply Antibody Block, select **DS Multimer AP Blocker**, then select Antibody Block reagent. Select your Multimer AP reagent of choice and desired incubation time, using previously optimized IHC-only conditions.

☐ Disable heat for 2nd Detection  
☐ DS Multimer HRP  
☒ DS Multimer AP  
☒ DS Multimer AP Blocker  
     [ Select Multimer blocker and Multimer species ]  
     [ Note: Recommended Multimer AP Reagent incubation time for ICW is 32min ]  
     ✓ Apply One Drop of [ [Antibody Block](#) ] ( DS Multimer AP Blocking #1 ), No Coverslip and Incubate for 32 Minutes  
     ✓ Apply One Drop of [ [UMap anti-Ms AP](#) ] ( DS Multimer AP ), and Incubate for [ [32 Minutes](#) ]  
     [ Select Multimer blocker and Multimer species ]  
☐ DS Proximity Detection  
☐ DS Enzyme conjugate  
☐ DS DISCOVERY Amplification  
☐ DS DAB

12. Select chromogens for IHC detection using the following recommendations:

|   |   |
|---|---|
| <input type="checkbox"/> DS Enzyme conjugate<br><input type="checkbox"/> DS DISCOVERY Amplification<br><input checked="" type="checkbox"/> DS DAB<br><input type="checkbox"/> DS Silver<br><input checked="" type="checkbox"/> DS DISCOVERY Purple<br><input type="checkbox"/> DS DISCOVERY Yellow HRP<br><input type="checkbox"/> DS DISCOVERY Teal HRP<br><input type="checkbox"/> DS DISCOVERY Green HRP<br><input type="checkbox"/> DS DISCOVERY Red HRP<br><input type="checkbox"/> DS DISCOVERY Blue HRP<br><input type="checkbox"/> DS Red<br><input type="checkbox"/> DS DISCOVERY Red<br><input type="checkbox"/> DS DISCOVERY Yellow<br><input type="checkbox"/> DS Blue<br><input type="checkbox"/> DS Rhodamine<br><input type="checkbox"/> DS FITC<br><input type="checkbox"/> DS Cy5<br><input type="checkbox"/> DS DCC<br><input type="checkbox"/> DS FAM<br><input type="checkbox"/> DS Red 610<br><input type="checkbox"/> DS Rhodamine 6G<br><input type="checkbox"/> DS Open Detection Kit | <div style="border-left: 1px solid black; padding-left: 10px; margin-left: 10px;"> <p><b>HRP Chromogens</b></p> </div> <div style="border-left: 1px solid black; padding-left: 10px; margin-left: 10px;"> <p><b>AP Chromogens</b></p> </div> <div style="border-left: 1px solid black; padding-left: 10px; margin-left: 10px;"> <p><b>DISCOVERY HRP Fluorescent Kits</b></p> </div> |
|---|---|

| IHC Chromogen Selection and Settings |                            |   |
|--------------------------------------|----------------------------|---|
| IHC Enzyme                           | IHC Chromogen              | Recommended Chromogen Settings*   |
| AP                                   | DISCOVERY Yellow           | 44 MIN – 2 HRS  |
|                                      | DISCOVERY Red              | 12 MIN†   |
|                                      | DS Red / ChromoMap Red‡    | Default   |
| HRP<br>(chromogenic)                 | DS DAB / ChromoMap DAB     | Default   |
|                                      | DISCOVERY Purple           | 40 MIN  |
|                                      | DISCOVERY Teal HRP         | DISCOVERY Teal H <sub>2</sub> O <sub>2</sub> – 16–32 MIN<br>DISCOVERY Teal Act – 16 MIN   |
|                                      | DISCOVERY Green HRP        | DISCOVERY Green H <sub>2</sub> O <sub>2</sub> – 16–32 MIN<br>DISCOVERY Green Act – 16 MIN |
| HRP<br>(fluorescent)                 | DISCOVERY DCC Kit          | 32 MIN  |
|                                      | DISCOVERY FAM Kit          | 20 MIN  |
|                                      | DISCOVERY FITC Kit         | 20 MIN  |
|                                      | DISCOVERY Rhodamine Kit    | 32 MIN  |
|                                      | DISCOVERY Rhodamine 6G Kit | 32 MIN  |
|                                      | DISCOVERY Red 610 Kit      | 32 MIN  |
|                                      | DISCOVERY Cy5 Kit          | 40 MIN  |

\*Protein detection with ICW can benefit from longer chromogen incubation or selection of a stronger chromogen relative to IHC alone


†Extending DISCOVERY Red incubation can result in a dot-like background

‡For stronger AP-based Red IHC detection, ChromoMap Red is recommended (select **DS Red**).

- Select your preferred Counterstain and Post-Counterstain settings.

Note: Hematoxylin staining can appear darker with ICW than with ISH-alone. Try a lighter counterstain if needed for optimal multiplex chromogenic co-detection.

- At the top of the Protocol Steps window, click **Save As**, then select a unique protocol number from the drop-down menu and choose a protocol name. Click **Active**, add relevant comments in the available field, and click **Save**.



15. Make a new protocol for each probe/antibody/chromogen combination. Save as a unique protocol number.

### Print the labels

1. Select the **Print Label** icon from the upper right corner of the home screen.
2. Select your preferred template or create a new template. To create a new template, refer to the *Ventana DISCOVERY ULTRA System User Manual* for details.
3. Click **Protocol**.
4. Select the protocols you created in the section above. Click the **Add** button. When the protocols for all slides have been assigned, click **Close/Print**.
5. Fill in the template for each slide. Click **Print** when completed


### Load the Reagents

1. Remove the nozzle caps from the filled dispensers and place each cap on the post located on the back of the dispenser.
2. Prime the user-fillable dispensers. For guidance, refer to the instructions provided by Roche Tissue Diagnostics.
3. If needed, remove any air bubbles at the nozzle tip by pushing down on the nozzle until the liquid reaches the tip of the nozzle or forms a small meniscus at the tip of the nozzle.
4. Remove the yellow locking ring from the dispensers in all the prefilled dispensers. Refer to the instructions provided by Roche Tissue Diagnostics.
5. Load the dispensers onto the reagent racks.
6. Load the reagent racks onto the reagent carousel.
7. Select the **Ready** button.

DISCOVERY ULTRA RNA Protein Co-Detection

MK 51-174/REV B/Effective date: 06/15/2022

For Research Use Only. Not for use in diagnostics procedures.




8. Open the slide drawers.
9. Load each slide onto a heater pad with the label facing upward and inward. Ensure that the slides sit securely on the pads.

---

**IMPORTANT!** Prior to loading the slides, ensure heater pads are completely dry. Wipe off any liquid using laboratory tissue paper.

---

10. Close the slide drawers.
11. Select the **Running** button.



12. The assay duration varies based on assay selections, approximately **15 – 20 HRS**.

---

**IMPORTANT!** Before leaving the instrument unattended, ensure all reagents and slides are successfully registered and the instrument is running.

---

### *Complete the run*

1. After the run is complete, remove the Dewax reagent, place nozzle cap on the dispenser, and store at room temperature.
2. For the remaining reagents, place nozzle caps back on the dispensers and place racks onto magnet locking tray.

---

**IMPORTANT!** Store reagent racks at **4°C** until next use. Store the Dewax dispenser at room temperature.

---

### *Wash and dry the slides*

1. Prepare 200 mL of diluted detergent by adding 1–2 drops detergent to 200 mL distilled water in a container with a cap.
2. Mix well by inverting the container 4–5 times.
3. Add diluted detergent to a Tissue-Tek Staining Dish.

**Note:** Store diluted detergent at **RT**.

4. Submerge a Tissue-Tek Slide Rack into the Tissue-Tek Staining Dish containing 200 mL diluted detergent.
5. Open the instrument slide drawers and unload slides.
6. Decant solution on the slides into the slide drawer, then *immediately* load slides into the Tissue-Tek Slide Rack submerged in detergent.
7. Rinse oil off the slides by moving the slide rack up and down in the dish 10 times.
8. Replace the detergent with distilled water and rinse slides by moving the slide rack up and down a minimum of **10** times.
9. Repeat Step 8 three to five times.
10. Transfer the slides into a Tissue-Tek Staining Dish containing **200 mL** distilled water.
11. Place slides in a drying oven at **60°C** for at least **30 MIN**.

*Mount the samples*

1. In a fume hood, fill two clearing agent dishes with ~200 mL fresh xylene.
2. Once slides are dry, move the Tissue-Tek Slide rack into the staining dish containing xylene for **1 MIN** with occasional agitation.
3. Move the Tissue-Tek Slide rack into the staining dish containing xylene for **1 MIN** with occasional agitation.
4. Lay each slide flat with the sections facing up in the fume hood then add 1–2 drops of EcoMount or other chromogen-compatible xylene-based mounting medium. Carefully place a 24 mm x 50 mm coverslip over the section and avoid trapping air bubbles.
5. Air dry slides for at least **15 MIN** before evaluation.

---

**IMPORTANT!** mRNA Teal, mRNA Green, DISCOVERY Teal HRP and DISCOVERY Green HRP chromogens are light sensitive and can fade over time. For best results, protect stored slides from the light and image within one week of staining.

---



## Appendix A. Sequential ISH-IHC Staining Optimization and Troubleshooting

Prior to using the sequential workflow for RNA-Protein Co-Detection, we recommend establishing working protocols for both ISH-only and IHC-only staining in the mRNA Universal procedure, which uses a different deparaffinization and cell conditioning method than the RUO Universal procedure. Perform protein target qualification to confirm that your protein epitope target is protease-tolerant by following the protocol setup outlined in **Chapter 1. Protein Target Qualification**. Many epitopes are sensitive to RNAscope protease and are not compatible with sequential staining. If your epitope is protease-sensitive, refer to **Chapter 3. Integrated Co-Detection Workflow**.

For Sequential Co-Detection, begin by applying the same baking, deparaffinization, cell conditioning, and protease conditions as the ISH-only protocol. Apply the same primary and secondary antibody, and IHC chromogen conditions as the IHC-only protocol. IHC optimization varies by tissue and primary antibody clone. If needed, the following parameters can be adjusted:

| Reagent                           | User Selection                                  | Recommended Incubation Temperature | Recommended Incubation Time | Optimization Range   |
|-----------------------------------|---|------------------------------------|-----------------------------|--|
| Target Retrieval                  | Cell Conditioning                               | 97°C                               | 24 MIN                      | 85–97°C<br>8–56 MIN  |
| RNAscope VS Protease              | 3 <sup>rd</sup> Pretreatment                    | 37°C                               | 16 MIN                      | 37–50°C<br>8–32 MIN  |
| Primary Antibody                  | DS Antibody                                     | 37°C                               | 32 MIN                      | Range available by user selection:<br>DS Antibody Default:<br><ul style="list-style-type: none"> <li>35–42°C, 16–60 MIN</li> </ul> DS Disable Heat:<br><ul style="list-style-type: none"> <li>Room Temperature, 4m–24h</li> </ul> DS High Temp Ab Incubation:<br><ul style="list-style-type: none"> <li>60–75°C, 4–60 MIN</li> </ul> DS Extended Antibody Incubation:<br><ul style="list-style-type: none"> <li>4–120 MIN</li> </ul> |
| Manually Applied Primary Antibody | DS Antibody –<br>DS Antibody Manual Application | 37°C                               | 32 MIN                      | Range available by user selection:<br>Default:<br><ul style="list-style-type: none"> <li>35–42°C, 16–60 MIN</li> </ul> DS Disable Heat:<br><ul style="list-style-type: none"> <li>Room Temperature, 4m–24h</li> </ul> DS High Temp Ab Incubation:<br><ul style="list-style-type: none"> <li>60–75°C, 4–60 MIN</li> </ul>   |
| Secondary Antibody / Multimer*    | DS Multimer AP<br>-or-<br>DS Multimer HRP       | Not Selectable                     | 32 MIN                      | 4–32 MIN   |
| IHC Chromogen                     |   | See recommendations on page 39     |                             |  |

\* For stronger detection, try an alternative Roche secondary system.

## Appendix B. ICW Staining Optimization and Troubleshooting

For optimal RNA and protein detection, we recommend establishing ISH-only and IHC-only protocols before proceeding to RNA-Protein Co-Detection. Begin by applying the same baking, deparaffinization, and cell conditioning conditions as the ISH-only protocol. Since a stronger protease treatment is required with Co-Detection than with ISH alone, we recommend using VS ICW Protease in most cases.

To optimize protein detection, titrate the antibody concentration for the co-detection workflow. You may need to use a higher primary antibody concentration for the co-detection workflow than you would normally use for IHC alone. For best preservation of RNA in the Integrated Co-Detection Workflow, we recommend diluting primary antibody concentrate in the ACD VS Co-Detection Antibody diluent. If further optimization is needed for a specific sample or target of interest, the following parameters can be adjusted:

| Reagent                               | User Selection   | Recommended Incubation Temperature | Recommended Incubation Time | Optimization Range  |
|---------------------------------------|--|------------------------------------|-----------------------------|---------------------|
| Target Retrieval                      | Cell Conditioning  | 97°C                               | 24 MIN                      | 85-97°C<br>8-56 MIN |
| Primary Antibody                      | Integrated Co-Detection Workflow                               | Not Selectable                     | 36 MIN                      | 16-60 MIN           |
| Primary Antibody (Manual Application) | Integrated Co-Detection Workflow / Manual Antibody Application | Not Selectable                     | 32 MIN                      | 16-60 MIN           |
| 10% NBF                               | Integrated Co-Detection Workflow                               | Not Selectable                     | 32 MIN                      | 16-60 MIN           |
| VS Co-Detection Protease              | ICW Protease   | 37°C                               | 16 MIN                      | 37°C*<br>8-32 MIN   |
| Secondary Antibody / Multimer**       | DS Multimer AP<br>-or-<br>DS Multimer HRP                      | Not Selectable                     | 32 MIN                      | 4-32 MIN            |
| IHC Chromogen                         |  | See recommendations on page 62     |                             |                     |

\*For best ISH sensitivity with ICW, do not increase VS Co-Detection Protease beyond 37°C

\*\* For stronger detection, try an alternative Roche secondary system.

† Some samples, such as cell pellets, can require reduced protease treatment. For these samples, use RNAscope VS Protease by selecting **3<sup>rd</sup> Pretreatment** rather than **ICW Protease**. If a sample requires stronger protease treatment than the maximum 32 minutes of VS Co-Detection Protease, we recommend two protease reagent applications. Select both **3<sup>rd</sup> Pretreatment** and **ICW Protease**. Heating either protease beyond 37°C during the integrated workflow is not recommended because this can result in decreased ISH sensitivity.



**For Research Use Only. Not for use in diagnostic procedures.**

Advanced Cell Diagnostics, Inc. and/or its affiliate(s) warrant their products as set forth in the ACD General Terms and Conditions of Sale found on the ACD website. Advanced Cell Diagnostics, Inc. reserves the right to change its products and services at any time to incorporate technological developments. This document is subject to change without notice. Although this document has been prepared with every precaution to ensure accuracy, Advanced Cell Diagnostics, Inc. assumes no liability for any errors, omissions, or for any damages resulting from the use of this information.

© 2022 **Advanced Cell Diagnostics**. All rights reserved. RNAscope is a registered trademark of Advanced Cell Diagnostics, Inc. All other trademarks belong to their respective owners.

**Headquarters**

7707 Gateway Blvd Suite 200, Newark, CA 94545 Phone 1-510-576-8800 Toll Free 1-877-576-3636 For support, email [support.acd@bio-technie.com](mailto:support.acd@bio-technie.com)  
[www.acdbio.com](http://www.acdbio.com)

