

miRNAscope™ LS Reagent Kit – RED

Description	The miRNAscope™ LS Reagent Kit - RED (Cat.No.324600) is designed to work with the miRNAscope™ LS Target Probe – New Target Probe, for fully automated <i>in situ</i> detection of RNA in tissue specimens on slides using the Leica BioSystems BOND Rx instruments. Each kit contains sufficient reagents in a convenient, Ready-To-Use (RTU) format to perform 60 slide/tests. The results of the assay can be viewed using a standard brightfield microscope.				
Contents and Storage	The components of the miRNAscope™ LS Reagent Kit - RED and quantities supplied are listed in the following table. The contents have a shelf life until the date of expiration indicated on the reagent label. It is important to follow the product storage conditions as specified on the label so that the reagents maintain their reliability.				
	Cat. No.	Sub-Kit	Component	Quantity	Storage
	324600	miRNAscope™ LS Reagent Kit – RED	miRNAscope™ LS AMP 1	21 mL x 1 bottle	2–8°C
			miRNAscope™ LS AMP 2	21 mL x 1 bottle	2–8°C
			miRNAscope™ LS AMP 3	21 mL x 1 bottle	2–8°C
			miRNAscope™ LS AMP 4 - RED	21 mL x 1 bottle	2–8°C
			miRNAscope™ LS AMP 5 - RED	21 mL x 1 bottle	2–8°C
			miRNAscope™ LS AMP 6 - RED	21 mL x 1 bottle	2–8°C
			RNAscope® 2.5 LS Protease III	21 mL x 1 bottle	2–8°C
			RNAscope® 2.5 LS Rinse	29 mL x 1 bottle	15–30°C
RNAscope® 2.5 LS Bluing Reagent			21 mL x 1 bottle	15–30°C	
RNAscope® 2.5 LS Hydrogen Peroxide			21 mL x 1 bottle	15–30°C	
Precautions	IMPORTANT! Before running the assay, please download and review the <i>miRNAscope™ LS Reagent Kit-RED User Manual</i> (Doc. No. UM-324600) at https://acdbio.com/technical-support/user-manuals .				
	Handle biological specimens and other materials according to local laboratory guidelines and regulations. Avoid cross-contamination of samples and reagents to ensure quality results. Use equipment and reagents as specified by the manufacturer. For safety and biohazard guidelines, refer to the “Safety” section in the user manual. For every chemical, read the Material Safety Data Sheets (MSDSs) and follow the handling instructions available at https://acdbio.com/documents/product-documents . Wear appropriate protective eyewear, clothing, and gloves.				
Contact Information	Address: 7707 Gateway Blvd CA 94545 USA Phone: 1-877-576-3636 option 1 Fax: 1-510-576-8801 Web: www.acdbio.com Email: support.acd@bio-techne.com				

Copyright
©2020. Advanced Cell Diagnostics, Inc. All rights reserved.

Trademarks
miRNAscope and RNAscope are trademarks or registered trademarks of Advanced Cell Diagnostics, Inc. Leica BioSystems and BOND Rx are trademarks of Leica BioSystems. All other trademarks belong to their respective owners. All other trademarks belong to their respective owners. Advanced Cell Diagnostics, Inc. reserves the right to change its products and services at any time to incorporate technological developments. This insert is subject to change without notice. Although this insert has been prepared with every precaution to ensure accuracy, Advanced Cell Diagnostics, Inc. assumes no liability for any errors, omissions, or for any damages resulting from the use of this information.