

READY, SET, GO! GETTING STARTED WITH RNASCOPE

Presented by:
Jacqueline Akech, Ph.D.
April 14th, 2015

Senior Scientist

Advanced Cell Diagnostics

©2014 Advanced Cell Diagnostics, Inc. | Confidential and Proprietary | For Research Use Only (RUO), not intended for diagnosis.



TOPICS


- How Does RNAscope® Work?
- Getting Started with RNAscope in your Laboratory
- Tips and Tricks on Running the Assay
- Frequently Asked Questions
- Time for Q&A



RNASCOPE OVERVIEW

RNASCOPE WORKFLOW


PERMEABILIZE
cells or tissue




HYBRIDIZE
to target RNA



AMPLIFY
signal




VISUALIZE
with morphology



QUANTIFY
single-cell expression




RNASCOPE[®] TECHNOLOGY: ZZ PROBE DESIGN






mRNA transcript detection: Highly specific & robust signal amplification

RNASCOPE PUBLICATION AND LITERATURE REFERENCES



Wang, F. et al. 2012







Visit www.acdbio.com and download a publication of your interest




RNASCOPE WORKFLOW

RNASCOPE ASSAY SELECTION




RNAscope Assays	RNAscope 2.0 HD BROWN	RNAscope 2.0 HD RED	RNAscope 2-plex	RNAscope Multiplex – Fluorescence
Assay type	Chromogenic	Chromogenic	Chromogenic	Fluorescent
Dye Used	Diaminobenzene (DAB)-HRP	Fast Red -ALP	HRP-Green, Fast Red -ALP	FITC, Cy3, Cy5,
Channel	Channel 1	Channel 1	Channel 1, 2	Channel 1, 2, 3
Probes channel	C1 Probes	C1 Probes	C1, C2 Probes	C1, C2, C3 Probes

RNASCOPE WORKFLOW: CHROMOGENIC ASSAY



TIP : Detection protocols will vary based on the chromogenic assay used
Download manuals: <http://www.acdbio.com/technical-support/downloads>

RNASCOPE WORKFLOW: FLUORESCENT ASSAY



TIP : Pretreatment conditions will vary based on sample type

Download manuals: <http://www.acdbio.com/technical-support/downloads>



ONE DAY OR TWO DAY ASSAY?

ONE DAY ASSAY

Sample preparation



Sample pretreatment



RNAscope assay

TWO DAY ASSAY

Sample preparation



Sample pretreatment




RNAscope assay

TIP : Review the User Manuals PART 1 and PART 2 for optional stopping points



GETTING STARTED WITH RNASCOPE IN YOUR LAB

GET STARTED BY FOLLOWING 3 EASY STEPS



TIP : Visit www.acdbio.com/go for more information on getting started

THE CHECKLIST: WHAT YOU NEED



<input checked="" type="checkbox"/>	RNAscope target probes (single-plex and multiplex)
<input checked="" type="checkbox"/>	Positive and Negative control probe
<input checked="" type="checkbox"/>	Hot-Plate for pretreatment/ target retrieval step
<input checked="" type="checkbox"/>	RNAscope reagent kit (Detection Kit, Pretreatment kits & Wash Buffer)
<input checked="" type="checkbox"/>	HybEZ Hybridization system
<input checked="" type="checkbox"/>	EZ Batch slide processing system or Tissue-Tek system
<input checked="" type="checkbox"/>	Ecomount for 2.0 HD Red & 2-plex chromogenic assay
<input checked="" type="checkbox"/>	RNAscope control slides
<input checked="" type="checkbox"/>	Immedge hydrophobic barrier pen
<input checked="" type="checkbox"/>	User supplied reagents (refer to user manual)
<input checked="" type="checkbox"/>	Read the user manual (Part 1 – Sample Prep & Part2 – Detection Assay)

***TIP : Visit www.acdbio.com/go for more information on getting started.
Checklist is available on the website and in the manual***

USING A HOT PLATE

✓ **Hotplate for retrieval/boiling**



TIP : When using a hot plate for pre-treatment step – pay close attention to the TIME and boiling TEMPERATURE

RNASCOPE REAGENT KIT CONTENTS



OLD



NEW



Contents of the reagent kit

1. Pretreatment reagents
2. RNAscope detection kit
3. Wash buffer

TIP : Warm probes at 40 °C for 10 minutes before use

TIP : Warm 50x wash buffer at 40 °C for 20 minutes if you notice a precipitation



USING A HYBEZ HYBRIDIZATION OVEN



✓ **HyBEZ hybridization system**




TIP: HybEZ oven is required as it provides both temperature and humidity control, necessary to obtain optimal RNAscope results

ACCESSORIES FOR WASHING STEPS



Tissue Tek washing tray

EZ Batch for slide processing





TIP : ACD EZ Batch slide processing tray is easy and convenient for loading multiple slides for washing steps.




QUALIFY YOUR SAMPLES USING CONTROLS

RNASCOPE® SEMI-QUANTITATIVE SCORING



Score	Criteria
0	No staining or <1 dot/ 10 cells*
1	1-3 dots/cell
2	4-9 dots/cell. None or very few dot clusters
3	10-15 dots/cell and <10% dots are in clusters
4	>15 dots/cell and >10% dots are in clusters




Score = 3





Score = 3


QUALIFY YOUR SAMPLES USING CONTROLS



QC check
PPIB > 2
DapB < 1



QC check
PPIB > 2
DapB < 1



PASS



Run your target probes

FAIL




- *Verify technique*
- *Check RNA quality with new samples*
- *Perform Assay optimization*

Technique check

Sample/RNA quality check


TIP : Always start with standard conditions

OPTIMIZE YOUR ASSAY



Technique check

Sample/
RNA quality check



OPTIMIZE YOUR ASSAY →


TIP : Refer to the Troubleshooting Guide





OPTIMIZE YOUR ASSAY

OPTIMIZE YOUR SAMPLE IN 3 EASY STEPS



STEP 1 START WITH STANDARD CONDITIONS



Observe Staining Pattern -


High background, over-digested? = **underfixed**

No signal/weak signal, under-digested? = **overfixed**

STEP 2 ADJUST PRETREATMENT 2, BOILING TIME



STEP 3 ADJUST PRETREATMENT 3/4, PROTEASE TIME*





TIP: For cultured cells, protease is diluted 1:15 in 1X PBS

* For fresh frozen samples, only protease pretreatment is required and is performed at room temperature



TROUBLESHOOTING: OVERDIGESTION



Sample: FFPE Xenograft

Assay: RNAscope 2.0 HD Red

Issue: Destroyed morphology, ghost nuclei, high nuclear background, weak hematoxylin staining


Optimization: Decrease pretreatment 2 conditions.

Result: Strong staining for positive controls with no/negligible background, intact morphology, strong hematoxylin staining


TIP: Visit <http://www.acdbio.com/technical-support/downloads/rnascope-ish-guide-troubleshooting/>

EXAMPLE OF SUCCESSFUL RNASCOPE RESULTS


Negative control, DapB



Positive control, PPIB



Target probe



RNAscope 2.0 HD Red

Human breast cancer tissue



TROUBLESHOOTING TIPS

FACTORS AFFECTING RNASCOPE ASSAY PERFORMANCE

<input checked="" type="checkbox"/>	Fixation conditions are not optimal
<input checked="" type="checkbox"/>	RNA is degraded
<input checked="" type="checkbox"/>	Hybridization conditions not optimal
<input checked="" type="checkbox"/>	Samples drying during assay




THE SOLUTIONS

<input checked="" type="checkbox"/>	Fix samples as recommended. E.g., for FFPE use 10% NBF RT, 16-32 hrs
<input checked="" type="checkbox"/>	Acquire new samples and assess RNA quality
<input checked="" type="checkbox"/>	Use the HybEZ hybridization oven only
<input checked="" type="checkbox"/>	Use Immedge pen and add adequate reagents to avoid drying


IMPACT OF FIXATION CONDITIONS

24 hours fixation/**Optimal**



Sample: FFPE brain sample

3 weeks fixation/**Over fixed**




Assay: RNAscope 2.0 Brown

Synaptophysin


TIP: Sample fixation has a great effect on the success of your assay
Solution: Increase pretreatment for better target accessibility

TROUBLESHOOTING: UNDER FIXATION

Positive control, Rn PPIB



Positive control, Rn PPIB



Sample: Flash Frozen followed by FFPE sample preparation (fixation), Rat intestines

Assay: RNAscope 2.0 HD Brown

Issue: Weak staining, destroyed morphology, FFPE sample is under fixed

Optimization: Fixation according to recommended guidelines for FFPE samples


Result: Strong staining for positive control, PPIB, intact morphology

TIP : Refer to the Troubleshooting Guide

<http://www.acdbio.com/technical-support/downloads/rnascope-ish-guide-troubleshooting/>

TROUBLESHOOTING: ASSAY WORKFLOW

Negative control, dapB



Sample: FFPE Hela pellet


Assay: RNAscope 2.0 HD Brown

Issue: Tissue dried out, high background

Optimization: Do not allow drying between amplification steps. Use the Immedge hydrophobic barrier pen

Result: Clean background

Negative control, dapB



TIP: Refer to the Troubleshooting Guide;

<http://www.acdbio.com/technical-support/downloads/rnascope-ish-guide-troubleshooting/>

REFER TO SAMPLE PRETREATMENT GUIDELINE

Tissue Pretreatment Guidelines



Follow the recommended pretreatment conditions based on your tissue type for:

- ❖ Any new or previously untested FFPE tissues
- ❖ Samples prepared suboptimally

Guidelines for Optimal Tissue Pretreatment

- ❑ Test representative samples with positive and negative control probes. [Controls should be: Positive = uniform signal; negative = blank].

- ❑ Fix sample in **FRESH 10% NBF for 16–32 HOURS at ROOM TEMPERATURE.**

NOTE: Do not fix at 4°C. DO NOT fix for < 16 hrs or >32 hrs. Refer to Table 1 for under/over-fixed tissue pretreatment guidelines.

- ❑ Vary **PRETREAT 2** and/or **PRETREAT 3 TIME** based on your tissue type (see Table 2).

NOTE: Certain Xenografts and Cell Pellets, require very mild pretreatment (**PRETREAT 2** for 8 min, **PRETREAT 3** for 15 min).

Table 1. Tissue Pretreatment Guidelines

Reagent	Mild	Standard	Extended
Pretreat 2	15 min	15 min	30 min
Pretreat 3	15 min	30 min	30 min

Table 2. Tissue Pretreatment Table

Species	Tissue type	Pathology	Pretreat Condition	Species	Tissue type	Pathology	Pretreat Condition	
Mouse / Rat	Intestine	Normal	Standard	Human	Cervical	Normal	Standard	
	Intestine	Tumor	Standard		Cervical dysplasia	Abnormal	Standard	
	Embryo	Normal	Standard		Brain	Tumor	Standard	
	Brain	Normal	Standard		Brain	Normal	Standard	
	Spleen	Normal	Mild		Head	Cancer	Standard	
	Eye/Retina	Normal	Standard		Neck	Cancer	Standard	
	Liver	Normal	Extended		Liver	Cancer	Standard	
	Kidney	Normal	Standard		Kidney	Normal	Standard	
	Human	Breast	Tumor		Standard	Skin	Normal	Standard
		Colon	Tumor		Standard	Melanoma	Tumor	Standard
Colon		Normal	Standard	Nevus	Benign	Standard		
Lung		Tumor	Standard	Placenta	Normal	Standard		
Lung		Normal	Standard	Skin (TMA*)	Normal	Standard		
Prostate		Tumor	Standard	Breast (TMA)	Normal	Standard		
Prostate		Normal	Standard	Melanoma (TMA)	Normal	Standard		
Lymph node		Tumor	Mild	Nevus (TMA)	Benign	Standard		
Lymph node		Normal	Mild	Stomach (TMA)	Normal	Standard		
Tonsil		Normal	Mild	Stomach (TMA)	Tumor	Standard		
Pancreas		Normal	Standard	Cell pellets**	—	Mild		
Cervical		Cancer	Standard	HeLa cells† (ACD control)	—	Standard		

* Tissue Microarray
 ** Fixed with 10% NBF
 † Fixed with 10% Formaldehyde/PBS

For information about species or tissue type not listed here, contact support at support@acdbio.com.

TIP : Refer to the user manual for tissue specific pretreatment guidelines





FREQUENTLY ASKED QUESTIONS

FREQUENTLY ASKED QUESTIONS

- **RNAscope assay compatibility with different tissues**

RNAscope manual assay can be used with FFPE, fresh-frozen, fixed-frozen and cultured cells. RNAscope automated assays are primarily supported with the FFPE tissue. Please refer to the User Manual Selection Guide: <http://www.acdbio.com/technical-support/downloads>

- **Key differences between RNAscope ISH assay and IHC**

No cooling is required during Epitope retrieval, users should directly put the slides in water at room temperature, dehydrate and proceed to Pretreatment 3 step as per the manual Part 1

TIP: Visit www.acdbio.com/support for additional FAQs

GUIDELINES TO FOLLOW WITH RNASCOPE ASSAY


- ✓ Read the manual and perform RNAscope exactly
- ✓ Utilize optional stopping points
- ✓ Flicking/tapping the slides for adequate drying of slides
- ✓ Storage in desiccant (FFPE) for RNA integrity
- ✓ Always use fresh reagents
- ✓ Warm probes and wash buffer at 40°C, precipitation occurs during storage
- ✓ Remember to optimize pretreatment conditions, when you switch tissues

TIP: Refer to the Troubleshooting Guide


<http://www.acdbio.com/technical-support/downloads/rnascope-ish-guide-troubleshooting/>

SUMMARY

1. Reviewed RNAscope technology
 - ZZ probe design and the workflow for manual assays
2. Tips for getting started with RNAscope in your lab
 - 3 easy steps to getting started in your lab



3. Tips for success with RNAscope assay, when you switch tissues
 - Pretreatment 2 and 3 optimization



VISIT THE SUPPORT PAGE TO LEARN MORE

ACD
Advanced Cell Diagnostics

Provider of the most sensitive and specific RNA *in situ* hybridization technology

Contact Us Quote Builder Search entire website

Technology Applications Science Diagnostics Products Services **Support** Company

**OPTIMIZE YOUR ASSAY,
OPTIMIZE YOUR RESULTS.**

Your formula for success in a few easy steps.

SUPPORT OVERVIEW

Technical Support Overview
Training Webinars
Downloads
Getting Started
Frequently Asked Questions (FAQs)
Online Training Videos
Product Literature

Technical Support Overview
Training Webinars
Download Manuals

We value your RNA research and we recognize you need the best to get superior results from our RNAscope® assays. Please start with our training videos.

For any product inquiries, please fill out our contact form to the right. For technical assistance with an RNAscope assay, please contact us.

Support you in getting the most out of your assays thanks to product tutorial

Contact Us

TIP: Visit www.acdbio.com/technical-support/support-overview



CONTACT ACD SUPPORT

- Support via email [_support@acdbio.com](mailto:support@acdbio.com)
- Support via phone-1-877-376-3636, option 3
 - Time 8:00am-6:00pm PST
- Support Resources available on website www.acdbio.com



 Manuals	 Getting Started	 FAQs	 Videos	 Product Literature
Download manuals, technical notes and MSDS.	Simple tips & tricks for you to get the best RNAscope result from day1.	Browse through our product frequently asked questions or add one of your own.	View our product and workflow videos on our Video page.	Find RNAscope publication lists, gene lists and download product brochures.
Go →	Go →	Go →	Go →	Go →

QUESTIONS?





Jacqueline Akech
ACD Technical Support
Advanced Cell Diagnostics, Inc.
3960 Point Eden Way, Hayward, CA 94545

Advanced Cell Diagnostics

©2013 Advanced Cell Diagnostics, Inc. | Confidential and Proprietary

