



ADVANCED CELL DIAGNOSTICS


TARGET PROBE: HS-HPV HR18

RNASCOPE® PROBE INFORMATION

Catalog number:	312591 (Manual)*, 312599 (2.5 VS assay)*, 312298 (2.5 LS assay)*
Species:	Virus-HPV
Target Gene:	HPV-HR18 E6/E7 mRNA
Full Name:	Human Papilloma Virus
Description:	18 high-risk types: HPV 16, 18, 26, 31, 33, 35, 39, 45, 51, 52, 53, 56, 58, 59, 66, 68, 73 and 82, E6/E7 mRNA

* Target region is the same across RNAscope® manual and automation assays

PROBE DATA



RNAscope® Manual Brown Assay performed on FFPE head and neck squamous carcinoma tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and HPV-HR18 probe (Bottom; Cat. No. 312591), 40X. Red arrow indicates representative target probe detection

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/1489078>

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/1489079>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker. Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.