


ADVANCED CELL DIAGNOSTICS

RNASCOPE® PROBE INFORMATION


Catalog number:	426581 (Manual)*, 426589 (2.5 VS assay)*, 426588 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	THBS1
Alternate names:	TSP; THBS; TSP1; TSP-1; THBS-1; THBS1
Entrez Gene ID:	7057
Accession#:	NM_003246.2
Unigene:	Hs.164226
Target Sequence:	Homo sapiens thrombospondin 1 (THBS1) mRNA
Probe Region:	3636 - 4657
ZZ pair number:	20

* Target region is the same across RNAscope® manual and automation assays

PROBE DATA



RNAscope® Manual Brown Assay performed on FFPE liver cancer tissue under standard condition, using the dapB Negative Control probe (Left; Cat. No. 310043) and Hs-THBS1 probe (Right; Cat. No. 426581), 40X



RNAscope® Manual Brown Assay performed on FFPE esophageal cancer tissue under standard condition, using the dapB Negative Control probe (Left; Cat. No. 310043) and Hs-THBS1 probe (Right; Cat. No. 426581), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/7057>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.

Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.